

# **Laboratory Water Analysis**

Photometric and Electrochemical Instruments, Reagents and Services



2 INTRODUCTION www.hach-lange.com

# HACH LANGE Photometric and Electrochemical Instruments, Reagents & Services

With high quality products, consumables, accessories and comprehensive services, HACH LANGE is your ideal partner for water analysis. Our laboratory solutions ensure accurate and reliable analysis for all key parameters in the municipal, regulatory and industrial sectors, in the lab and in the field.



#### **Everything from a single supplier**

From instrumentation out in the field or in the lab, to sampling, reagents, accessories and consumables. HACH LANGE provides all you need for your water analysis.

#### For every application

HACH LANGE water analysis is the result of decades of practical experience. We supply you with tailor-made solutions for reliable monitoring of wastewater, drinking water and process water.

#### Parameters from A-Z

From Ammonium to Zirconium.
Consistently user friendly solutions; proven in daily practice. Our customers know they can rely on HACH LANGE for their water analysis; from sample preparation to quality control. In the interest of continuous product improvement, technical specifications are subject to change without prior notice.

#### CONTENT

- 4 Photometric System
- 6 Photometric Tests
- 7 LCK Cuvette Tests
- 12 Powder Pillows
- 14 SWIFTEST
- 15 ACCUVAC
- 16 Liquid Reagent Tests
- 18 TEST'N'TUBE
- 20 Standard Solutions
- 22 Sample Preparation
- 23 Accessories
- 25 Photometers
- 26 Photometer Quick Reference Guide
- 28 AP 3900
- 30 DR 6000
- 32 DR 3900
- 34 DR 900
- 36 POCKET Colorimeter II
- 38 HT 200S Thermostat
- 39 LT 200 Thermostat
- 40 Cuvette Quick Reference Guide
- 41 Photometer Accessories
- 42 Turbidity
- 43 Electrochemistry
- 44 HQD and SENSION+ Meters and Probes Introduction
- 46 HQD Benchtop and Portable Meters
- 48 HQD INTELLICAL Probes
- 50 HQD Meter & Probe Accessories
- 52 SENSION+ Portable and Benchtop Meters
- 55 SENSION+ Portable, Benchtop and ISE Probes
- 58 SENSION+ Accessories
- 59 Pocket Pro and Pocket Pro+ Testers
- 61 pH Buffer and Conductivity Standard Solutions
- 62 HACH LANGE Services
- 65 Parameter Index
- 66 How to Order



## **Photometric System**

HACH LANGE offers a perfectly coordinated system of photometers and reagents, required accessories and services. For all key parameters from Ammonium to Zirconium. From a fast screening test to standard comparable analysis, with sample preparation and Quality Assurance.



#### Systematic quality and efficiency

Only a perfect interaction guarantees highest efficiency and accuracy — starting with the individual components of the spectrophotometer and the ready-to-use chemistry up to the interaction with you and your laboratory equipment. HACH LANGE delivers to you a perfectly coordinated system — as a developer, manufacturer and sales & service partner.

#### Easy & safe handling

By means of bar-coded cuvettes, the DR spectrophotometer automatically identifies test parameter, range, method, lot number and expiry date. Colour coded cuvettes, packaging, pictograms, and instructions in multiple languages simplify testing. DOSICAP ZIP reagent delivery provides ease of use and eliminates hazardous chemistry handling.

#### Sustainable & environmentally friendly

Continuous environmental investment is a high priority in the development of the HACH LANGE cuvette tests. Since 1978, we have collected used reagents for proper disposal. Thanks to the special reagent processing techniques applied in the HACH LANGE Environment Centre, more than 75% of all returned test components are fed back into the production and material cycles.

## Analysis made simple

#### LCK cuvette tests - unrivalled analysis

- ► Safe Maximum safety for users, thanks to the closed cuvette system and low amounts of reagents.

  Complete labelling of the individual cuvettes, including barcode label for automatic recognition in the photometer.
- ► Easy Convenient and error-free dosing of the reagents without pipetting and reagent contact, thanks to DOSICAP and DOSICAP ZIP: cuvette caps containing an exactly pre-dispensed amount of freeze-dried reagent.
- ► Approved HACH LANGE cuvette tests are officially approved for legally required consent limits.

  With the help of standard solutions and round-robin test solutions, they provide the assurance you need.
- ► Versatile 50 parameters and more than 100 measuring ranges for all applications in water analysis from extremely polluted industrial wastewater to trace analysis in drinking water.



#### IBR+ increases reliability

During the rotating ten times measurement process using the IBR+ Integrated Barcode Reader, the DR spectrophotometer immediately picks up all the information on the cuvette, also including lot number and expiry date. Both are documented with the measurement value. In case of exceeding the expiry date you automatically get an alarm.



#### RFID for traceability and rapid data updates

By means of RFID (Radio-Frequency Identification) technology you can trace your samples completely back to the sampling point. All important data like sampling point, operator, date and time are stored on an RFID tag on the sample bottle. In addition, Certificates of Analysis (CoA) can be retrieved immediately from the RFID tag on the packaging. In the laboratory, all this information will be transferred via RFID identification to the DR spectrophotometer in a matter of seconds.

Never before has updating or programming of methods into the spectrophotometer been so easy and quick. You simply hold the cuvette test box in front of the DR's RFID module, wait for the acoustic signal and that's it. The measurement starts instantly - with the correct calibration data leading to the right result.



#### **Analytical Quality Assurance (AQA)**

Quality assurance and analysis are completely interlinked. QA procedures can be easily defined and documented within the instrument without additional software. Results are only dependable in conjunction with AQA. HACH LANGE offers classic single standard solutions as well as practical multi-standard solutions in application-oriented combinations. In addition the comprehensive ADDISTA AQA system for cuvette tests contains two round-robin solutions which entitle you to participate in external round-robin tests free of charge.



#### Alignment of laboratory and process analysis

Compare your online value with your laboratory reference value directly in the spectrophotometer - via LINK2SC connection between SC controller and DR 3900/DR 6000. The exchange of data works bidirectional, which means that you can do a matrix correction of your process probe straight from the laboratory.



## **Photometric Tests**

HACH LANGE offers a family of prepared reagents designed to streamline your testing and save your time. Whether you are using Cuvette Tests, Powder Pillows, TEST'N'TUBEs, Liquid Reagents, or ACCUVACs, our reagents deliver the quality and convenience you expect from us.



#### **Focus on quality**

To make sure you receive high-quality reagents, HACH LANGE conducts stringent quality control checks throughout the entire manufacturing process. Several in-process control steps take place as well as comprehensive inspections of the finished products. Documentation of the final quality checks is carried out on our Certificates of Analysis (CoA).

#### Reagents tailored to your needs

Lot after lot, HACH LANGE prepares reagents providing accurate results, enhanced stability, repeatable results, and extended shelf life. Especially suitable for photometric measurements under difficult conditions are Powder Pillows and ACCUVACs. For demanding analysis the cuvette tests with their outstanding precision are just what is needed.

#### **Everything for your security**

It goes without saying that safety comes first. Hazard code information can be found on product labels, packages, Safety Data Sheets AND, what's more, in this catalogue! Safety at a glance: in accordance with Article 48 of the CLP regulation (EC) No. 1272/2008, the following tables provide EU and/or GHS hazard codes of our laboratory chemicals.

## Which LCK cuvette tests for my photometer?

### Quick reference guide

#### LCK - Outstanding precision and handling



Our cuvette tests cover all water analysis applications. They satisfy the most demanding tasks, e.g. monitoring consent limits as an equivalent alternative to time-consuming reference methods. The 2D barcode also details the lot number and the expiry date of the reagents. The Certificate of Analysis (CoA) is directly available via RFID tag on the packaging.

| Article<br>number | Parameter        | Measuring<br>range                     | Method                               | According to standard      | Quality<br>control | Number<br>of tests | PC II | DR<br>3900 | DR<br>6000 | EU<br>hazard<br>code | GHS<br>hazard<br>code                         |
|-------------------|------------------|--|--------------------------------------|----------------------------|--------------------|--------------------|-------|------------|------------|----------------------|---|
| LCK362            | Acid capacity    | 0.5 - 8.0 mmol/L                       | HACH LANGE Method                    |                            |                    | 25                 |       |            | -          | -                    | -   |
| LCK300            | Alcohol          | 0.01 - 0.12 g/L                        | Alcohol Oxidase (Enzymatic)          |                            |                    | 24                 |       |            |            | -                    | -   |
| LCK301            | Aluminium        | 0.02 - 0.5<br>mg/L Al                  | Chromazurol S                        |                            | LCA702             | 24                 |       | •          | •          | T                    | GHS02,<br>GHS05,<br>GHS07,<br>GHS08           |
| LCK302            | Ammonium         | 47 - 130 mg/L<br>NH <sub>4</sub> -N    | Indophenol Blue                      | ISO 7150-1, DIN 38406 E5-1 | LCA705             | 25                 |       | •          | •          | Xn, N                | GHS05,<br>GHS07,<br>GHS09                     |
| LCK303            | Ammonium         | 2 - 47 mg/L<br>NH <sub>4</sub> -N      | Indophenol Blue                      | ISO 7150-1, DIN 38406 E5-1 | LCA703             | 25                 | •     | •          | •          | Xn, N                | GHS05,<br>GHS07,<br>GHS09                     |
| LCK304            | Ammonium         | 0.015 - 2.0 mg/L<br>NH <sub>4</sub> -N | Indophenol Blue                      | ISO 7150-1, DIN 38406 E5-1 | LCA700             | 25                 | •     | •          | •          | Xn, N                | GHS05,<br>GHS07,<br>GHS09                     |
| LCK305            | Ammonium         | 1 - 12 mg/L<br>NH <sub>4</sub> -N      | Indophenol Blue                      | ISO 7150-1, DIN 38406 E5-1 | LCA704             | 25                 | •     | •          | •          | Xn, N                | GHS05,<br>GHS07,<br>GHS09                     |
| LCK390            | AOX              | 0.05 - 3.0 mg/L<br>AOX                 | Digestion +<br>Iron(III)-Thiocyanate | DIN EN ISO 9562            | LCA390             | 24                 |       | •          | •          | T, C, F              | GHS02,<br>GHS06,<br>GHS08                     |
| LCK391            | AOX              | 0.005 - 0.5 mg/L<br>AOX                | Digestion +<br>Iron(III)-Thiocyanate | DIN EN ISO 9562            | LCA390             | 12                 |       | •          | •          | T, C, F              | GHS02,<br>GHS05,<br>GHS06,<br>GHS08           |
| LCK241            | Bitter units     | ≥ 2 BU                                 | Analogous MEBAK-Method               | MEBAK II                   |                    | 25                 |       |            | •          | Xn, F, N             | GHS02,<br>GHS05,<br>GHS07,<br>GHS08,<br>GHS09 |
| LCK554            | BOD <sub>5</sub> | 0.5 - 12 mg/L O <sub>2</sub>           | Dilution Method                      | EN 1899-1                  |                    | 20                 |       | •          | •          | Xi                   | GHS05.<br>GHS07                               |
| LCK555            | BOD <sub>5</sub> | 4 - 1650 mg/L<br>O <sub>2</sub>        | Dilution Method                      | EN 1899-1                  | LCA555             | 39                 |       |            |            | Xi                   | GHS05.<br>GHS07                               |

PC II: Single Parameter Colorimeter, DR 3900: VIS Spectrophotometer, DR 6000: UV-VIS Spectrophotometer

Please note: Some methods require reagent blanks. For these, the number of tests varies.

- \*: GHS hazard code will be available in the future
- -: product is not subject to classification







# Which LCK cuvette tests for my photometer?

| Article number | Parameter                                  | Measuring<br>range                                       | Method                 | According to standard                | Quality<br>control                      | Number<br>of tests | PC II | DR<br>3900 | DR<br>6000 | EU<br>hazard | GHS<br>hazard   |
|----------------|--|--|------------------------|--------------------------------------|---|--------------------|-------|------------|------------|--------------|---|
| Hullibel       |  | Tallyo   |                        |                                      | Control                                 | OI tests           |       | 3300       | 0000       | code         | code  |
| LCK307         | Boron                                      | 0.05 - 2.50<br>mg/L B                                    | Azomethine-H           | DIN 38405-D17                        | 191442                                  | 25                 |       | -          | •          | -            | GHS07   |
| LCK308         | Cadmium                                    | 0.02 - 0.3 mg/L<br>Cd                                    | Cadion                 |                                      | LCA702                                  | 25                 |       | •          | •          | T, N         | GHS02,<br>GHS05,<br>GHS06,<br>GHS07,<br>GHS08,<br>GHS09 |
| LCK388         | Carbonate/<br>carbon<br>dioxide            | 55 - 550 mg/L<br>CO <sub>2</sub>                         | pH Indicator           |                                      |   | 25                 |       | •          | •          | -            | -   |
| LCK311         | Chloride<br>Chloride                       | 1 - 70 mg/L Cl<br>70 - 1000<br>mg/L Cl                   | Iron(III)-Thiocyanate  |                                      | LCA700,<br>LCA703,<br>LCA704,<br>LCA705 | 24                 |       |            | •          | T, C         | GHS02,<br>GHS05,<br>GHS06                               |
| LCK410         | Chlorine, free                             | 0.05 - 2.0 mg/L<br>Cl <sub>2</sub> free/ClO <sub>2</sub> | DPD                    | ISO 7393-1-2-1985,<br>DIN 38408 G4-2 | LCA310                                  | 24                 |       |            | •          | -            | GHS07   |
| LCK310         | Chlorine/<br>Ozone/<br>Chlorine<br>dioxide | 0.05 - 2.0 mg/L<br>Cl <sub>2</sub>                       | DPD                    | ISO 7393-1-2-1985,<br>DIN 38408 G4-2 | LCA310                                  | 24                 | •     | •          | •          | -            | GHS07   |
| LCK213         | Chromic acid                               | 0.5 - 5.0 g/L<br>CrO <sub>3</sub>                        | Intrinsic Baths Colour |                                      |   | 25                 |       | •          | •          | Xi           | GHS07   |
| LCK313         | Chromium                                   | 0.03 - 1.0 mg/L<br>Cr (VI)                               | Diphenylcarbazide      | EN ISO 11885,<br>DIN 38405-D24       | LCA702                                  | 25                 |       | •          | •          | Xi, Xn       | GHS05,<br>GHS07,<br>GHS08                               |
| LCS313         | Chromium,<br>trace                         | 0.005 - 0.25<br>mg/L Cr (VI)                             | Diphenylcarbazide      | EN ISO 11885,<br>DIN 38405-D24       | LCA702                                  | 25                 |       | •          | •          | Xn           | GHS05,<br>GHS07,<br>GHS08                               |
| LCI400         | COD  | 0 - 1000 mg/L<br>0 <sub>2</sub>                          | Dichromate             | ISO 15705                            | LCA703                                  | 24                 |       | •          | •          | T, C         | GHS05,<br>GHS06,<br>GHS08,<br>GHS09                     |
| LCI500         | COD  | 0 - 150 mg/L 0 <sub>2</sub>                              | Dichromate             | ISO 15705                            | LCA704                                  | 24                 |       | •          | •          | T, C         | GHS05,<br>GHS06,<br>GHS08,<br>GHS09                     |
| LCK014         | COD  | 1000 - 10000<br>mg/L 0 <sub>2</sub>                      | Dichromate             | ISO 6060-1989,<br>DIN 38409-H41-H44  | LCA705                                  | 25                 | •     | •          | •          | T, C         | GHS05,<br>GHS06,<br>GHS08,<br>GHS09                     |
| LCK1014        | COD  | 100 - 2000<br>mg/L 0 <sub>2</sub>                        | Dichromate             | ISO 6060-1989,<br>DIN 38409-H41-H44  | LCA708                                  | 25                 |       | •          | •          | T, C         | GHS05,<br>GHS06,<br>GHS08,<br>GHS09                     |
| LCK114         | COD  | 150 - 1000<br>mg/L 0 <sub>2</sub>                        | Dichromate             | ISO 6060-1989,<br>DIN 38409-H41-H44  | LCA703                                  | 25                 | •     | •          | •          | T, C         | GHS05,<br>GHS06,<br>GHS08,<br>GHS09                     |
| LCK314         | COD  | 15 - 150 mg/L<br>0 <sub>2</sub>                          | Dichromate             | ISO 6060-1989,<br>DIN 38409-H41-H44  | LCA704                                  | 25                 | •     | -          | •          | T, C         | GHS05,<br>GHS06,<br>GHS08,<br>GHS09                     |
| LCK414         | COD  | 5.0 - 60 mg/L 0 <sub>2</sub>                             | Dichromate             | ISO 6060-1989,<br>DIN 38409-H41-H44  | LCA700                                  | 24                 |       | •          | •          | T, C         | GHS05,<br>GHS06,<br>GHS08,<br>GHS09                     |

PC II: Single Parameter Colorimeter, DR 3900: VIS Spectrophotometer, DR 6000: UV-VIS Spectrophotometer Please note: Some methods require reagent blanks. For these, the number of tests varies.





<sup>\*:</sup> GHS hazard code will be available in the future

 $<sup>\</sup>boldsymbol{\cdot} \boldsymbol{\cdot}$  product is not subject to classification

# Which LCK cuvette tests for my photometer?

| Article<br>number | Parameter           | Measuring<br>range                     | Method                            | According to standard                 | Quality<br>control | Number<br>of tests | PC II | DR<br>3900 | DR<br>6000 | EU<br>hazard<br>code | GHS<br>hazard<br>code               |
|-------------------|---------------------|--|-----------------------------------|---------------------------------------|--------------------|--------------------|-------|------------|------------|----------------------|-------------------------------------|
| LCK514            | COD                 | 100 - 2000<br>mg/L 0 <sub>2</sub>      | Dichromate                        | ISO 6060-1989,<br>DIN 38409-H41-H44   | LCA708             | 25                 |       | •          | Ť          | T, C                 | GHS05,<br>GHS06,<br>GHS08,<br>GHS09 |
| LCK614            | COD                 | 50 - 300 mg/L<br>O <sub>2</sub>        | Dichromate                        | ISO 6060-1989,<br>DIN 38409-H41-H44   | LCA709             | 25                 |       | -          | •          | T, C                 | GHS05,<br>GHS06,<br>GHS08,<br>GHS09 |
| LCK714            | COD                 | 100 - 600<br>mg/L 0 <sub>2</sub>       | Dichromate                        | ISO 6060-1989,<br>DIN 38409-H41-H44   | 1218629            | 25                 |       | -          | •          | T, C                 | GHS05,<br>GHS06,<br>GHS08,<br>GHS09 |
| LCK914            | COD                 | 5 - 60 g/L 0 <sub>2</sub>              | Dichromate                        | ISO 6060-1989,<br>DIN 38409-H41-H44   |                    | 25                 |       | -          | •          | T, C                 | GHS05,<br>GHS06,<br>GHS08,<br>GHS09 |
| LCK214            | COD, mercury free   | 100 - 1000<br>mg/L 0 <sub>2</sub>      | Dichromate                        | ISO 6060-1989,<br>DIN 38409-H41       | 1218629            | 25                 |       | -          | •          | T, C                 | GHS05,<br>GHS08,<br>GHS09           |
| LCK329            | Copper              | 0.1 - 8.0 mg/L<br>Cu                   | Bathocuproine Disulphonic Acid    |                                       | LCA701             | 25                 |       | •          |            | -                    | -                                   |
| LCK229            | Copper              | 2 - 100 g/L Cu                         | Intrinsic Bath Colour             |                                       |                    | 25                 |       |            | -          | Xi                   | GHS05                               |
| LCK529            | Copper, trace       | 0.01 - 1.0 mg/L<br>Cu                  | Bathocuproine Disulphonic<br>Acid |                                       | LCA706             | 20                 |       | -          |            | -                    | -                                   |
| LCK315            | Cyanide             | 0.01 - 0.6 mg/L<br>CN                  | Babituric Acid-Pyridine           | ISO 6703-1-2-3-1984,<br>DIN 38405 D13 |                    | 25                 |       | •          | •          | C, Xn                | GHS05,<br>GHS07,<br>GHS08           |
| LCK319            | Cyanide             | 0.03 - 0.35<br>mg/L CN                 | HACH LANGE Method                 |                                       |                    | 24                 |       | -          | •          | Xi, N                | GHS05,<br>GHS07,<br>GHS09           |
| LCK323            | Fluoride            | 0.1 - 2.5 mg/L F                       | SPADNS                            |                                       | 29153              | 25                 |       | -          |            | -                    | GHS05                               |
| LCK325            | Formaldehyde        | 0.5 - 10 mg/L<br>H <sub>2</sub> C0     | Acetylacetone                     |                                       |                    | 24                 | •     | -          | •          | -                    | GHS07                               |
| LCS325            | Formaldehyde, trace | 0.01 - 1.0 mg/L<br>H <sub>2</sub> CO   | Acetylacetone                     |                                       |                    | 24                 |       | •          | •          | -                    | GHS07                               |
| LCK320            | Iron                | 0.2 - 6.0 mg/L<br>Fe                   | 1.10-Phenanthroline               | DIN 38405-D17                         | 2833649            | 24                 |       | •          | •          | Xi                   | GHS07                               |
| LCK321            | Iron                | 0.2 - 6.0 mg/L<br>Fe                   | 1.10-Phenanthroline               | ISO 6332-1988,<br>DIN 38406 E1-1      | LCA701             | 25                 |       | •          | •          | -                    | -                                   |
| LCK521            | Iron, trace         | 0.01 - 1.0<br>mg/L Fe                  | 1.10-Phenanthroline               | ISO 6332-1988,<br>DIN 38406 E1-1      | LCA706             | 20                 |       | •          |            | -                    |                                     |
| LCK306            | Lead                | 0.1 - 2.0 mg/L<br>Pb                   | PAR                               |                                       | LCA701             | 25                 |       | •          | •          | T+, N,<br>Xn         | GHS06,<br>GHS07,<br>GHS09           |
| LCK326            | Magnesium           | 0.5 - 50 mg/L<br>Mg                    | Metalphthalein                    |                                       | 1479442            | 25                 |       | -          | •          | -                    | -                                   |
| LYW185            | Menthol             | 0.5 - 15 mg/100<br>mL Menthol          | p-Dimethylaminobenzaldehyde       |                                       |                    | 25                 |       | •          | •          | С                    | GHS05                               |
| LCK330            | Molybdenum          | 3 - 300 mg/L Mo                        | Thioglycolic Acid                 |                                       |                    | 24                 |       | •          |            | Т                    | GHS05,<br>GHS06                     |
| LCK337            | Nickel              | 0.1 - 6.0 mg/L Ni                      | Dimethylglyoxime                  | DIN 38406-E11                         | LCA701             | 25                 |       | -          | •          | С                    | GHS05,<br>GHS07,<br>GHS08           |
| LCK237            | Nickel              | 5 - 120 g/L Ni                         | Intrinsic Baths Colour            |                                       |                    | 25                 |       | -          | -          | Xi                   | GHS05                               |
| LCK537            | Nickel, trace       | 0.05 - 1.0<br>mg/L Ni                  | Dimethylglyoxime                  |                                       | LCA706             | 20                 |       | =          | •          | C, 0                 | GHS05,<br>GHS07,<br>GHS08           |
| LCK339            | Nitrate             | 0.23 - 13.5 mg/L<br>NO <sub>3</sub> -N | 2.6-Dimethylphenol                | ISO 7890-1-2-1986,<br>DIN 38405 D9-2  | LCA703             | 25                 |       | •          | •          | С                    | GHS02,<br>GHS05,<br>GHS07           |



# Which LCK cuvette tests for my photometer?

| Article<br>number | Parameter                             | Measuring<br>range                       | Method  | According to standard                  | Quality<br>control | Number<br>of tests | PC II | DR<br>3900 | DR<br>6000 | EU<br>hazard<br>code | GHS<br>hazard<br>code               |
|-------------------|---------------------------------------|--|---|--|--------------------|--------------------|-------|------------|------------|----------------------|-------------------------------------|
| LCK340            | Nitrate                               | 5 - 35 mg/L<br>NO <sub>3</sub> -N        | 2.6-Dimethylphenol  | ISO 7890-1-2-1986, DIN 38405<br>D9-2   | LCA704             | 25                 |       | •          |            | С                    | GHS02,<br>GHS05                     |
| LCK341            | Nitrite                               | 0.015 - 0.6 mg/L<br>NO <sub>2</sub> -N   | Diazotisation   | EN ISO 26777, DIN 38405 D10            | LCA707             | 25                 |       | •          |            | Xi                   | GHS07                               |
| LCK342            | Nitrite                               | 0.6 - 6.0 mg/L<br>NO <sub>2</sub> -N     | Diazotisation   | EN ISO 26777, DIN 38405 D10            | LCA 709            | 25                 |       | •          | •          | Xi                   | GHS07                               |
| LCK541            | Nitrite, trace                        | 0.0015 - 0.03<br>mg/L NO <sub>2</sub> -N | Diazotisation   | EN ISO 26777, DIN 38405 D10            | 2340249            | 50                 |       | •          | •          | -                    | GHS07                               |
| LCK138            | Nitrogen total<br>(LATON)             | 1 - 16 mg/L TN <sub>b</sub>              | Koroleff Digestion<br>(Peroxodisulphate), and<br>Photometric Detection with<br>2.6-Dimethylphenol | EN ISO 11905-1                         | LCA709             | 25                 |       | •          | •          | T, C, Xi,<br>Xn      | GHS02,<br>GHS05,<br>GHS07,<br>GHS08 |
| LCK238            | Nitrogen total<br>(LATON)             | 5 - 40 mg/L TN <sub>b</sub>              | Koroleff Digestion<br>(Peroxodisulphate), and<br>Photometric Detection with<br>2.6-Dimethylphenol | EN ISO 11905-1                         | LCA700             | 25                 |       | •          | •          | T, C, Xi,<br>Xn      | GHS02,<br>GHS05,<br>GHS07,<br>GHS08 |
| LCK338            | Nitrogen total<br>(LATON)             | 20 - 100 mg/L<br>TN <sub>b</sub>         | Koroleff Digestion<br>(Peroxodisulphate), and<br>Photometric Detection with<br>2.6-Dimethylphenol | EN ISO 11905-1                         | LCA708             | 25                 |       | •          | •          | T, C, Xi,<br>Xn      | GHS02,<br>GHS05,<br>GHS07,<br>GHS08 |
| LCK365            | Organic acids                         | 50 - 2500 mg/L<br>as Acetic Acid         | Esterification  |  |                    | 25                 |       | •          | •          | С                    | GHS05,<br>GHS07,<br>GHS08,<br>GHS09 |
| LCK345            | Phenols                               | 0.05 - 5 mg/L<br>Phenols                 | 4-Nitroaniline  |  |                    | 24                 |       | •          | •          | Xn                   | GHS05,<br>GHS07,<br>GHS09           |
| LCK346            | Phenols                               | 5 - 200 mg/L<br>Phenols                  | 4-Aminoantipyrine   | ISO 6439-1990, DIN 38409 H16           |                    | 24                 |       | •          | •          | Xn, 0                | GHS03,<br>GHS07,<br>GHS08           |
| LCK049            | Phosphate, ortho                      | 1.6 - 30 mg/L<br>PO <sub>4</sub> -P      | Vanadate-Molybdate  |  | LCA703             | 25                 |       | •          | •          | С                    | GHS05                               |
| LCK348            | Phosphate,<br>ortho + total           | 0.5 - 5.0 mg/L<br>PO <sub>4</sub> -P     | Phosphormolybdenum Blue   | EN ISO 6878-1-1986, DIN 38405<br>D11-4 | LCA700,<br>LCA707  | 25                 |       |            | Ť.         | С                    | GHS05,<br>GHS07,<br>GHS08           |
| LCK349            | Phosphate,<br>ortho + total           | 0.05 - 1.5 mg/L<br>PO <sub>4</sub> -P    | Phosphormolybdenum Blue   | ISO 6878-1-1986, DIN 38405<br>D11-4    | LCA704,<br>LCA709  | 25                 | •     | •          | •          | С                    | GHS05,<br>GHS07,<br>GHS08           |
| LCK350            | Phosphate,<br>ortho + total           | 2 - 20 mg/L<br>PO <sub>4</sub> -P        | Phosphormolybdenum Blue   | ISO 6878-1-1986, DIN 38405<br>D11-4    | LCA703,<br>LCA708  | 25                 | •     | •          |            | C, Xn                | GHS05,<br>GHS07,<br>GHS08           |
| LCS349            | Phosphate,<br>ortho + total           | 0.01 - 0.5 mg/L<br>PO <sub>4</sub> -P    | Phosphormolybdenum Blue   | ISO 6878-1-1986, DIN 38405<br>D11-4    | LCA704,<br>LCA709  | 25                 |       | •          |            | С                    | GHS05,<br>GHS07,<br>GHS08           |
| LCK240            | Photometric<br>lodine sample<br>(PIS) | > 0.2 lodine<br>value                    | MEBAK Method  | Mebak II                               |                    | 25                 |       |            | •          | F                    | GHS02,<br>GHS05                     |
| LCK228            | Potassium                             | 5 - 50 mg/L K                            | Kalignost   |  | LCA700             | 25                 |       | -          | •          | -                    | GHS05,<br>GHS06,<br>GHS07,<br>GHS08 |
| LCK328            | Potassium                             | 8 - 50 mg/L K                            | Kalignost   |  | LCA700             | 24                 |       |            |            | Xn                   | GHS06                               |
| LCK354            | Silver                                | 0.04 - 0.8 mg/L<br>Ag                    | HACH LANGE Method   |  | 1461342            | 25                 |       | •          | •          | F                    | GHS02,<br>GHS07,<br>GHS08           |
| LCK355            | Silver                                | 5 - 400 mg/L<br>Ag (l)                   | HACH LANGE Method   |  | 1461342            | 24                 |       | •          | •          | С                    | GHS05                               |

PC II: Single Parameter Colorimeter, DR 3900: VIS Spectrophotometer, DR 6000: UV-VIS Spectrophotometer Please note: Some methods require reagent blanks. For these, the number of tests varies.





<sup>\*:</sup> GHS hazard code will be available in the future

 $<sup>\</sup>boldsymbol{\cdot} \boldsymbol{\cdot}$  product is not subject to classification

# Which LCK cuvette tests for my photometer?

| Article<br>number | Parameter                      | Measuring<br>range                | Method  | According to standard                  | Quality<br>control           | Number<br>of tests | PC II | DR<br>3900 | DR<br>6000 | EU<br>hazard<br>code | GHS<br>hazard<br>code               |
|-------------------|--------------------------------|-----------------------------------|---|--|------------------------------|--------------------|-------|------------|------------|----------------------|-------------------------------------|
| LCK318            | Sludge activity                | 5 - 200 μg<br>Formazan (SA)       | Colorimetric  | DIN 38409-H7                           |                              |                    |       |            |            | F                    | GHS02                               |
| LCK357            | Starch                         | 2 - 150 mg/L<br>Starch            | HACH LANGE Method   |  |                              | 25                 |       |            |            | -                    | -                                   |
| LCK153            | Sulphate                       | 40 - 150 mg/L<br>SO <sub>4</sub>  | Barium Sulphate   |  | LCA704                       | 25                 |       | •          | •          | Т                    | GHS06                               |
| LCK353            | Sulphate                       | 150 - 900 mg/L<br>SO <sub>4</sub> | Barium Sulphate   |  | LCA701,<br>LCA702,<br>LCA703 | 25                 |       | •          | •          | T                    | GHS06                               |
| LCK653            | Sulphide                       | 0.1 - 2.0 mg/L<br>S <sup>2-</sup> | Dimethyl-p-phenylenediamine   | ISO 10530-1991, DIN 38405-D26          |                              | 25                 |       | •          | •          | С                    | GHS05                               |
| LCK654            | Sulphite                       | 0.1 - 5.0 mg/L<br>SO <sub>3</sub> | HACH LANGE Method   |  |                              | 25                 |       |            |            | -                    | -                                   |
| LCK332            | Surfactants,<br>anionic        | 0.05 - 2.0 mg/L                   | Methylene Blue (MBA)  | ISO 7875-1-2-1984, DIN 38409-H<br>23-1 |                              | 25                 |       | •          | •          | Xn                   | GHS07,<br>GHS08                     |
| LCK331            | Surfactants, cationic          | 0.2 - 2.0 mg/L                    | Bromophenol Blue  |  |                              | 25                 |       | •          | •          | Xn, F                | GHS02,<br>GHS07,<br>GHS08           |
| LCK333            | Surfactants, nonionic          | 0.2 - 6.0 mg/L<br>as TRITON x 100 | TBPE  | DIN 38409-H23-2                        | LCA333                       | 25                 |       | •          | •          | Xn                   | GHS02,<br>GHS08                     |
| LCK334            | Surfactants, nonionic          | 0.1 - 20 g/L                      | CTAS  | DIN 38409-H23-2                        |                              | 25                 |       |            |            | Xn                   | GHS07,<br>GHS08                     |
| LCK433            | Surfactants, nonionic          | 6 - 200 mg/L<br>as TRITON x 100   | TBPE  | DIN 38409-H23-2                        |                              | 25                 |       | •          | •          | Xn                   | GHS02,<br>GHS08                     |
| LCK359            | Tin                            | 0.1 - 2.0 mg/L<br>Sn              | Pyridinfluoron (PYF)  |  |                              | 24                 |       | •          | •          | Т, О                 | GHS02,<br>GHS03,<br>GHS07,<br>GHS08 |
| LCK380            | TOC                            | 2 - 65 mg/L C                     | Difference Method (TOC is<br>determined as the difference<br>between the TC and TIC<br>values), Persulphate Digestion | DIN 38409-H3                           | 2833249                      | 25                 |       | •          | •          | Xn, Xi, O            | GHS03,<br>GHS07,<br>GHS08           |
| LCK381            | TOC                            | 60 - 735 mg/L C                   | Difference Method (TOC is<br>determined as the difference<br>between the TC and TIC<br>values), Persulphate Digestion | DIN 38409-H3                           | 2833149                      | 25                 |       | •          | •          | Xn, 0                | GHS03,<br>GHS07,<br>GHS08           |
| LCK385            | TOC                            | 3 - 30 mg/L C                     | Purging Method, Persulphate<br>Digestion  | EN 1484 ,DIN 38409-H3                  | LCA704                       | 25                 |       | •          | •          | Xn                   | GHS07,<br>GHS08                     |
| LCK386            | TOC                            | 30 - 300 mg/L C                   | Purging Method, Persulphate<br>Digestion  | EN 1484, DIN 38409-H3                  | LCA703                       | 25                 |       |            |            | Xn                   | GHS07,<br>GHS08                     |
| LCK387            | TOC                            | 300 - 3000<br>mg/L C              | Purging Method, Persulphate<br>Digestion  | EN 1484, DIN 38409-H3                  | LCA705                       | 20                 |       | •          | •          | Xn                   | GHS07,<br>GHS08                     |
| LCK242            | Vicinal<br>diketones<br>(VDK)  | 0.015 - 0.5<br>mg/kg Diacetyl     | Analogous MEBAK-Method  | MEBAK II                               |                              | 25                 |       |            | •          | T, N                 | GHS05,<br>GHS06,<br>GHS08,<br>GHS09 |
| LCK327            | Water<br>hardness              | 1 - 20 °dH                        | Metalphthalein  |  | 2833449                      | 25                 |       | •          | •          | -                    | -                                   |
| LCK427            | Water<br>hardness,<br>residual | 0.02 - 0.6 °dH                    | Metalphthalein  |  | 2833449                      | 24                 |       | •          | •          | Xi                   | -                                   |
| LCK360            | Zinc                           | 0.2 - 6.0 mg/L<br>Zn              | PAR   |  | LCA701                       | 24                 | -     | •          | •          | Xn                   | GHS07                               |
| LCS360            | Zinc, trace                    | 0.02 - 0.8<br>mg/L Zn             | PAR   |  | LCA701                       | 24                 |       | •          | •          | Xn                   | GHS07                               |
| LCK364            | Zirconium                      | 10 - 60 mg/L Zr                   | SurTec/HACH LANGE Method  |  |                              | 12 - 24            | •     | •          |            |                      | GHS05                               |



## Which Powder Pillows for my photometer?

### Quick reference guide

#### Powder Pillows - low-price methods with long shelf life



Powder Pillows are available for a large number of parameters and measuring ranges. Hermetically sealed in aluminium foil pillows, the PERMACHEM reagents have a shelf life of many years. The reagent is simply poured into the measuring cuvette together with the sample. The evaluation can be carried out visually, e.g. with a colour disk, or with a HACH LANGE photometer.

| Article<br>number | Parameter  | Measuring range  | Method                         | Quality<br>control  | Number<br>of tests | PC II | DR<br>900 | DR<br>3900 | DR<br>6000 | EU<br>hazard<br>code | GHS<br>hazard<br>code |
|-------------------|--|--|--------------------------------|---------------------|--------------------|-------|-----------|------------|------------|----------------------|-----------------------|
| 2242000           | Aluminium  | 0.008 - 0.800 mg/L Al  | Aluminon                       | 1417442             | 100                |       |           |            |            | Xi                   | *                     |
| 2603700           | Aluminium  | 0.002 - 0.250 mg/L Al  | Eriochrome Cyanine R           | 1417442             | 100                |       |           |            |            | F, Xn                | *                     |
| 2653299           | Ammonia  | 0.01 - 0.50 mg/L NH <sub>3</sub> -N                            | Salicylate                     | 189149              | 100                |       |           |            |            | Xn                   | *                     |
| 2668000           | Ammonia  | 0.01 - 0.50 mg/L NH <sub>3</sub> -N                            | Salicylate                     | 15349               | 100                |       |           |            |            | Xn                   | *                     |
| 2459200           | Ammonium compounds, quaternary                               | 0.2 - 5.0 mg/L as CTAB   | Direct Binary Complex          |                     | 100                |       |           | •          | •          | Xi                   | GHS07                 |
| 1206499           | Barium   | 2 - 100 mg/L Ba  | Turbidimetric                  | 1461142             | 100                |       |           |            |            | Xi                   | *                     |
| 2141299           | Benzotriazole, Tolyltriazole<br>Benzotriazole, Tolyltriazole | 1.0 - 20.0 mg/L Tolyltriazole<br>1.0 - 16.0 mg/L Benzotriazole | UV Photolysis                  |                     | 100                |       |           |            | •          | Xn                   | *                     |
| 1417099           | Boron  | 0.2 - 14.0 mg/L B  | Carmine                        |                     | 100                |       |           |            |            | -                    | *                     |
| 2802246           | Chloramine, mono   | 0.04 - 4.50 mg/L Cl <sub>2</sub>                               | Indophenol                     |                     | 50                 |       |           |            |            | C, Xn                | *                     |
| 2105569           | Chlorine, free   | 0.02 - 2.00 mg/L Cl <sub>2</sub>                               | DPD                            | 1426810,<br>2630020 | 100                | •     | •         |            |            | -                    | *                     |
| 1407099           | Chlorine, free   | 0.1 - 10.0 mg/L Cl <sub>2</sub>                                | DPD                            |                     | 100                |       |           |            |            | -                    | *                     |
| 2105528           | Chlorine, free,<br>Chlorine dioxide                          | 0.02 - 2.00 mg/L Cl <sub>2</sub>                               | DPD                            | 1426810,<br>2630020 | 1000               | -     |           |            |            | -                    | *                     |
| 2105628           | Chlorine, total  | 0.02 - 2.00 mg/L Cl <sub>2</sub>                               | DPD                            | 1426810,<br>2630020 | 1000               | -     | •         |            | •          | -                    | *                     |
| 2105669           | Chlorine, total,<br>Bromine, lodine                          | 0.02 - 2.00 mg/L Cl <sub>2</sub>                               | DPD                            | 1426810,<br>2630020 | 100                | •     |           |            |            | -                    | *                     |
| 2770900           | Chlorine dioxide   | 0.04 - 5.00 mg/L ClO <sub>2</sub>                              | DPD/Glycine                    |                     | 100                | •     |           |            | •          | -                    | *                     |
| 1271099           | Chromium   | 0.010 - 0.700 mg/L Cr (VI)                                     | 1,5-Diphenylcarbohydrazide     | 1425610             | 100                |       |           |            |            | Xi                   | *                     |
| 2242500           | Chromium, total  | 0.01 - 0.70 mg/L Cr  | Alkaline Hypobromite Oxidation | 1425610             | 100                |       | •         | •          | •          | T, C                 | *                     |
| 2651600           | Cobalt, Nickel   | 0.01 - 2.00 mg/L Co  | PAN                            | 2150342,<br>1417642 | 100                | •     |           |            |            | T, Xi                | *                     |
| 2105869           | Copper   | 0.04 - 5.00 mg/L Cu  | Bicinchoninate                 | 12842               | 100                |       |           |            |            | -                    | GHS07                 |
| 2603300           | Copper   | 2 - 210 μg/L Cu  | Porphyrin                      | 12842               | 100                |       |           |            |            | Xn                   | *                     |
| 2430200           | Cyanide  | 0.002 - 0.240 mg/L CN  | Pyridine-Pyrazalone            |                     | 100                |       |           |            |            | -                    | *                     |
| 246066            | Cyanuric acid  | 5 - 50 mg/L  | Turbidimetric                  |                     | 50                 |       |           |            |            | -                    | *                     |
| 2544800           | Iron   | 0.01 - 1.80 mg/L Fe  | FerroMo                        | 1417542             | 100                |       |           |            |            | -                    | *                     |
| 2105769           | Iron   | 0.02 - 3.00 mg/L Fe  | FerroVer                       | 1417542             | 100                |       |           |            |            | Xn                   | *                     |

PC II: Single Parameter Colorimeter, DR 900: Multi-Parameter Colorimeter, DR 3900: VIS Spectrophotometer, DR 6000: UV-VIS Spectrophotometer

<sup>\*:</sup> GHS hazard code will be available in the future

<sup>-:</sup> product is not subject to classification

# Which Powder Pillows for my photometer?

| Article<br>number | Parameter                     | Measuring range                       | Method                    | Quality<br>control | Number<br>of tests | PC II | DR<br>900 | DR<br>3900 | DR<br>6000 | EU<br>hazard<br>code | GHS<br>hazard<br>code |
|-------------------|-------------------------------|---------------------------------------|---------------------------|--------------------|--------------------|-------|-----------|------------|------------|----------------------|-----------------------|
| 2608799           | Iron                          | 0.012 - 1.800 mg/L Fe                 | TPTZ                      | 1417542            | 100                |       |           |            |            | Xn                   | *                     |
| 230166            | Iron                          | 0.009 - 1.400 mg/L Fe                 | FerroZine                 | 1417542            | 50                 |       |           |            |            | T                    | *                     |
| 103769            | Iron, ferrous                 | 0.02 - 3.00 mg/L Fe (II)              | 1,10 Phenanthroline       | 1417542            | 100                |       |           |            |            | Xn, N                | *                     |
| 2430000           | Manganese                     | 0.1 - 20.0 mg/L Mn                    | Periodate Oxidation       | 1279142            | 100                |       |           |            |            | Xi, O                | *                     |
| 2604100           | Molybdenum                    | 0.3 - 40.0 mg/L Mo                    | Mercaptoacetic Acid       | 1426510            | 100                |       |           |            |            | Xn                   | *                     |
| 2449400           | Molybdenum, Molybdate         | 0.02 - 3.00 mg/L Mo                   | Ternary Complex           | 1426510            | 100                |       |           |            |            | -                    | *                     |
| 2243500           | Nickel                        | 0.02 - 1.80 mg/L Ni                   | Heptoxime                 | 1417642            | 50                 |       |           |            |            | Xn                   | *                     |
| 2106169           | Nitrate                       | 0.3 - 30.0 mg/L NO <sub>3</sub> -N    | Cadmium Reduction         | 30749              | 100                |       |           |            |            | T, N                 | *                     |
| 2429800           | Nitrate                       | 0.01 - 0.50 mg/L NO <sub>3</sub> -N   | Cadmium Reduction         | 30749              | 100                |       |           |            |            | T, Xi, N             | *                     |
| 2107169           | Nitrite                       | 0.002 - 0.300 mg/L NO <sub>2</sub> -N | Diazotisation             | 2340249            | 100                |       |           |            |            | Xi                   | *                     |
| 2107569           | Nitrite                       | $2$ - $250~\mathrm{mg/L~NO}_2$        | Ferrous Sulphate          |                    | 100                |       |           |            |            | Xi                   | *                     |
| 2446600           | Oxygen scavengers             | 5 - 600 g/L Carbohydrazide            | Iron Reduction            |                    | 100                |       |           |            |            | С                    | *                     |
| 2243900           | Phenols                       | 0.002 - 0.200 mg/L Phenol             | 4-Aminoantipyrine         |                    | 100                |       |           |            |            | Xn                   | *                     |
| 2106069           | Phosphate, ortho              | 0.02 - 2.50 mg/L PO <sub>4</sub>      | Ascorbic Acid             | 256949             | 100                |       |           |            |            | Xi                   | *                     |
| 212528            | Phosphate, ortho              | 0.02 - 2.50 mg/L PO <sub>4</sub>      | Ascorbic Acid             | 256949             | 1000               |       |           |            |            | Xi                   | *                     |
| 2106028           | Phosphate, ortho              | 0.02 - 2.50 mg/L PO <sub>4</sub>      | Ascorbic Acid             | 256949             | 1000               |       |           |            |            | Xi                   | *                     |
| 2429700           | Phosphonates                  | 0.02 - 2.50 mg/L                      | Persulfate UV Oxidation   |                    | 100                |       |           |            |            | 0, Xi, N             | *                     |
| 2459100           | Potassium                     | 0.1 - 7.0 mg/L K                      | Tetraphenylborate         | 2240442            | 100                |       |           |            |            | F, T, Xi             | *                     |
| 2429600           | Silica                        | 1 - 100 mg/L SiO <sub>2</sub>         | Silicomolybdate           | 110649             | 100                |       |           |            |            | Xn                   | GHS07                 |
| 2459300           | Silica                        | 0.010 - 1.600 mg/L SiO <sub>2</sub>   | Heteropoly Blue           | 110649             | 100                |       |           |            |            | Xn                   | *                     |
| 2296600           | Silver                        | 0.02 - 0.70 mg/L Ag                   | Colorimetric              | 1461342            | 50                 |       |           | -          |            | T, Xi                | *                     |
| 2106769           | Sulphate                      | 2 - 70 mg/L SO <sub>4</sub>           | SulfaVer 4, turbidimetric | 257849             | 100                |       |           |            |            | T                    | *                     |
| 2495300           | Total Kjeldahl Nitrogen (TKN) | 1 - 150 mg/L TKN                      | Nessler                   |                    | 250                |       | -         | -          | -          | T+, C, N             | *                     |
| 2429300           | Zinc                          | 0.01 - 3.00 mg/L Zn                   | Zincon                    | 237842             | 100                |       |           |            |            | T, N                 | *                     |



## Which SWIFTESTs for my photometer?

## Quick reference guide

#### The right amout of DPD with the SWIFTEST



The SWIFTEST is a powder dispenser that releases the correct amount of DPD (N,N-diethyl-p-phenylenediamine) at the press of a button. It contains enough reagent for 250 chlorine tests (free or total chlorine). As a practical, attractively priced alternative, the SWIFTEST is ideal for laboratories with a high sample throughput, and for analysis in the field.

| Article<br>number | Product description  | Measuring range                  | Method | Quality<br>control  | Number<br>of tests | PC II | DR<br>900 | DR<br>3900 | DR<br>6000 | EU<br>hazard<br>code | GHS<br>hazard<br>code |
|-------------------|--|----------------------------------|--------|---------------------|--------------------|-------|-----------|------------|------------|----------------------|-----------------------|
| 2802400           | SWIFTEST DPD Total chlorine<br>reagent dispenser and<br>reagent vial | 0.02 - 2.00 mg/L Cl <sub>2</sub> | DPD    | 1426810,<br>2630020 | 250                | •     | •         | •          | •          | -                    | GHS07                 |
| 2802300           | SWIFTEST DPD Free chlorine<br>reagent dispenser and<br>reagent vial  | 0.02 - 2.00 mg/L Cl <sub>2</sub> | DPD    | 1426810,<br>2630020 | 250                | •     | Ť         | •          | •          | -                    | *                     |
| 2105660           | DPD Total chlorine, SWIFTEST dispenser reagent (refill)              | 0.02 - 2.00 mg/L Cl <sub>2</sub> | DPD    | 1426810,<br>2630020 | 250                |       |           | •          |            | -                    | GHS07                 |
| 2105560           | DPD Free Chlorine, SWIFTEST dispenser reagent (refill)               | 0.02 - 2.00 mg/L Cl <sub>2</sub> | DPD    | 1426810,<br>2630020 | 250                | •     | •         | -          | •          | ٠                    | -                     |

PC II: Single Parameter Colorimeter, DR 900: Multi-Parameter Colorimeter, DR 3900: VIS Spectrophotometer, DR 6000: UV-VIS Spectrophotometer

<sup>\*:</sup> GHS hazard code will be available in the future

<sup>-:</sup> product is not subject to classification

# Which ACCUVACs for my photometer?

## Quick reference guide

#### ACCUVAC - analysing without pipetting



The secret of the ACCUVAC is the vacuum in the sealed glass cuvette containing a measured amount of reagent. The test is carried out by immersing the tip of the ACCUVAC in the sample, then breaking it by applying moderate pressure. The vacuum draws the sample into the cuvette, whilst ensuring thorough mixing. The resulting colour is measured visually or photometrically.

| Article<br>number | Parameter                           | Measuring range                         | Method                     | Quality<br>control  | Number<br>of tests | PC II | DR<br>900 | DR<br>3900 | DR<br>6000 | EU<br>hazard<br>code | GHS<br>hazard<br>code |
|-------------------|-------------------------------------|---|----------------------------|---------------------|--------------------|-------|-----------|------------|------------|----------------------|-----------------------|
| 2502025           | Chlorine, free,<br>Chlorine dioxide | 0.02 - 2.00 mg/L Cl <sub>2</sub>        | DPD                        | 1426810,<br>2630020 | 25                 | •     |           | •          | •          | -                    | *                     |
| 2503025           | Chlorine, total,<br>Bromine, lodine | 0.05 - 4.50 mg/L Br <sub>2</sub>        | DPD                        | 2630020             | 25                 | •     | •         | •          | •          | -                    | *                     |
| 2503025           | Chlorine, total,<br>Bromine, lodine | 0.07 - 7.00 mg/L I <sub>2</sub>         | DPD                        | 2630020             | 25                 |       |           | •          |            | -                    | *                     |
| 2503025           | Chlorine, total,<br>Bromine, lodine | 0.02 - 2.00 mg/L Cl <sub>2</sub>        | DPD                        | 2630020             | 25                 |       |           | •          | •          | -                    | *                     |
| 2505025           | Chromium                            | 0.010 - 0.700 mg/L Cr (VI)              | 1,5-Diphenylcarbohydrazide | 1425610             | 25                 |       |           |            |            | Xi                   | *                     |
| 2504025           | Copper                              | 0.04 - 5.00 mg/L Cu                     | Bicinchoninate             | 2833649             | 25                 |       |           |            |            | Xn                   | *                     |
| 2506025           | Fluoride                            | 0.02 - 2.00 mg/L F                      | SPADNS                     | 29153               | 25                 |       |           |            |            | С                    | *                     |
| 2507025           | Iron                                | 0.02 - 3.00 mg/L Fe                     | FerroVer                   | 1417542             | 25                 |       |           |            |            | Xn                   | *                     |
| 2510025           | Iron                                | 0.012 - 1.800 mg/L Fe                   | TPTZ                       | 1417542             | 25                 |       |           |            |            | Xn                   | *                     |
| 2514025           | Iron                                | 0.02 - 3.00 mg/L Fe (II)                | 1,10 Phenanthroline        | 2833649             | 25                 |       |           | •          |            | Xn, N                | *                     |
| 2511025           | Nitrate                             | 0.3 - 30.0 mg/L mg/L NO <sub>3</sub> -N | Cadmium Reduction          | 30749               | 25                 |       |           |            |            | T+, N                | *                     |
| 2512025           | Nitrite                             | 0.002 - 0.300 mg/L NO <sub>2</sub> -N   | Diazotisation              | 2340249             | 25                 |       |           |            |            | Xi                   | *                     |
| 2501025           | Oxygen, dissolved                   | 6 - 800 μg/L 0 <sub>2</sub>             | Indigo Carmine             |                     | 25                 |       |           | -          |            | -                    | *                     |
| 2515025           | Oxygen, dissolved                   | 0.3 - 15.0 mg/L O <sub>2</sub>          | HRDO                       |                     | 25                 |       |           |            |            | Xn, N                | *                     |
| 2516025           | Ozone                               | 0.01 - 0.25 mg/L 0 <sub>3</sub>         | Indigo                     |                     | 25                 |       |           |            |            | Xn                   | *                     |
| 2517025           | Ozone                               | 0.01 - 0.75 mg/L 0 <sub>3</sub>         | Indigo                     |                     | 25                 |       |           |            | -          | Xn                   | *                     |
| 2518025           | Ozone                               | 0.01 - 1.50 mg/L O <sub>3</sub>         | Indigo                     |                     | 25                 |       |           |            |            | Xn                   | *                     |
| 2508025           | Phosphate                           | 0.02 - 2.50 mg/L PO <sub>4</sub>        | Ascorbic Acid              | 256949              | 25                 |       |           |            |            | Xi                   | *                     |
| 2525025           | Phosphate, ortho                    | 0.3 - 45.0 mg/L PO <sub>4</sub>         | Molybdovanadate            | 256949              | 25                 |       |           | -          |            | С                    | *                     |
| 2509025           | Sulphate                            | 2 - 70 mg/L SO <sub>4</sub>             | SulfaVer 4                 | 257849              | 25                 |       |           | -          |            | T                    | *                     |

PC II: Single Parameter Colorimeter, DR 900: Multi-Parameter Colorimeter, DR 3900: VIS Spectrophotometer, DR 6000: UV-VIS Spectrophotometer



<sup>\*:</sup> GHS hazard code will be available in the future

<sup>-:</sup> product is not subject to classification

## Which liquid reagent tests for my photometer?

## Quick reference guide

#### Reagent solutions, economic liquid reagent tests and rapid liquid systems



Reagent tests for the determination of numerous parameters required in drinking, waste and process water applications as well as product control and monitoring. A cost-effective solution for your high-volume testing and serial analysis.

| Article<br>number | Parameter           | Measuring range   | Method   | Quality<br>control  | Number<br>of tests | PC II | DR<br>900 | DR<br>3900 | DR<br>6000 | EU<br>hazard<br>code | GHS<br>hazard<br>code |
|-------------------|---------------------|---|--|---------------------|--------------------|-------|-----------|------------|------------|----------------------|-----------------------|
| 2458200           | Ammonia             | 0.02 - 2.50 mg/L NH <sub>3</sub> -N                         | Nessler  |                     | 250                |       |           |            | -          | T+, N                | *                     |
| 2242200           | Cadmium             | 0.7 - 80 g/L Cd   | Dithizone  | 1402442             | 60 - 100           |       |           |            |            | T+, C, N             | *                     |
| 2556900           | Chlorine            | 0.02 - 2.00 mg/L $\mathrm{Cl}_2$                            | DPD  | 1426810,<br>2630020 | 450                |       |           | •          | •          | Xn, N                | *                     |
| 2557000           | Chlorine            | 0.02 - 2.00 mg/L $\mathrm{Cl}_2$                            | DPD  | 2630020,<br>1426810 | 450                |       |           | •          | •          | C, Xn, N             | *                     |
| HPT310            | Chlorine            | 0.02 - 2.00 mg/L $\mathrm{Cl}_2$                            | DPD  | 2630020,<br>1426810 | 100                | •     | •         | •          | •          | Xi                   | *                     |
| HPT210            | Chlorine, free      | 0.02 - 2.00 mg/L $\mathrm{Cl}_2$                            | DPD  | 2630020,<br>1426810 | 100                |       | •         | •          | •          | Xi                   | *                     |
| LCW510            | Chlorine/Ozone      | $0.1 - 1.5 \text{ mg/L Cl}_2 / \text{ O}_3$ (round cuvette) | DPD  |                     |                    |       |           | •          | •          | -                    | GHS07                 |
| 2242300           | Chlorine dioxide    | 0.01 - 1.00 mg/L ClO <sub>2</sub>                           | Chlorophenol Red   |                     | 100                |       |           |            |            | Xi                   | *                     |
| HPT240            | Chlorine dioxide    | $0.02$ - $0.50~\mathrm{mg/L}~\mathrm{ClO}_{_2}$             | Amaranth Method  |                     | 100                |       |           |            |            | -                    | *                     |
| 2651600           | Cobalt, Nickel      | 0.01 - 2.00 mg/L Co   | PAN  | 2150342,<br>1417642 | 100                | •     |           | •          | •          | T, Xi                | *                     |
| 44449             | Fluoride            | 0.02 - 2.00 mg/L F  | SPADNS   | 29153               | 125                |       |           |            |            | С                    | *                     |
| 2257700           | Formaldehyde        | 3 - 500 μg/L CH <sub>2</sub> O                              | MBTH   |                     | 100                |       |           |            | •          | Xn                   | *                     |
| 2603100           | Hardness            | 8 - 1000 μg/L CaCO <sub>3</sub>                             | Chlorophosphonazo  | 2833449             | 100                |       |           |            |            | Xn,                  | *                     |
| 2319900           | Hardness, Ca and Mg | $0.05$ - $4.00$ mg/L Ca as $\mathrm{CaCO_3}$                | Calmagite Colorimetric   | 218710              | 100                |       | •         | •          | •          | С                    | *                     |
| 179032            | Hydrazine           | 4 - 600 μg/L N <sub>2</sub> H <sub>4</sub>                  | p-Dimethylaminobenzaldehyde  |                     | 100                |       |           |            |            | С                    | *                     |
| LCW025            | Hydrazine           | 0.01 - 2.0 mg/L N <sub>2</sub> H <sub>4</sub>               | 4-Dimethylaminobenzaldehyde  |                     | 60                 |       |           |            |            | -                    | GHS05                 |
| LCW058            | Hydrogen peroxide   | 1 - 10 g/L H <sub>2</sub> O <sub>2</sub>                    | Peroxomolybdate  |                     | 40                 |       |           |            |            | -                    | GHS05                 |
| 230149            | Iron                | 0.009 - 1.400 mg/L Fe                                       | FerroZine  | 1417542             | 500 -<br>1000      |       |           | •          | •          | T                    | *                     |
| LCW021            | Iron                | 0.005 - 0.25 mg/L Fe  | Iron(II) ions react with<br>FerroZine to form a violet<br>complex compound |                     | 50                 |       |           | •          | •          | С                    | GHS05                 |
| 2375000           | Lead                | 5 - 150 μg/L Pb   | LeadTrak   | 1426210             | 20                 |       |           |            |            | С                    | *                     |

PC II: Single Parameter Colorimeter, DR 900: Multi-Parameter Colorimeter, DR 3900: VIS Spectrophotometer, DR 6000: UV-VIS Spectrophotometer Please note: Some methods require reagent blanks. For these, the number of tests varies.

 $<sup>\</sup>ensuremath{^{\star}}\xspace$  GHS hazard code will be available in the future

<sup>-:</sup> product is not subject to classification

# Which liquid reagent tests for my photometer?

| Article<br>number | Parameter       | Measuring range  | Method                               | Quality<br>control | Number<br>of tests | PC II | DR<br>900 | DR<br>3900 | DR<br>6000 | EU<br>hazard<br>code | GHS<br>hazard<br>code                         |
|-------------------|-----------------|--|--------------------------------------|--------------------|--------------------|-------|-----------|------------|------------|----------------------|---|
| 2651700           | Manganese       | 0.006 - 0.700 mg/L Mn  | PAN                                  | 1279142            | 50                 |       |           |            |            | T, N                 | *   |
| LCW532            | Manganese       | 0.005 - 0.5 mg/L Mn  | 1-(2-pyridylazo)-2-naphthol<br>(PAN) |                    | 50                 |       |           | •          | •          | T, C,<br>N, Xi       | GHS02,<br>GHS05,<br>GHS06,<br>GHS08,<br>GHS09 |
| LCW032            | Manganese       | 0.2 - 5 mg/L Mn<br>(round cuvette or 10 mm<br>rectangular cuvette) | Formaldoxime                         | LCA706             | 50                 |       |           |            |            | Т                    | GHS05,<br>GHS06,<br>GHS07,<br>GHS08,<br>GHS09 |
| 2658300           | Mercury         | 0.1 - 2.5 μg/L Hg  | Cold Vapour Concentration            | 1419542            | 25                 |       |           |            |            | 0, T, C              | *   |
| 2076032           | ortho Phosphate | 0.3 - 45.0 mg/L PO <sub>4</sub>                                    | Molybdovanadate                      | 2109210            | 50                 |       |           |            |            | С                    | *   |
| 2657512           | рН              | 6.5 - 8.5 units  | Colorimetric Phenol Red              |                    | 50                 |       |           |            |            | -                    | *   |
| 2076049           | Phosphate       | 0.3 - 45.0 mg/L PO <sub>4</sub>                                    | Molybdovanadate                      | 2109210            | 250                |       |           |            |            | С                    | *   |
| 2244100           | Phosphate       | 0.23 - 30.00 mg/L PO <sub>4</sub>                                  | Amino Acid                           | 2109210            | 100                |       |           |            |            | T                    | *   |
| LCW250            | Reducing agent  | 0.02 - 1.0 mg/L DEHA   | Iron Reduction Method                |                    | 100                |       |           | -          |            | С                    |   |
| 2553500           | Silica          | 3 - 1000 μg/L SiO <sub>2</sub>                                     | Heteropoly Blue                      | 110649             | 100                |       |           | •          | •          | Xi                   | *   |
| 2581400           | Silica          | 3 - 1000 μg/L SiO <sub>2</sub>                                     | Heteropoly Blue                      | 110649             | 40                 |       |           |            | •          | Xi                   | *   |
| 2678500           | Silica          | 3 - 1000 μg/L SiO <sub>2</sub>                                     | Heteropoly Blue                      | 110649             | 250                |       |           | -          | •          | Xn                   | *   |
| LCW028            | Silica          | 0.01 - 0.8 mg/L SiO <sub>2</sub>                                   | Molybdenum Blue                      |                    | 50                 |       |           | •          |            | Xi, Xn               |   |
| 2244500           | Sulphide        | 5 - 800 μg/L S <sup>2-</sup>                                       | Methylene Blue                       |                    | 100                |       |           | •          |            | C, T                 | *   |
| LCW053            | Sulphide        | 0.1 - 2.0 mg/L S <sup>2-</sup>                                     | Dimethyl-p-phenylenediamine          |                    | 25 - 49            |       |           |            |            | С                    |   |
| HPT430            | Sulphite        | 0.1 - 5.0 mg/L SO <sub>3</sub>                                     | HACH LANGE Method                    | 2267410            | 100                |       |           |            |            | -                    | *   |
| LCW054            | Sulphite        | 0.1 - 5.0 mg/L SO <sub>3</sub>                                     | HACH LANGE Method                    | 2267410            | 100                |       |           |            |            | -                    | GHS07   |
| 2244600           | Tannin & Lignin | 0.1 - 9.0 mg/L as Tannic Acid                                      | Tyrosine                             |                    | 100                |       |           |            |            | -                    | *   |
| 2790800           | Trihalomethanes | 10 - 600 μg/L CHCl <sub>3</sub>                                    | THM Plus                             |                    | 50 - 99            |       |           |            |            | T, C                 | *   |
| 2244700           | Volatile acids  | 27 - 2800 mg/L HOAc  | Esterification                       |                    | 100                |       |           |            |            | C, Xn, N             | *   |



# Which TEST'N'TUBEs for my photometer?

## Quick reference guide

#### TEST'N'TUBEs - Safe and convenient testing



TEST'N'TUBE cuvette tests are completely equipped with all premeasured reagents, optimised for reliable measuring results and easy handling. Capped 16 mm vials provide a self-contained package for mixing and measurement. All necessary reagents and vials are contained in the package.

| Article<br>number | Parameter | Measuring range                       | Method                     | Quality<br>control              | Number<br>of tests | DR<br>900 | EU<br>hazard<br>code | GHS<br>hazard<br>code               |
|-------------------|-----------|---------------------------------------|----------------------------|---------------------------------|--------------------|-----------|----------------------|-------------------------------------|
| 2604545           | Ammonia   | 0.02 - 2.50 mg/L NH <sub>3</sub> -N   | Salicylate                 | 189149,<br>15349                | 25 - 50            |           | C, Xn                | *                                   |
| 2606945           | Ammonia   | 0.4 - 50.0 mg/L NH <sub>3</sub> -N    | Salicylate                 | 189149,<br>15349                | 25 - 50            |           | C, Xn                | *                                   |
| 2105545           | Chlorine  | 0.09 - 5.00 mg/L $\mathrm{Cl}_2$      | DPD                        | 1426810,<br>2630020             | 50                 |           | -                    | *                                   |
| 2125851           | COD       | 3 - 150 mg/L O <sub>2</sub>           | Dichromate                 | 1218629,<br>1218649,<br>2253929 | 25                 | •         | С                    | GHS05,<br>GHS06,<br>GHS08,<br>GHS09 |
| 2125951           | COD       | 20 - 1500 mg/L $\mathrm{O}_2$         | Dichromate                 | 1218629,<br>1218649,<br>2253929 | 25                 | •         | T, C                 | GHS05,<br>GHS06,<br>GHS08,<br>GHS09 |
| 2345852           | COD       | 25 - 150 mg/L 0 <sub>2</sub>          | Dichromate without mercury | 1218629,<br>1218649,<br>2253929 | 25                 | •         | С                    | GHS05,<br>GHS09                     |
| 2345952           | COD       | 0 - 1500 mg/L 0 <sub>2</sub>          | Dichromate without mercury | 1218629,<br>1218649,<br>2253929 | 25                 | •         | T, C                 | GHS05,<br>GHS06,<br>GHS08,<br>GHS09 |
| 2415851           | COD       | 0.7 - 40 mg/L 0 <sub>2</sub>          | Dichromate                 | 1218629,<br>1218649,<br>2253929 | 25                 | •         | С                    | GHS05,<br>GHS08,<br>GHS09           |
| 2415951           | COD       | 200 - 15000 mg/L 0 <sub>2</sub>       | Dichromate                 | 1218629,<br>1218649,<br>2253929 | 25                 | •         | T, C                 | GHS05,<br>GHS06,<br>GHS08,<br>GHS09 |
| 2623451           | COD       | 20 - 1000 mg/L O <sub>2</sub>         | Manganese (III)            | 1218629,<br>1218649,<br>2253929 | 25                 | •         | С                    | *                                   |
| 2605345           | Nitrate   | 0.2 - 30.0 mg/L NO <sub>3</sub> -N    | Chromotropic Acid          | 30749                           | 50                 | •         | С                    | GHS05,<br>GHS07                     |
| 2608345           | Nitrite   | 0.003 - 0.500 mg/L NO <sub>2</sub> -N | Diazotization              | 2340249                         | 50                 |           | Xi                   | *                                   |

DR 900: Multi-Parameter Colorimeter


Please note: Some methods require reagent blanks. For these, the number of tests varies.

<sup>\*:</sup> GHS hazard code will be available in the future

<sup>-:</sup> product is not subject to classification

# Which TEST'N'TUBEs for my photometer?

| Article<br>number | Parameter        | Measuring range                  | Method   | Quality<br>control           | Number<br>of tests | DR<br>900 | EU<br>hazard<br>code | GHS<br>hazard<br>code               |
|-------------------|------------------|----------------------------------|--|------------------------------|--------------------|-----------|----------------------|-------------------------------------|
| 2672245           | Nitrogen, total  | 0.5 - 25.0 mg/L N                | Persulphate Digestion                          | 189149,<br>15349,<br>2406549 | 25 - 50            | •         | C, 0                 | GHS03,<br>GHS05,<br>GHS07,<br>GHS08 |
| 2714100           | Nitrogen, total  | 10 - 150 mg/L N                  | Persulphate Digestion                          | 15349,<br>2406549            | 25 - 50            | •         | C, 0                 | GHS03,<br>GHS05,<br>GHS07,<br>GHS08 |
| 2742545           | Phosphate        | 0.06 - 5.00 mg/L PO <sub>4</sub> | Ascorbic Acid                                  | 2109210                      | 25 - 50            |           | Xi                   | GHS07                               |
| 2742745           | Phosphate        | 0.06 - 5.00 mg/L PO <sub>4</sub> | Ascorbic Acid                                  | 2109210                      | 25 - 50            | •         | C, 0                 | GHS03,<br>GHS05,<br>GHS07,<br>GHS08 |
| 2767345           | Phosphate, ortho | 1.0 - 100.0 mg/L PO <sub>4</sub> | Molybdovanadate                                | 256949                       | 25 - 50            |           | Xi                   | *                                   |
| 2742645           | Phosphate, total | 0.06 - 3.50 mg/L PO <sub>4</sub> | PhosVer 3 with Acid Persulfate<br>Digestion    | 2109210                      | 25 - 50            | •         | C, 0                 | GHS03,<br>GHS05,<br>GHS07,<br>GHS08 |
| 2767245           | Phosphate, total | 1.0 - 100 mg/L PO <sub>4</sub>   | Molybdovanadate with Acid Persulfate Digestion | 256949                       | 25 - 50            | •         | C, 0                 | GHS03,<br>GHS05,<br>GHS07,<br>GHS08 |
| 2760345           | тос              | 0.3 - 20.0 mg/L C                | Direct   |                              | 25 - 50            | •         | 0, C, Xn             | GHS03,<br>GHS05,<br>GHS07,<br>GHS08 |
| 2760445           | тос              | 100 - 700 mg/L C                 | Direct   |                              | 25 - 50            | •         | O, C, Xn             | GHS03,<br>GHS05,<br>GHS07,<br>GHS08 |
| 2815945           | TOC              | 15 - 150 mg/L C                  | Direct   |                              | 25 - 50            | •         | 0, C, Xn             | GHS03,<br>GHS05,<br>GHS07,<br>GHS08 |





# Standard Solutions - Multi-parameter for Analytical Quality Assurance



The comprehensive ADDISTA AQA system for HACH LANGE cuvette tests contains a standard solution plus two round-robin solutions which allow the user to participate in analysis checking free of charge.

Lot number, expiry date and target values by parameter are delivered via RFID tag on the packaging.

| Article number | For the following cuvette tests / parameters  |
|----------------|---|
| LCA700         | $ \begin{array}{l} - \text{LCK304 Ammonium, 0.015-2.0 mg/L NH}_4\text{-N} \\ - \text{LCK311 Chloride, 1-70 mg/L Cl} \\ - \text{LCK228 Potassium, 5-50 mg/L K} \\ - \text{LCK328 Potassium, 8-50 mg/L K} \\ - \text{LCK348 Phosphate (ortho), 0.5-5.0 mg/L PO}_4\text{-P} \\ - \text{LCK414 COD, 5-60 mg/L O}_2 \\ - \text{LCK238 LATON, 5-40 mg/L TN}_6 \end{array} $   |
| LCA701         | - LCK306 Lead, 0.1-2.0 mg/L Pb<br>- LCK321 Iron, 0.2-6.0 mg/L Fe<br>- LCK329 Copper, 0.1-8.0 mg/L Cu<br>- LCK337 Nickel, 0.1-6.0 mg/L Ni<br>- LCK353 Sulphate, 150-900 mg/L SO <sub>4</sub><br>- LCK360 Zinc, 0.2-6.0 mg/L Zn   |
| LCA702         | - LCK301 Aluminium, 0.02-0.5 mg/L Al - LCK308 Cadmium, 0.02-0.3 mg/L Cd - LCK313 Chromium (VI), 0.03-1.0 mg/L Cr - LCK313 Chromium (total), 0.03-1.0 mg/L Cr - LCS313 Chromium trace, 0.005-0.25 mg/L Cr - LCK353 Sulphate, 150-900 mg/L SO <sub>4</sub>  |
| LCA703         | $ \begin{array}{l} - \text{LCK049 Orthophosphate, 1.6-30 mg/L PO}_4\text{-P} \\ - \text{LCK114 COD, 150-1000 mg/L O}_2 \\ - \text{LCI400 COD, 0-1000 mg/L O}_2 \\ - \text{LCK303 Ammonium, 2-47 mg/L NH}_4\text{-N} \\ - \text{LCK311 Chloride, 1-70 mg/L Cl} \\ - \text{LCK339 Nitrate, 0.23-13.5 mg/L NO}_3\text{-N} \\ - \text{LCK350 Phosphate (ortho), 2-20 mg/L PO}_4\text{-P} \\ - \text{LCK353 Sulphate, 150-900 mg/L SO}_4 \\ - \text{LCK386 TOC, 30-300 mg/L C} \end{array} $ |
| LCA704         | - LCK153 Sulphate, 40-150 mg/L $\rm SO_4$<br>- LCK305 Ammonium, 1-12 mg/L $\rm NH_4$ -N<br>- LCK311 Chloride, 1-70 mg/L $\rm Cl$<br>- LCK314 COD, 15-150 mg/L $\rm O_2$<br>- LCK340 Nitrate, 5-35 mg/L $\rm NO_3$ -N<br>- LCK349 Phosphate (ortho), 0.05-1.5 mg/L $\rm PO_4$ -P<br>- LCK385 TOC, 3-30 mg/L $\rm C$  |

| Article number | For the following cuvette tests / parameters  |
|----------------|---|
| LCA705         | - LCK014 COD, 1000-10000 mg/L O <sub>2</sub>  |
|                | - LCK302 Ammonium, 47-130 mg/L NH <sub>4</sub> -N<br>- LCK311 Chloride, 1-70 mg/L Cl                    |
|                | - LCK387 TOC, 300-3000 mg/L C   |
| LCA706         | - LCK521 Iron trace, 0.01-1.0 mg/L Fe   |
|                | - LCK529 Copper trace, 0.01-1.0 mg/L Cu   |
|                | - LCK537 Nickel trace, 0.05-1.0 mg/L Ni   |
|                | - LCW032 Manganese, 0.02-5.0 mg/L Mn  |
| LCA707         | - LCK341 Nitrite, 0.015-0.6 mg/L NO <sub>2</sub> -N   |
|                | - LCK614 COD, 50-300 mg/L O <sub>2</sub><br>- LCK348 Phosphate (total), 0.5-5.0 mg/L PO <sub>4</sub> -P |
| 1.04700        | 7   |
| LCA708         | - LCK338 LATON, 20-100 mg/L TN <sub>b</sub><br>- LCK514 COD, 100-2000 mg/L O <sub>2</sub>               |
|                | - LCK350 Phosphate (total), 2-20 mg/L PO <sub>4</sub> -P  |
| LCA709         | - LCK138 LATON, 1-16 mg/L TN  |
|                | - LCK614 COD, 50-300 mg/L 0   |
|                | - LCK349 Phosphate (total), 0.05-1.5 mg/L PO <sub>4</sub> -P  |
|                | - LCK342 Nitrite, 0.6-6.0 mg/L NO <sub>2</sub> -N   |
| 2833149        | - Ammonia 15 mg/L NH <sub>3</sub> -N  |
|                | - Nitrate 10 mg/L NO <sub>3</sub> -N  |
|                | - COD 500 mg/L O <sub>2</sub>   |
|                | - Phosphate 100 mg/L PO <sub>4</sub>  |
|                | - Sulphate 400 mg/L SO <sub>4</sub><br>- TOC 161 mg/L C   |
| 2833249        | - Ammonia 2.0 mg/L NH <sub>3</sub> -N / 2.1 mg/L NH <sub>4</sub> -N                                     |
| 2033243        | - Nitrate 4.0 mg/L NO <sub>3</sub> -N   |
|                | - Phosphate 2.0 mg/L PO <sub>4</sub>  |
|                | - COD 25 mg/L O <sub>2</sub>  |
|                | - Sulphate 50 mg/L SO <sub>4</sub>  |
|                | - TOC 8 mg/L C  |
|                |   |






# Standard Solutions - Single parameter for **Analytical Quality Assurance**



Regular use of standard solutions can ensure laboratory process control, increase your confidence, and help provide evidence of performance to inspectors, regulators, and clients. Single parameters are available in a variety of analytes and concentrations for proof of accuracy.

| Parameter              | Article<br>number | Product description  | Concentration                | EU<br>hazard<br>code | GHS<br>hazard<br>code |
|------------------------|-------------------|--|------------------------------|----------------------|-----------------------|
| Alkalinity             | 2349732           | Sulphuric acid standard solution, 0.035 N, 100 mL MDB  | 0.035 N                      | -                    | *                     |
| Alkalinity             | 20353             | Sulphuric acid standard solution, 0.020 N, 1 L   | 0.020 N                      | -                    | *                     |
| Ammonia                | 15349             | Ammonia standard solution, 10 mg/L NH <sub>2</sub> -N, 500 mL                                  | 10 mg/L NH <sub>2</sub> -N   | -                    | *                     |
| Ammonia                | 189149            | Ammonia standard solution, 1mg/L NH <sub>2</sub> -N, 500 mL                                    | 1 mg/L NH <sub>3</sub> -N    | -                    | *                     |
| Ammonia                | 2406549           | Ammonia standard solution, 100 mg/L NH <sub>3</sub> -N, 500 mL                                 | 100 mg/L NH <sub>3</sub> -N  | -                    | *                     |
| AOX                    | LCA390            | ADDISTA Mono standard for AOX cuvette test LCK390  | Lot specific concentration   | -                    | *                     |
| BOD                    | LCA555            | ADDISTA Mono standard for BOD cuvette test LCK555  | 200 mg/L O <sub>2</sub>      | Xn, 0                | GHS03,<br>GHS07       |
| BOD                    | 1486510           | BOD standard solution, 300 mg/L O <sub>2</sub> (NIST), 10 mL, 16 pcs.                          | 300 mg/L 0 <sub>2</sub>      | -                    | *                     |
| BOD                    | 1486610           | BOD standard solution, 3000 mg/L 0 <sub>2</sub> (NIST), 10 mL, 16 pcs.                         | 3000 mg/L O <sub>2</sub>     | -                    | *                     |
| Chlorine               | LCA310            | ADDISTA Mono standard for chlorine cuvette test LCK310   | 25 - 30 mg/L Cl <sub>2</sub> | -                    | *                     |
| Chlorine               | 1426810           | Chlorine standard solution, 50-75 mg/L Cl <sub>2</sub> (NIST)                                  | 50 - 75 mg/L Cl <sub>2</sub> | -                    | *                     |
| Chlorine               | 2630020           | Chlorine standard solution, 25-30 mg/L Cl <sub>2</sub> (NIST), 20 pcs.                         | 25 - 30 mg/L Cl <sub>2</sub> | -                    | *                     |
| Chlorine               | 2635300           | SpecCheck Gel secondary standard kit-LR chlorine, DPD  | 0 - 2.0 mg/L Cl <sub>2</sub> | -                    | *                     |
| COD                    | 1218629           | COD standard solution, 300 mg/L 0 <sub>2</sub> (NIST), 200 mL                                  | 300 mg/L 0 <sub>2</sub>      | -                    | *                     |
| COD                    | 2253929           | COD standard solution, 1000 mg/L 0 <sub>2</sub> (NIST), 200 mL                                 | 1000 mg/L O <sub>2</sub>     | -                    | *                     |
| COD                    | 1218649           | COD standard solution, 300 mg/L O <sub>2</sub> (NIST), 500 mL                                  | 300 mg/L 0 <sub>2</sub>      | -                    | *                     |
| Colour                 | 141453            | Colour standard solution, 500 Pt Co Units, 1 L   | 500 Pt Co units              | T                    | *                     |
| Colour                 | 2602853           | Colour standard solution, 15 Pt Co Units, 1L   | 15 Pt Co units               | T                    | *                     |
| Conductivity           | 1440042           | Sodium chloride standard solution, 1000 µS/cm (NIST), 100 mL                                   | 1000 μS/cm                   | -                    | *                     |
| Conductivity           | 1440049           | Sodium chloride standard solution, 1000 µS/cm (NIST), 500 mL                                   | 1000 μS/cm                   | -                    | *                     |
| Conductivity           | 210553            | Sodium chloride standard solution, 1990 µS/cm (NIST), 1 L                                      | 1990 μS/cm                   | -                    | *                     |
| Conductivity           | 2971849           | Sodium chloride standard solution, 100 µS/cm (NIST), 500 mL                                    | 100 μS/cm                    | -                    | *                     |
| Conductivity           | 2972249           | Sodium chloride standard solution, 10000 µS/cm (NIST), 500 mL                                  | 10000 μS/cm                  | -                    | *                     |
| Iron                   | 1417542           | Iron standard solution 100.0 mg/L Fe (NIST), 100 mL  | 100 mg/L Fe                  | -                    | *                     |
| Nitrite                | 2340249           | Nitrite standard solution, 250 µg/mL NO <sub>2</sub> -N, APHA, 500 mL                          | 250 mg/L NO <sub>2</sub> -N  | -                    | *                     |
| Phosphate              | 1424342           | Phosphate standard solution, 15 mg/L PO <sub>4</sub> , 100 mL                                  | 15 mg/L PO <sub>4</sub>      | -                    | *                     |
| Phosphate              | 17149             | Phosphate standard solution, 50 mg/L PO <sub>4</sub> (NIST), 500 mL                            | 50 mg/L PO <sub>4</sub>      | -                    | *                     |
| Phosphate              | 256949            | Phosphate standard solution, 1 mg/L PO <sub>4</sub> , 500 mL                                   | 1 mg/L PO <sub>4</sub>       | -                    | *                     |
| Silica                 | 110649            | Silica standard solution, 1 mg/L SiO <sub>2</sub> (NIST), 500 mL                               | 1 mg/L SiO <sub>2</sub>      | -                    | *                     |
| Sulphate               | 2175749           | Sulphate standard solution, 1000 mg/L SO <sub>4</sub> (NIST), 500 mL 1000 mg/L SO <sub>4</sub> |                              | -                    | *                     |
| Sulphate               | 257849            | Sulphate standard solution, 50 mg/L SO <sub>4</sub> (NIST), 500 mL                             | 50 mg/L SO <sub>4</sub>      | -                    | *                     |
| Surfactants, non-ionic | LCA333            | ADDISTA Surfactants standard for LCK333 1g/L TRITON x 100                                      | 1 g/L TRITON x 100           | -                    | *                     |
| Varies                 | 244932            | Sulphuric acid standard solution, 5.25 N, 100 mL   | 5.25 N                       | С                    | *                     |
| Varies                 | 20253             | Sulphuric acid standard solution, 0.100 N, 1 L   | 0.100 N                      | -                    | *                     |
| Varies                 | 2332453           | Sodium hydroxide standard solution, 6 N, 1 L   | 6.0 N                        | С                    | *                     |
| Varies                 | 2339349           | Sulphuric acid 0.04 N, 500 mL  | 0.04 N                       | -                    | *                     |
| Varies                 | 28249             | Potassium hydroxide standard solution, 8.00 N, 500 mL  | 8.00 N                       | С                    | *                     |

<sup>\*:</sup> GHS hazard code will be available in the future







# **Sample Preparation**




A selection of sample preparation accessories to photometric analysis for the purpose of digestion, filtration, homogenisation, and dilution.

| Product description   | Article number | EU hazard code | GHS hazard code            |
|---|----------------|----------------|----------------------------|
| Dilution water, organic free, 500 mL                                  | 2641549        | -              | *                          |
| Chloride test strips, low range, 30 - 600 mg/L, 40 pcs                | 2744940        | -              | *                          |
| Chloride test strips, 300 - 6000 mg/L, 0.05 - 1.0 % NaCl , 40 tests   | 2751340        | -              | *                          |
| CRACK SET Reagent set for metal digestions                            | LCW902         | Xn, C, O       | GHS03, GHS05, GHS07, GHS08 |
| Calcium separation set  | LCW903         | -              | GHS07                      |
| Membrane filtration set with 50 membrane filters 1.2 μm               | LCW904         |                |                            |
| Screening test for organic complexing agents                          | LCW907         | -              | GHS05                      |
| Digestion solution for chloride in concrete                           | LCW908         | С              | GHS05                      |
| Total Kjeldahl nitrogen, reagents for digestion                       | LCW909         | С              | GHS05                      |
| Nitrification inhibitor for BOD <sub>5</sub> , acc. DIN 38409-51, 35g | LCW910         | -              | *                          |
| Powder dispenser  | LCW912         |                |                            |
| Membrane filtration set with 50 membrane filters 0.45 μm              | LCW916         |                |                            |
| CleanUp set for cyanide cuvette test LCK319                           | LCW923         | Xn             | GHS08                      |
| Chloride elimination set  | LCW925         | C, O           | GHS03, GHS05               |
| Set for digestion of total silver                                     | LCW954         | Xn, O          | GHS03, GHS07, GHS08        |
| Magnetic stirrer rods, 3 pieces                                       | LYW064         |                |                            |
| Chromium digestion for highly loaded samples                          | LYW513         | Xn, O          | GHS03, GHS07, GHS08        |
| Magnetic stirrer, 0 - 1500 rpm  | LYW854         |                |                            |
| Timer clock   | LZC902         |                |                            |

 $<sup>\</sup>ensuremath{^{\star}}\xspace$  GHS hazard code will be available in the future

<sup>-:</sup> product is not subject to classification

# **Accessories for Cuvette Tests**



A selection of accessories for determination of AOX, BOD, cyanide, organic acids, phenols, sludge activity, surfactants and TOC.

 ${\tt TOC\ Shaker\ for\ purging\ the\ inorganic\ carbon\ (TIC)\ to\ determine\ TOC\ with\ LCK385,\ LCK386\ and\ LCK387.}$ 

| Product description  | Article number  |
|--|-----------------|
| AOX  |                 |
| CARBODISK Active carbon disks for the AOX reference analysis                                   | LZC910          |
| Magnetic stirrer, 0 - 1500 rpm   | LYW854          |
| BOD  |                 |
| $\operatorname{BioKit}$ for $\operatorname{BOD}_5$ cuvette test, as inoculation mat., 20 tests | LZC555          |
| BOD <sub>5</sub> dilution water set  | LZC901          |
| Set of reaction glasses with caps, 60 pieces   | LZC924          |
| $\label{eq:AquaKit} AquaKit \ for \ BOD_{\scriptscriptstyle{5}} \ dilution \ water \ set$      | LZC955          |
| Reaction vessels with screw caps, 20 mm diameter, 5 pieces                                     | LZP065          |
| Funnel   | EBT006          |
| Beaker 150 mL  | HBG011          |
| Magnetic stirrer, 0 - 1500 rpm   | LYW854          |
| LT 20 BOD <sub>5</sub> thermostat  | LTV073          |
| Cyanide, organic acids and phenols   |                 |
| MICRO DIST Block, digital, complete  | MDI001          |
| MICRO DIST Digestion tubes, user-fill, 10 pieces   | A17017          |
| MICRO DIST Digestion tubes, user-fill, 100 pieces  | A17117          |
| MICRO DIST Digestion tubes, user-fill, 50 pieces   | A17517          |
| MICRO DIST Cap press   | 17023L          |
| Rack collector 24 positions  | 17012           |
| Sample rack, 60 positions for 16 mm tubes - Al & ASX-500 samplers                              | 21302           |
| Sludge activity  |                 |
| Accessory kit sludge activity  | LZC918          |
| Membrane filtration set with 50 membrane filters 1.2 $\mu m$                                   | LCW904          |
| Surfactants  |                 |
| LS 120 Shaker for surfactant analysis  | LQV148.99.10001 |
| TOC  |                 |
| Powder dispenser   | LCW912          |
| Membrane filtration set with 50 membrane filters 0.45 $\mu m$                                  | LCW916          |
| TOC-X5 TOC Shaker for purging method   | LQV148.99.00001 |
|  |                 |



## **Accessories**



A selection of cuvette consumables, glass vessels, pipettes, pipette tips, and safety accessories.

Just comfortable: TENSETTE plus electronic pipette complete with rack, storage battery and power supply, 0.2 - 5.0 mL volume.

| Product description   | Article number |
|---|----------------|
| Cuvette consumables   |                |
| Rack for 16 HACH LANGE round cuvettes or rectangular cuvettes 10 mm                   | LYW915         |
| Rack for 8 reaction vessels   | LYW918         |
| Rack for 7 cuvettes with layer thickness of 50 mm                                     | ETS016         |
| Cooling rack for 8 tubes (COD tubes, 16 mm)   | 1864100        |
| Rack, test tube, Polyethylene, 30 mm (O.D.), 21 holes                                 | 2497904        |
| Stopper, Neoprene, solid, size 2, 12 pieces   | 1480802        |
| Stopper for 18 mm glass viewing tube, 6 pieces  | 173106         |
| Disposable paper tissues, white, 200 pieces   | EZZ073         |
| Blank value cuvette set   | LCW919         |
| Glass vessels   |                |
| Beaker 150 mL   | HBG011         |
| Reaction vessels with screw caps, 20 mm diameter, 5 pieces                            | LZP065         |
| Volumetric flask 50 mL, class A wide neck, NS12/21 PP-stopper, 2 pieces               | LZP141         |
| Volumetric flask 100 mL, class A, NS 14/23, PP-stopper, 2 pieces                      | LZP142         |
| Graduated cylinder 50:1 mL, tall form, class B, 2 pieces                              | LZP143         |
| Graduated cylinder 100 mL: 1 mL, tall form, class B, 2 pieces                         | LZP144         |
| Set of reaction glasses with caps, 60 pieces  | LZC924         |
| Pipettes, pipette tips  | LEGGET         |
| Pipettes, pipette tips Pipette tips 1.0-5.0 mL for variable volume pipette, 75 pieces | BBP068         |
| Rack for 5 pipettes   | LYW964         |
| TENSETTE plus electronic pipette  | BBP087         |
| Pipette, fixed volume, 1.0 mL   | BBP163         |
| Pipette, fixed volume, 2.0 mL   | BBP164         |
| Pipette, fixed volume, 0.1 mL   | LYW785         |
| Pipette, fixed volume, 0.2 mL   | LYW790         |
| Pipette tips 0.2-1.0 mL   | BBP079         |
| Pipette tips 0.2-5.0 mL for electronic pipette, 10 pieces                             | LYW250         |
| Pipette tips 0.1 and 0.2 mL   | LYW786         |
| Pipette tips 1.0-5.0 mL   | LYW787         |
| Pipette tips 0.2-1.0 mL   | LYW788         |
| Pipette, variable, volume 1.0-5.0 mL  | BBP065         |
| Pipette, variable, volume 0.2-1.0 mL  | BBP078         |
| Set of 2 pipettes, variable volume, incl. tips  | LZP320         |
|   |                |
| Pipette validation kit  | LCA722         |
| Safety accessories  | F77004         |
| Safety goggles, transparent, DIN 582, suitable for spectacle wearers                  | EZZ031         |
| Safety goggles Uvex according DIN 58211, green / purple                               | EZZO42         |
| Adhesive tape width 75 mm, for transportation of hazardous materials                  | HYB008         |
| Protective gloves size L, blue, nitrile, powder-free, 50 pieces                       | SM743L         |
| Protective gloves size 7 (M), blue, nitrile, powder-free, 50 pieces                   | SM743M         |
| Single use latex gloves size 7 (M), powder-free, green, 100 pieces                    | SM995417       |
| Single use latex gloves size L, powder-free, green, 100 pieces                        | SM995418       |

## **Photometers**

HACH LANGE photometers consistently break innovation barriers to deliver top performance for both routine laboratory tasks and demanding photometry applications. Our photometers are engineered to achieve exceptional analytical accuracy in a simple testing format.



#### Fully automated water analysis

Including sample preparation, the laboratory robot AP 3900 even processes the critical parameters of COD, total P and total N in parallel.

#### Spectrophotometers for the lab

The high-performance VIS and UV-VIS spectrophotometers DR 3900 & DR 6000 deliver reliable and traceable measurement results for both routine laboratory analysis and user applications. With RFID technology, IBR+, AQA+ and LINK2SC. Detailed explanations can be found on page 5.

#### Portable colorimeters for field testing

Designed for use in the field, the handheld outdoor colorimeters DR 900 and POCKETs with battery operation are pre-programmed, easy to use and have a rugged construction.



26 PHOTOMETERS www.hach-lange.com

# **Benchtop Photometer Finder**

## Quick reference guide





## NEW

|                             | DR 6000<br>UV-VIS spectrophotometer with RFID technology   | DR 3900<br>Spectrophotometer with RFID technology                 |  |  |
|-----------------------------|--|---|--|--|
| IBR+                        | Automatic test recognition, lot  | control and expiry date check                                     |  |  |
| Specific technology         | RFID for easy method update, sample ID and Certificate of Analysis   |   |  |  |
| LINK2SC                     | Data exchange with   | SC 1000 controller  |  |  |
| Quality assurance           | Function to schedule and docur   | nent QA with pass/fall indication                                 |  |  |
| Pre-programmed methods      | > 240  | > 220   |  |  |
| Cuvette compatibility       | Rectangular: 10, 20, 30, 50 mm, 1 inch; round: 13 mm, 1 inch<br>Optional 100 mm rectangular cell with additional adapter | Rectangular: 10, 20, 30, 50 mm, 1 inch; round: 13 mm, 1 inch      |  |  |
| Display                     | 7" TFT WVGA co   | lour touchscreen  |  |  |
| Operating mode              | Transmittance (%), Absorbanc   | e and Concentration, Scanning                                     |  |  |
| Wavelength range            | 190 - 1100 nm  | 320 - 1100 nm   |  |  |
| Photometric measuring range | ± 3 Abs (wavelength range 340 to 900 nm)   |   |  |  |
| Photometric accuracy        | 5 mAbs at 0  | .0 to 0.5 Abs   |  |  |
|                             | 1 % at 0.50 to 2.0 Abs   |   |  |  |
| Wavelength resolution       | 0.1 nm   | 1 nm  |  |  |
| Spectral bandwidth          | 2 nm   | 5 nm  |  |  |
| Optical system              | Reference be   | e beam, spectral  |  |  |
| Source lamp                 | Tungsten (VIS), Deuterium lamp (UV)  | Tungsten (VIS)  |  |  |
| Printer compatibility       | Supports most offi   | ce deskjet printers   |  |  |
| Data storage                | 5000 measured values (Result, Date, Time, Sample ID, Operator ID)  | 2000 measured values (Result, Date, Time, Sample ID, Operator ID) |  |  |
| Interfaces                  | USB type A (2), USB type   | B, Ethernet, RFID module  |  |  |
| Power supply                | Power cord, 100 - 240 V, 50 - 60 Hz  | External power supply, 100 - 240 V, 50 - 60 Hz                    |  |  |
| Dimensions (H x W x D)      | 215 mm x 500 mm x 460 mm   | 151 mm x 350 mm x 255 mm  |  |  |
| Weight                      | 11 kg  | 4.2 kg  |  |  |
|                             | Detailed description on pages 30 and 31  | Detailed description on pages 32 and 33                           |  |  |

## Portable Photometer Finder

## Quick reference guide





## NEW

|   | IVEVV  |   |  |  |
|---|--|---|--|--|
|   | DR 900   | POCKET Colorimeter II   |  |  |
|   | Multi-parameter colorimeter  | Single parameter colorimeter  |  |  |
| Supported chemistry                         | HACH tests   | HACH tests and HACH LANGE cuvette tests   |  |  |
| Operating mode                              | Transmittance (%), Absorbance and Concentration  |   |  |  |
| Source lamp                                 | Light Emitti   | ng Diode (LED)  |  |  |
| Wavelength range                            | 420, 520, 560, 610 nm  | varies with model   |  |  |
| Photometric measuring range                 | 0 -  | 2 Abs   |  |  |
| Wavelength accuracy                         | $\pm$ 1 nm (fixed, varies with model)  | Fixed wavelength $\pm 2 \text{ nm}$ varies with model   |  |  |
| Spectral bandwidth                          | 15 nm filt   | ter bandwidth   |  |  |
| Photometric accuracy                        | ± 0.005 Abs @1.0 ABS Nominal   |   |  |  |
| Photometric linearity                       | ± 0.002 Abs (0 - 1 Abs)  |   |  |  |
| Wavelength selection                        | Automatic  | Fixed wavelength  |  |  |
| Stray light                                 | < 1.0 % at 400 nm  |   |  |  |
| Display                                     | Graphical display 240 x 160 pixel (Backlit)  | LCD, backlit  |  |  |
| User programmes                             | 10   | Custom programming 1  |  |  |
| Data storage                                | 500 measured values (result, date, time, sample ID, user ID according to GLP)  | 10 measured values + time   |  |  |
| Cuvette compatibility                       | 1 inch round / 16 mm round (with adapter)  | 1 inch round / 13 mm round /1cm square (with optional adapter)  |  |  |
| Dimensions (H x W x D)                      | 231 mm x 96 mm x 48 mm   | 155 mm x 61 mm x 35 mm  |  |  |
| Weight                                      | 0.6 kg with battery  | 0.23 kg   |  |  |
| Environmental conditions: temperature       | 10   | - 40 °C   |  |  |
| Environmental conditions: relative humidity | max. 90 % relative humidity (non-condensing)   |   |  |  |
| Battery requirements                        | 4 AA size batteries  | 4 AAA size batteries  |  |  |
| Battery life                                | 6 months (typical) at 5 readings a day / 5 day week without backlight (backlight usage will decrease battery life)   | 2000 tests * backlight will decrease battery life   |  |  |
| Interface                                   | USB type Mini IP67   | None  |  |  |
| Enclosure waterproof rating                 | IP67   | IP67 (excludes battery compartment)   |  |  |
| User interface                              | English, French, German, Italian, Spanish, Portuguese, Bulgarian,<br>Chinese, Czech, Danish, Dutch, Finnish, Greek, Hungarian, Japanese,<br>Korean, Polish, Romanian, Russian, Slovenian, Swedish, Turkish   | Numeric   |  |  |
| Includes                                    | DR 900 Colorimeter, two 1-inch glass sample cells marked at 10, 20 and 25 mL, two 1 cm plastic sample cells, 1 x 16-mm COD/ TEST'N'TUBE adapter, 4 AA alkaline batteries, instrument manual multilingual printed, instrument and procedure manuals on CD; USB Mini to USB cable, European CE mark. | PC II, sample cells, manual, carry case. Reagents included: 100 tests each for low range or 50 tests each for high range. |  |  |
|   | Detailed description on pages 34 and 35  | Detailed description on pages 36 and 37   |  |  |



28 PHOTOMETERS www.hach-lange.com

# AP 3900: Laboratory robot for fully automated water analysis

Laboratory robot for water analysis including sample preparation. Modular concept. Basic version contains COD, total P, total N, Ammonium, Nitrate and Nitrite.

#### **NEW**



- ► Saves time and costs
- ► Increases productivity and flexibility
- ► Highest precision and accuracy due to automated procedures
- ► Parallel execution of different samples and methods
- ► Reliable by complete traceability of results
- ▶ Uses DR 3900 as detector
- ► Cost effective for 30 or more tests per day

This unique product processes the critical parameters of COD, total P and total N in parallel using our well established, pre-programmed cuvette tests. The control software ensures the optimal sequence for processing all samples to minimise total time to results through sample preparation, digestion of complex samples, waiting times and measurement. Additional samples can be added at any time, even when the sequence is running and the current status of the analysis is accessible any time with a simple mouse click. Rapid yet simple even untrained users are able to enter all necessary information to the system due to the easy-to-use software.





AllPhred demonstrates laboratory automation with AP 3900: The QR code leads you to an animation.

#### Technical data

| Number of cuvette positions<br>160  | <b>Dosing system (reagent)</b><br>Reagent - pipette tips  | <b>Compressed air pressure</b> 5 bar  |  |  |
|---|---|---|--|--|
| Number of heating positions<br>2 x 24 (optional 2 x 48)                                       | <b>Calibration</b><br>Range 0.2 - 2.0 mL  | <b>Power requirements (Hz)</b> 50/60 Hz   |  |  |
| Number of reagent positions 12  | Measurement method Automatic HACH LANGE cuvette test (13 mm test tube); 10 times measurement and 2D barcode | Power requirements (Voltage) 230 V AC  Dimensions (H x W x D) 950 mm x 1290 mm x 840 mm  Temperature Selectable 40°C, 100°C, 110°C, 148°C and 150°C |  |  |
| Number of sample positions<br>24 (optional 48); 50 mL sample volume                           | Detector DR 3900  |   |  |  |
| <b>Dispenser</b> Calibrated Hamilton Dispenser 2.5 mL   | Photometric accuracy<br>1% at 0.5 - 2.0 E   |   |  |  |
| <b>Dosing system (sample)</b> Sample - PTFE sheathed needle, ID 2 mm Stirrer with 9 mm paddle | Photometric linearity < 0.5 % - 2 E   | Subject to change without notice.   |  |  |

## Specially designed for AP 3900

### APC chemistry - exclusively suitable for the laboratory robot

| Article<br>number | Parameter                 | Measuring range                     | Method   | According to standard                  | Quality<br>control | Number<br>of tests | EU<br>hazard<br>code | GHS<br>hazard<br>code               |
|-------------------|---------------------------|-------------------------------------|--|--|--------------------|--------------------|----------------------|-------------------------------------|
| APC303            | Ammonium                  | 2 - 47 mg/L NH <sub>4</sub> -N      | Indophenol Blue  | ISO 7150-1,<br>DIN 38406 E5-1          | LCA703             | 100                | Xn, N                | GHS05,<br>GHS07,<br>GHS09           |
| APC304            | Ammonium                  | 0.015 - 2.0 mg/L NH <sub>4</sub> -N | Indophenol Blue  | ISO 7150-1,<br>DIN 38406 E5-1          | LCA700             | 100                | Xn, N                | GHS05,<br>GHS07,<br>GHS09           |
| APC114            | COD                       | 150 - 1000 mg/L 0 <sub>2</sub>      | Dichromate   | ISO 6060-1989,<br>DIN 38409-H41-H44    | LCA703             | 100                | T, C                 | GHS05,<br>GHS06,<br>GHS08,<br>GHS09 |
| APC314            | COD                       | 15 - 150 mg/L 0 <sub>2</sub>        | Dichromate   | ISO 6060-1989,<br>DIN 38409-H41-H44    | LCA704             | 100                | T, C                 | GHS05,<br>GHS06,<br>GHS08,<br>GHS09 |
| APC500            | COD                       | 0 - 150 mg/L 0 <sub>2</sub>         | Dichromate   | ISO 15705                              | LCA704             | 100                | T, C                 | GHS05,<br>GHS06,<br>GHS08,<br>GHS09 |
| APC400            | COD                       | 0 - 1000 mg/L 0 <sub>2</sub>        | Dichromate   | ISO 15705                              | LCA703             | 100                | T, C                 | GHS05,<br>GHS06,<br>GHS08,<br>GHS09 |
| APC339            | Nitrate                   | 0.23 - 13.5 mg/L NO <sub>3</sub> -N | 2.6-Dimethylphenol   | ISO 7890-1-2-1986,<br>DIN 38405 D9-2   | LCA703             | 100                | С                    | GHS02,<br>GHS05,<br>GHS07           |
| APC340            | Nitrate                   | 5 - 35 mg/L NO <sub>3</sub> -N      | 2.6-Dimethylphenol   | ISO 7890-1-2-1986,<br>DIN 38405 D9-2   | LCA704             | 100                | С                    | GHS02,<br>GHS05                     |
| APC341            | Nitrite                   | 0.015 - 0.6 mg/L NO <sub>2</sub> -N | Diazotisation  | EN ISO 26777, DIN 38405 D10            | LCA707             | 100                | Xi                   | GHS07                               |
| APC342            | Nitrite                   | 0.6 - 6.0 mg/L NO <sub>2</sub> -N   | Diazotisation  | EN ISO 26777, DIN 38405 D10            | LCA709             | 100                | Xi                   | GHS07                               |
| APC138            | Nitrogen total<br>(LATON) | 1 - 16 mg/L TN <sub>b</sub>         | Koroleff Digestion<br>(Peroxodisulphate), and Photometric<br>Detection with 2.6-Dimethylphenol | EN ISO 11905-1                         | LCA709             | 50                 | T, C                 | GHS02,<br>GHS05,<br>GHS07,<br>GHS08 |
| APC238            | Nitrogen total<br>(LATON) | 5 - 40 mg/L TN <sub>b</sub>         | Koroleff Digestion<br>(Peroxodisulphate), and Photometric<br>Detection with 2.6-Dimethylphenol | EN ISO 11905-1                         | LCA700             | 50                 | T, C                 | GHS02,<br>GHS05,<br>GHS07,<br>GHS08 |
| APC338            | Nitrogen total<br>(LATON) | 20 - 100 mg/L TN <sub>b</sub>       | Koroleff Digestion<br>(Peroxodisulphate), and Photometric<br>Detection with 2.6-Dimethylphenol | EN ISO 11905-1                         | LCA708             | 50                 | T, C                 | GHS02,<br>GHS05,<br>GHS07,<br>GHS08 |
| APC348            | Phosphate                 | 0.5 - 5.0 mg/L PO <sub>4</sub> -P   | Phosphormolybdenum Blue  | EN ISO 6878-1-1986,<br>DIN 38405 D11-4 | LCA700,<br>LCA707  | 100                | С                    | GHS05,<br>GHS07,<br>GHS08           |
| APC349            | Phosphate                 | 0.05 - 1.5 mg/L PO <sub>4</sub> -P  | Phosphormolybdenum Blue  | ISO 6878-1-1986,<br>DIN 38405 D11-4    | LCA704,<br>LCA709  | 100                | С                    | GHS05,<br>GHS07,<br>GHS08           |
| APC350            | Phosphate                 | 2 - 20 mg/L PO <sub>4</sub> -P      | Phosphormolybdenum Blue  | ISO 6878-1-1986,<br>DIN 38405 D11-4    | LCA703,<br>LCA708  | 100                | С                    | GHS05,<br>GHS07,<br>GHS08           |

APC chemistry was specially developed for the AP 3900 laboratory robot and only runs on this instrument. Please note: APC400 and APC500 require reagent blanks. For these, the number of tests varies. Hazard code descriptions: see page 19

### **Order information**

| Article number | Product description   |   |
|----------------|---|---|
| SMAP3900-MULTI | AP 3900 Multi laboratory robot standard configuration               | Laboratory robot for water analysis including sample preparation.<br>Modular concept. Basic version contains COD, total P, total N,<br>Ammonium, Nitrate and Nitrite. |
| SMUPGRADE-24F  | Filtration upgrade incl. software for the AP 3900 (24 pos)          |   |
| SMUPGRADE-48S  | AP 3900 multi upgrade kit to 48 sample positions                    |   |
| SMUPGRADE-96H  | Heater upgrade incl. software for the AP 3900 (2 x 48 pos)          |   |
| LZX521         | Compressor with connection kit, 2 x fitting 1/4" outer , 5 m tubing |   |



30 PHOTOMETERS www.hach-lange.com

# DR 6000: Combining quality and efficiency in the professional laboratory

The UV-VIS spectrophotometer delivers top performance for both routine laboratory tasks and demanding photometry applications.

#### NEW



- ► Improved laboratory efficiency more than 240 pre-programmed methods directly available
- ► Comparable and reliable results with the approved HACH LANGE cuvette tests
- ► Transparent working processes in every situation with access to all raw data
- ► Integrated Quality Assurance with function for scheduling, evaluation, documentation
- ▶ Optimised data management, LIMS compatible
- ► Traceability back to the sampling point by means of RFID technology

The new UV-VIS spectrophotometer is engineered and manufactured in Germany to deliver in fourth generation exceptional analytical accuracy. The Czerny-Turner monochromator design reduces aberrations and guarantees a minimal spectral bandwidth. The output coupler mirror optimally aligns the measurement beam.

Four sequential band-pass filters reduce internal scattered light. The reference beam technology compensates for signal fluctuations in the instrument. Two low-noise silicon detectors ensure high selectivity and basic stability of the measurement signal.

The UV-VIS spectrophotometer unites reliable results with efficiency. The intuitive menu navigation with colour touch screen allows you to enter and calibrate your own methods in just a few simple steps. The instrument provides a wide range of pre-programmed methods. Application packages, e.g. for enzymology and colorimetry, open up further application opportunities including drinking water and brewery analysis.



Application solftware available for drinking water, breweries and food industry.











## Technical data

| IBR+ Automatic test recognition, lot control and expiry date check  | Operating mode Transmittance (%), Absorbance and Concentration, Scanning   | <b>Printer compatibility</b> Supports most office deskjet printers                          |
|---|--|---|
| Specific technology RFID for easy method update, sample ID and Certificate of Analysis  | Wavelength range<br>190 - 1100 nm  | <b>Data storage</b><br>5000 measured values (Result, Date, Time, Sample ID,<br>Operator ID) |
| LINK2SC Data exchange with SC 1000 controller   | Photometric measuring range ± 3 Abs (wavelength range 340 to 900 nm)       | Interfaces USB type A (2), USB type B, Ethernet, RFID module                                |
| Quality assurance Function to schedule and document QA with pass/fall indication  | Photometric accuracy<br>5 mAbs at 0.0 to 0.5 Abs<br>1 % at 0.50 to 2.0 Abs | Power supply Power cord, 100 - 240 V, 50 - 60 Hz  |
| Pre-programmed methods > 240  | <b>Wavelength resolution</b> 0.1 nm  | <b>Dimensions (H x W x D)</b> 215 mm x 500 mm x 460 mm                                      |
| Cuvette compatibility Rectangular: 10, 20, 30, 50 mm, 1 inch; round: 13 mm, 1 inch Optional 100 mm rectangular cell with additional adapter | Spectral bandwidth 2 nm  | Weight<br>11 kg   |
|   | <b>Optical system</b><br>Reference beam, spectral                          |   |
| <b>Display</b><br>7" TFT WVGA colour touchscreen  | Source lamp<br>Tungsten (VIS), Deuterium lamp (UV)                         | Subject to change without notice.   |

## **Order information**

| Article number  | Product description                                      |  |
|-----------------|--|--|
| LPV441.99.00011 | DR 6000 UV-VIS spectrophotometer with RFID technology    | The UV-VIS spectrophotometer delivers top performance for both routine laboratory tasks and demanding photometry applications.   |
| LQV156.99.10011 | LOC 100 RFID set for sample identification               | The set contains: 1 RFID locator LOC 100, 15 sample RFID tags in 5 colours, 5 location RFID tags and 2 operator RFID tags.   |
| LQV157.99.30001 | SIP 10 Sipper set for DR 6000 with 1 cm quartz cell      | Sipper set for flow through applications in UV range. 1 cm flow through cell, quartz.  |
| LQV157.99.20001 | SIP 10 Sipper set for DR 6000 with 1 inch round cell     | Sipper for pour through applications with spectrophotometer DR 6000. With dual path length 1 inch/cm round cell, USB cable and pump tubing.  |
| TSE-CC-DR6000   | Comfort contract for DR 6000                             | Service Package 2.2: Comfort maintenance contract for spectrophotometer DR 6000. Excludes travel and wear parts. Includes warranty extension.  |
| TSE-BC-DR6000   | Basic contract for DR 6000                               | Service Package 2.1: Basic maintenance contract for spectrophotometer DR 6000. Excludes travel and wear parts.   |
| TSE-MC-DR6000   | One-off inspection for DR 6000                           | Service Package 1: One-off inspection for spectrophotometer DR 6000. Excludes travel and wear parts.   |
| TSE-IQOQ-DR6000 | Equipment qualification IQOQ for DR 6000                 | Equipment qualification IQOQ for UV-VIS spectrophotometer DR 6000. Installation and operational qualification. Including required documents, excluding travel and standards.                           |
| LZV935          | DR 6000 application software for drinking water analysis | The drinking water analysis software LZV925 is a compilation of all spectrophotometric applications that are relevant for drinking water analysis.   |
| LZV936          | DR 6000 application software for brewery analysis        | The additional brewery analysis software LZV936 is a compilation of all spectrophotometric applications that are relevant for brewery analysis.  |
| LZV937          | DR 6000 application software for enzymatic food analysis | The additional software LZV937 is a compilation of enzymology tests manufactured by R-Biopharm AG, Darmstadt that can be performed with spectrophotometer DR 6000 and carousel insert LZV902.99.00001. |
| LZV938          | Remote operating software photometer                     | The operating software for the photometer allows the instrument to be controlled remotely using a PC.  |
| LZV902.99.00001 | Carousel holder 1 cm for DR 6000                         | Cell holder with 7 positions to measure mini series or ezymatic methods.   |



32 PHOTOMETERS www.hach-lange.com


# DR 3900: Accuracy from start to finish

High-performance VIS spectrophotometer with RFID technology for reliable and traceable measurement results of routine analysis and user applications.



- ► Traceability starts with sampling Samples are encoded and identified with RFID
- ► IBR+ increases the reliability of your measurement values A 2D barcode on the cuvette delivers lot number and expiry date
- ► Rapid data updates
  The RFID labeling allows a touchless data transfer
- ► Quality assurance made easy with AQA+ Definition and documentation of QA procedures, retrieval of Certificates of Analysis
- ► Alignment of laboratory and process analysis LINK2SC - adjustment of process on-line value and lab reference
- ► Data transfer is simple via USB or Ethernet

Compact and reliable VIS spectrophotometer with reference beam technology. Samples are traced back to sample location due to RFID technology. Lot number and expiry date information of reagents are now included on the 2D barcode. The RFID module reads out all batch specific information like factors, updated methods and the current batch certificate from cuvette test box. All information can be retrieved immediately on the spectrophotometer and printed out. Process results can be compared to laboratory reference values in the photometer via LINK2SC connection between SC controller and photometer. Data can be exchanged bi-directionally via Ethernet, i.e. matrix corrections of process probes can be done directly from the laboratory.



During the rotating ten times measurement process using the IBR+ barcode reader, the DR 3900 immediately picks up all the information on the cuvette.











## Technical data

| IBR+ Automatic test recognition, lot control and expiry date check                     | Operating mode Transmittance (%), Absorbance and Concentration, Scanning   | Source lamp<br>Tungsten (VIS)   |
|--|--|---|
| Specific technology RFID for easy method update, sample ID and Certificate of Analysis | Wavelength range<br>320 - 1100 nm  | <b>Printer compatibility</b> Supports most office deskjet printers                          |
| LINK2SC Data exchange with SC 1000 controller  | Photometric measuring range ± 3 Abs (wavelength range 340 to 900 nm)       | <b>Data storage</b><br>2000 measured values (Result, Date, Time, Sample ID,<br>Operator ID) |
| Quality assurance Function to schedule and document QA with pass/fail indication       | Photometric accuracy<br>5 mAbs at 0.0 to 0.5 Abs<br>1 % at 0.50 to 2.0 Abs | Interfaces USB type A (2), USB type B, Ethernet, RFID module                                |
| Pre-programmed methods > 220   | Wavelength resolution<br>1 nm  | <b>Power supply</b> External power supply, 100 - 240 V, 50 - 60 Hz                          |
| Cuvette compatibility Rectangular: 10, 20, 30, 50 mm, 1 inch; round: 13 mm, 1 inch     | Spectral bandwidth<br>5 nm   | <b>Dimensions (H x W x D)</b> 151 mm x 350 mm x 255 mm                                      |
|  | Optical system Reference beam, spectral                                    | <b>Weight</b><br>4.2 kg   |
| <b>Display</b> 7" TFT WVGA colour touchscreen  |  | Subject to change without notice.   |

## Order information

| Article number  | Product description                                  |   |
|-----------------|--|---|
| LPV440.99.00001 | DR 3900 Spectrophotometer with RFID technology       | High-performance VIS spectrophotometer with RFID technology for reliable and traceable measurement results of routine analysis and user applications.                                     |
| LPV440.99.10001 | DR 3900 RFID spectrophotometer / LOC 100 kit         | High-performance VIS spectrophotometer with RFID technology for reliable and traceable measurement results of routine analysis and user applications; sample identification set included. |
| LQV156.99.10011 | LOC 100 RFID set for sample identification           | The set contains: 1 RFID locator LOC 100, 15 sample RFID tags in 5 colours, 5 location RFID tags and 2 operator RFID tags.  |
| LQV157.99.10001 | SIP 10 Sipper set for DR 3900 with 1 inch round cell | Sipper for pour through applications with spectrophotometer DR 3900. With dual path length 1 inch/cm round cell, USB cable and pump tubing.   |
| TSE-CC-DR3900   | Comfort contract for DR 3900                         | Service Package 2.2: Comfort maintenance contract for spectrophotometer DR 3900. Excludes travel and wear parts. Includes warranty extension.   |
| TSE-BC-DR3900   | Basic contract for DR 3900                           | Service Package 2.1: Basic maintenance contract for spectrophotometer DR 3900. Excludes travel and wear parts.  |
| TSE-MC-DR3900   | One-off inspection for DR 3900                       | Service Package 1: One-off inspection for laboratory spectrophotometer DR 3900. Excludes travel and wear parts.   |
| TSE-IQOQ-DR3900 | Equipment qualification IQOQ for DR 3900             | Equipment qualification IQOQ for VIS spectrophotometer DR 3900. Installation and operational qualification. Including required documents, excluding travel and standards.                 |



34 PHOTOMETERS www.hach-lange.com

# DR 900: Robust outdoor handheld Colorimeter

Portable and robust, micro-processor-controlled colorimeter with power-saving LED technology. Pre-programmed with 90 HACH Methods.

#### NEW



- ► Rugged construction

  Dustproof, waterproof, shockproof
- ► Designed for use in the field

  True handheld analysis for use anywhere
- ► Easy to use Menu driven, step-by-step analysis
- ► Reliable results without a main connection Improved user interface allowing quick selection of tests
- ► Configured for immediate use Pre-programmed, ready to use out of the box

The handheld colorimeter saves time in the field by allowing quick and easy access to the most used testing methods in less than four clicks. This colorimeter is waterproof, dustproof, shock resistant, and has been drop tested for greater quality assurance.

This instrument comes with an intuitive user interface, a large data store and a built-in USB port for easy transferring of information. The handheld colorimeter also helps satisfy core testing needs by offering at least 90 of the most common testing methods.

Combining all of these features with a push button backlit display for use in low light areas, you have a handheld colorimeter which is field ready in every possible way, and makes testing in harsh field environments a little less challenging.



DR 900 with opened cover and sample cells



## Technical data

| Operating mode Transmittance (%), Absorbance and Concentration | <b>Display</b> Graphical display 240 x 160 pixel (Backlit)                                 | Battery life 6 months (typical) at 5 readings a day / 5 day week without backlight (backlight usage will decrease battery life)   |
|--|--|---|
| Source lamp<br>Light Emitting Diode (LED)                      | User programmes<br>10  | Interface USB type Mini IP67  |
| <b>Wavelength range</b> 420, 520, 560, 610 nm                  | Data storage 500 measured values (result, date, time, sample ID, user ID according to GLP) | Enclosure waterproof rating   |
| Photometric measuring range<br>0 - 2 Abs                       | Supported chemistry HACH tests   | User interface<br>English, French, German, Italian, Spanish, Portuguese,  |
| Wavelength accuracy<br>± 1 nm (fixed, varies with model)       | Cuvette compatibility 1 inch round / 16 mm round (with adapter)                            | Bulgarian, Chinese, Czech, Danish, Dutch, Finnish, Greek,<br>Hungarian, Japanese, Korean, Polish, Romanian, Russian,<br>Slovenian, Swedish, Turkish                     |
| Spectral bandwidth 15 nm filter bandwidth                      | <b>Dimensions (H x W x D)</b> 231 mm x 96 mm x 48 mm                                       | Includes  DR 900 Colorimeter, two 1-inch glass sample cells marked at 10, 20 and 25 mL, two 1 cm plastic sample cells, 1 x  |
| Photometric accuracy<br>± 0.005 Abs @1.0 ABS Nominal           | Weight 0.6 kg with battery   | 16-mm COD/TEST'N'TUBE adapter, 4 AA alkaline batteries, instrument manual multilingual printed, instrument and procedure manuals on CD: USB Mini to USB cable, European |
| Photometric linearity<br>± 0.002 Abs (0 - 1 Abs)               | Environmental conditions: temperature<br>10 - 40 °C  | CE mark.  |
| <b>Wavelength selection</b> Automatic                          | Environmental conditions: relative humidity max. 90 % relative humidity (non-condensing)   |   |
| <b>Stray light</b> < 1.0 % at 400 nm                           | Battery requirements 4 AA size batteries   | Subject to change without notice.   |

## Order information

| Article number | Product description                            |  |
|----------------|--|--|
| 9385200        | DR 900 Robust portable datalogging colorimeter | Portable and robust, micro-processor-controlled colorimeter with power-saving LED technology. Pre-programmed with 90 HACH Methods.           |
| 4942500        | Case assy, DR 800 and DR 900 colorimeter       | Suitcase for keeping and transport of colorimeter, accessories and reagents.   |
| 2722000        | Case, soft 11.5 H x 2.5 D x 11.5 W             | Instrument carrying case, soft-sided with shoulder strap   |
| 2763900        | DR/Check absorbance standard kit               | Standard set (4 pcs) for checking of the photometric accuracy  |
| TSE-MC-DR900   | One-off inspection for DR 900                  | Service Package 1: One-off inspection for spectrophotometer DR 900. Excludes travel and wear parts.  |
| TSE-CC-DR900   | Comfort contract for DR 900                    | Service Package 2.2: Comfort maintenance contract for spectrophotometer DR 900. Excludes travel and wear parts. Includes warranty extension. |
| TSE-BC-DR900   | Basic contract for DR 900                      | Service Package 2.1: Basic maintenance contract for spectrophotometer DR 900. Excludes travel and wear parts.                                |



36 PHOTOMETERS www.hach-lange.com

# POCKET Colorimeter II: Small in size, big on waterproof analysis

Portable colorimeter programmed for the determination of one or two parameters. In carrying case, complete with reagents, cuvettes and operating instructions.



- ➤ Simple: all functions are available via only four keys
- ► Powerful: battery operation for up to 2000 tests
- ► Clear display: a backlit display and large readout for difficult conditions
- ► Reliable results any place due to battery operation
- ► Rugged construction yet lightweight
- ➤ Waterproof to IP 67: for watertight results!

#### Technical data

| Source lamp<br>Light Emitting Diode (LED)   | <b>Data storage</b> 10 measured values + time  | Battery requirements 4 AAA size batteries                            |
|---|--|--|
| Wavelength range varies with model  | Supported chemistry HACH tests and HACH LANGE cuvette tests                              | Battery life 2000 tests * backlight will decrease battery life       |
| Photometric measuring range<br>0 - 2 Abs  | Cuvette compatibility 1 inch round / 13 mm round /1cm square (with optional adapter)     | Enclosure waterproof rating IP67 (excludes battery compartment)      |
| $\begin{tabular}{ll} \textbf{Wavelength accuracy} \\ \textbf{Fixed wavelength} \pm 2 \text{ nm varies with model} \\ \end{tabular}$ | <b>Dimensions (H x W x D)</b> 155 mm x 61 mm x 35 mm                                     | User interface<br>Numeric  |
| Spectral bandwidth 15 nm filter bandwidth   | Weight<br>0.23 kg  | Includes PC II, sample cells, manual, carry case. Reagents included: |
| <b>Wavelength selection</b> Fixed wavelength  | Environmental conditions: temperature 10 - 40 °C   | 100 tests each for low range or 50 tests each for high range.        |
| <b>Display</b><br>LCD, backlit  | Environmental conditions: relative humidity max. 90 % relative humidity (non-condensing) |  |
| User programmes<br>Custom programming 1   |  | Subject to change without notice.                                    |

### **Order information**

| Article number      | Parameter   | Method  | Range   | Number of tests |
|---------------------|---|---|---|-----------------|
| Kits for HACH LAN   | GE cuvette tests                                      |   |   |                 |
| 5953000V.01         | Ammonium (suitable for LCK303, 304, 305)              | m (suitable for LCK303, 304, 305) Indophenol Blue 2.0 - 47.0 mg/L / 0.015 - 2.0 mg/L / 1.0 - 12.0 mg/L NH <sub>4</sub> -N |   | 25              |
| 5953000V.02         | Chlorine (suitable for LCK310)                        | DPD   | 0.05 - 2.0 mg/L Cl <sub>2</sub>   | 24              |
| 5953000V.03         | COD LR (suitable for LCK314, 614)                     | Dichromate  | 15 - 150 mg/L O <sub>2</sub> /<br>50 - 300 mg/L O <sub>2</sub>                        | 25              |
| 5953000V.04         | COD HR (suitable for LCK014, 114)                     | Dichromate  | 1000 - 10000 mg/L O <sub>2</sub> /<br>150 - 1000 mg/L O <sub>2</sub>                  | 25              |
| 5953000V.05         | Phosphate (suitable for LCK348, 349, 350)             | Phosphormolybdenum Blue   | 0.5 - 5.0 mg/L / 0.05 - 1.5 mg/L PO <sub>4</sub> / 2.0 - 20.0 mg/L PO <sub>4</sub> -P | 25              |
| 5953000V.06         | Zinc (suitable for LCK360)                            | PAR   | 0.2 - 6.0 mg/L Zn   | 24              |
| 5953000V.07         | Chloride (suitable for LCK311)                        | Iron(III)-Thiocyanate   | 1 - 70 mg/L / 70 - 1000 mg/L Cl   | 24              |
| 5953000V.08         | Formaldehyde (suitable for LCK325)                    | Acetylacetone   | 0.5 - 10.0 mg/L H <sub>2</sub> CO   | 24              |
| Kits for HACH tests | 8   | ,   | J 2   |                 |
| 5870025             | Aluminum  | Aluminon  | 0.02 - 0.80 mg/L Al   | 100             |
| 5870026             | Chloramine, Mono and Free Ammonia                     | Indophenol  | 0.02 - 0.50 mg/L NH <sub>3</sub> -N /<br>0.04 - 4.50 mg/L Cl <sub>2</sub>             | 50 - 100        |
| 5870040             | Nitrogen, Ammonia, mid range                          | Salicylate  | 0.01 - 0.80 mg/L NH <sub>3</sub> -N   | 100             |
| 5870001             | Bromine   | DPD   | 0.05 - 4.50/0.2 - 70.0 mg/L Br <sub>2</sub>   | 50 - 100        |
| 5870051             | Chlorine Dioxide                                      | DPD   | 0.05 - 5.00 mg/L CIO <sub>2</sub>   | 100             |
| 5870000             | Chlorine, low range - as free & total Cl <sub>2</sub> | DPD   | 0.02 - 2.00 / 0.1 - 8.0 mg/L Cl <sub>2</sub>  | 50 - 100        |
| 5870012             | pH, and high range Chlorine                           | DPD   | 0.1 - 10.0 Cl <sub>2</sub> / 6.0 - 8.5 pH   | 100             |
| 5870023             | Chlorine, free, SWIFTEST dispenser                    | DPD   | 0.02 - 2.00 mg/L Cl <sub>2</sub> /<br>0.1 - 8.0 mg/L Cl <sub>2</sub>                  | 125 - 250       |
| 5870024             | Chlorine, total, SWIFTEST dispenser                   | DPD   | 0.02 - 2.00/0.1 - 8.0 mg /L Cl <sub>2</sub>   | 125 - 250       |
| 5870017             | Chromium, low range - as Cr (VI)                      | 1,5 Diphenylcarbohydrazide  | 0.01 - 0.70 mg/L Cr   | 100             |
| 5870019             | Copper - as free Cu                                   | Bicinchoninate  | 0.04 - 5.00 mg/L Cu   | 100             |
| 5870005             | Fluoride  | SPADNS  | 0.1 - 2.0 mg/L F  | 50              |
| 5870016             | Iron, low range                                       | TPTZ  | 0.01 - 1.70 mg/L Fe   | 50 - 100        |
| 5870022             | Iron, medium range                                    | FerroVer  | 0.02 - 5.00 mg/L Fe   | 100             |
| 5870021             | Lead  | LeadTrak Fast Column Extraction   | 5 - 150 μg/L Pb   | 20              |
| 5870015             | Manganese   | Periodate Oxidation   | 0.2 - 20.0 mg/L Mn  | 100             |
| 5870018             | Manganese   | PAN   | 0.01 - 0.70 mg/L Mn   | 50              |
| 5870010             | Molybdate - as Molybdenum                             | Ternary Complex   | 0.02 - 3.00/0.1 - 12.0 mg/L Mo  | 100             |
| 5870020             | Nickel & Cobalt                                       | PAN   | 0.01 - 1.00 mg/L Ni /<br>0.02 - 2.00 mg/L Co  | 100             |
| 5870002             | Nitrate-Nitrogen                                      | Cadmium Reduction   | 0.4 - 30.0 mg/L NO <sub>3</sub> -N  | 100             |
| 5870003             | Dissolved Oxygen (DO)                                 | HRD0  | 0.2 - 10.0 mg/L DO  | 25              |
| 5870004             | Ozone   | Indigo Trisulfonate   | 0.01 - 0.25 mg/L O <sub>3</sub>   | 25              |
| 5870006             | Phosphorus, Orthophosphate (reactive)                 | Ascorbic Acid   | 0.02 - 3.00 mg/L PO <sub>4</sub>  | 100             |
| 5870007             | Phosphonate   | Persulfate UV Oxidation   | 0.1 - 2.5/1 - 125 mg/L PO <sub>4</sub>  | 100             |
| 5870034             | Silica, high range                                    | Silicomolybdate  Turbidimetric  | 1 - 100 mg/L SiO <sub>2</sub>   | 100<br>100      |
| 5870029<br>5870009  | Sulphate<br>Zinc                                      | Zincon  | 2 - 70 mg/L SO <sub>4</sub><br>0.02 - 3.00 mg/L Zn                                    | 100             |
|                     |   | ZIIICOII  | 0.02 - 3.00 Hig/E Zii   | 100             |
| Kits for special ap | •   |   |   |                 |
| 5870042             | Colorimeter 420 nm                                    | _   |   |                 |
| 5870045             | Colorimeter 450 nm                                    |   |   |                 |
| 5870050             | Colorimeter 500 nm                                    |   |   |                 |
| 5870052<br>5870055  | Colorimeter 528 nm Colorimeter 550 nm                 | $\dashv$  |   |                 |
| 5870058             | Colorimeter 530 nm                                    | $\dashv$  |   |                 |
| 5870060             | Colorimeter 580 film                                  |   |   |                 |
| 5870065             | Colorimeter 655 nm                                    | _   |   |                 |

Test kits must be ordered separately. Subject to change without notice.



38 PHOTOMETERS www.hach-lange.com

# HT 200S: Fast and cost effective digestion in 15 minutes

Heating block with HSD technology (High Speed Digestion) for extremely fast digestion of samples



- ► Saves time in the analysis of COD, TN<sub>b</sub>, P<sub>tot</sub> and heavy metals
- ► Automatic fast cooling
- ► Variable digestion time and temperature for special digestions
- ► COD results in just 35 minutes

### Technical data

| Heating programmes Pre-programmed for 100°C, HT and COD mode and freely selectable 40-170°C, 5-240 min | Max. operating humidity 90 %   | <b>Power supply</b> 230 V +5%/-15%, 50 Hz, 1300 VA     |
|--|--|--|
| User programmes 9 free temperature/time  | Number of cuvettes<br>12 x 20 mm diameter  | <b>Dimensions (H x W x D)</b> 330 mm x 300 mm x 430 mm |
| Heating rate from 20 °C - 148 °C in 8 minutes  | <b>Display size</b> 2 x 16 characters  | Weight<br>12 kg  |
| Temperature stability ± 1 °C in conformity with EN, ISO, EPA methods                                   | Display type<br>LCD  |  |
| <b>Operating temperature range</b> 10 - 45 °C  | User interface<br>English, French, German, Italian, Spanish, Dutch,<br>Swedish, Polish, Danish | Subject to change without notice.                      |

### **Order information**

| Article number | Product description                 |   |
|----------------|-------------------------------------|---|
| LTV077         | HT 200S High temperature thermostat | Heating block with HSD technology (High Speed Digestion) for extremely fast digestion of samples  |
| OHA104         | Reduction insert for 13 mm cuvettes |   |
| TSE-MC-HT200   | One-off inspection for HT 200S      | Service Package 1: One-off inspection for laboratory high temp. thermostat HT 200S. Excludes travel and wear parts.                               |
| TSE-CC-HT200   | Comfort contract for HT 200S        | Service Package 2.2: Comfort maintenance contract for high temp. thermostat HT 200S. Excludes travel and wear parts. Includes warranty extension. |
| TSE-BC-HT200   | Basic contract for HT 200S          | Service Package 2.1: Basic maintenance contract for high temp. thermostat HT 200S. Excludes travel and wear parts.                                |

# LT 200: Thermostat for standard and special digestions

Pre-programmed for all standard digestions and freely programmable for user specific digestions



- ► Great flexibility
- ► Excellent reproducibility
- ► Simple to use

### Technical data, dual block version

| Toolin Joseph Grand Discourse Control of the Contro |  |   |  |  |
|--|--|---|--|--|
| Heating programmes Pre-programmed for 40°C, 100°C, 148°C and freely selectable from 37-150°C, 1-480 min  |  | <b>Power supply</b><br>115 V - 230 V +5%/-15%, 50 - 60 Hz, 900 VA max |  |  |
| User programmes 6 free temperature/time  | Number of cuvettes 21 x 13 mm diameter 4 x 20 mm diameter                              | <b>Dimensions (H x W x D)</b> 145 mm x 250 mm x 310 mm                |  |  |
| <b>Heating rate</b> from 20 - 148 °C in 10 minutes   | Display size 2 x 16 characters   | Weight<br>2.8 kg  |  |  |
| Temperature stability $\pm$ 1 °C in conformity with EN, ISO, EPA methods   | Display type<br>LCD  |   |  |  |
| Operating temperature range $10$ - $45\ ^{\circ}\text{C}$  | User interface<br>English, French, German, Italian, Spanish, Dutch,<br>Swedish, Polish | Subject to change without notice.                                     |  |  |

### **Order information**

| Article number  | Product description  |   |
|-----------------|--|---|
| LTV082.99.21002 | LT 200 Dry thermostat with 2 blocks, 15 x 13 mm, 6 x 13 mm / 4 x 20 mm | - Illuminated digital display for remaining time /  |
| LTV082.99.10002 | LT 200 Dry thermostat with 1 block, 9 x 13 mm / 2 x 20 mm              | temperature and operator guidance   |
| LTV082.99.23002 | LT 200 Dry thermostat with 2 blocks, 15 x 13 mm, 15 x 13 mm            | - 3 pre-programmed and 6 freely programmable storage spaces                                       |
| LTV082.99.51002 | LT 200 Dry thermostat with 2 blocks, 6 x 13 mm and 4 x 20 mm per block | - Temperature settings between 37 °C and 150 °C in 1 °C steps, time setting between 1 and 480 min |
| LZT144          | Adapter sleeve LT 200 for MICRO DIST                                   |   |



40 PHOTOMETERS www.hach-lange.com

# Which cuvettes for my photometer?

### Quick reference guide

| Article number | Optical path length | Material      | Volume / Package size                                  | PC II | DR 900 | DR 3900 | DR 6000 |
|----------------|---------------------|---------------|--|-------|--------|---------|---------|
| LCW906         | 13 mm round         | Glass         | 25 pcs 7 mL, with rubber caps                          |       |        | •       | -       |
| LZP045         | 10 mm rectangular   | Glass         | 1 pcs 3.5 mL   |       |        |         | -       |
| LZP169         | 50 mm rectangular   | Glass         | 1 pcs 7 mL, semi-micro                                 |       |        |         | •       |
| LZP269         | 50 mm rectangular   | Glass         | 1 pcs 7 mL, semi-micro                                 |       |        |         | •       |
| 2095100        | 10 mm rectangular   | Glass         | 2 pcs 3.5 mL, matched pair, with caps                  |       |        |         | -       |
| 2122800        | 1 inch round        | Glass         | 1 pcs 10 mL, with cap                                  |       |        |         | •       |
| 2401906        | 25 mm round         | Glass         | 6 pcs 25 mL, with caps                                 |       | -      |         |         |
| 2427606        | 1 inch round        | Glass         | 6 pcs 10 mL, with caps                                 |       | •      |         | •       |
| 2495402        | 1 inch square       | Glass         | 2 pcs matched pair                                     |       |        |         |         |
| 2612602        | 1 inch square       | Glass         | 2 pcs 25 mL, matched pair, with caps                   |       |        |         |         |
| 2629250        | 50 mm rectangular   | Glass         | 1 pcs 17,5 mL, with cap                                |       |        |         |         |
| 2665902        | 1 inch square       | Glass         | 2 pcs 25 mL, matched pair                              |       |        |         |         |
| LCW919         | 11 mm round         | Glass         | 5 pcs blank value cuvette set, 7 mL, rubber caps       |       |        |         |         |
| LZP167         | 50 mm rectangular   | Optical glass | 1 pcs 20 mL  |       |        |         | •       |
| LZP331         | 20 mm rectangular   | Optical glass | 1 pcs 7 mL   |       |        |         |         |
| 5940506        | 1 inch round        | Plastic       | 6 pcs 25 mL, 10 mm & 1 dual pathlength, with cap       |       | -      |         | -       |
| LZP341         | 50 mm rectangular   | PMMA          | 10 pcs 7 mL, semi-micro, with caps                     |       |        |         | •       |
| EBK019         | 10 mm rectangular   | Polystyrene   | 1000 pcs 3.5 mL  |       |        |         |         |
| 2410212        | 1 inch rectangular  | Polystyrene   | 12 pcs 25 mL, with caps                                |       |        |         |         |
| 2629500        | 10 mm rectangular   | Polystyrene   | 100 pcs 1.5 mL,  |       |        |         | •       |
| 4864302        | 1 inch round        | Polystyrene   | 2 pcs 10 mL, with caps                                 | •     | -      |         |         |
| LZP332         | 10 mm rectangular   | Quartz glass  | 1 pcs 3.5 mL   |       |        |         | -       |
| LZP333         | 50 mm rectangular   | Quartz glass  | 1 pcs 17.5 mL  |       |        |         | •       |
| A24209         | 10 mm rectangular   | Quartz glass  | 1 pcs 160 $\mu$ L, pour through cell, CH = 10 mm       |       |        |         | -       |
| LZV510         | 10 mm rectangular   | Quartz glass  | 1 pcs 450 $\mu\text{L},$ pour through cell, CH = 10 mm |       |        |         |         |
| LZV649         | 50 mm rectangular   | Quartz glass  | 1 pcs 370 μL, flow through cell                        |       |        |         | -       |
| 2624450        | 50 mm rectangular   | Quartz glass  | 1 pcs 17.5 mL, with cap                                |       |        |         | -       |
| 2624410        | 10 mm square        | Quartz glass  | 1 pcs 3.5 mL, with cap                                 |       |        |         |         |

PC II: Single Parameter Colorimeter, DR 900: Multi-Parameter Colorimeter, DR 3900: VIS Spectrophotometer, DR 6000: UV-VIS Spectrophotometer

# Portable water quality laboratories, **BOD** laboratory measurement stations



**CELs - Factory-configured** laboratories for several applications: Colorimeter, reagent sets, required apparatus, additional instruments and all items needed to conduct field

| Article number | Product description   |
|----------------|---|
| 251231         | Portable colorimeter laboratory for water conditioning      |
| 251232         | Portable colorimeter lab for environmental water quality    |
| 251233         | Portable colorimeter laboratory for aquaculture             |
| 251234         | Basic portable colorimeter laboratory for drinking water    |
| 251235         | Advanced portable colorimeter laboratory for drinking water |
| 251236         | Basic portable colorimeter laboratory for wastewater        |
| 251237         | Advanced portable colorimeter laboratory for wastewater     |
| 251238         | Portable colorimeter lab for professional water treatment   |
| 251239         | Advanced portable colorimeter/pH/conductivity laboratory    |



**DRELs - Factory**configured laboratories for professional water analysis: Photometer with batteries and cuvettes, instrument case, reagent & apparatus case, reagent set, and instrument set.

| Article number | Product description       |
|----------------|---------------------------|
| LZV729         | DREL Complete water lab   |
| LZV735         | DREL Industrial water lab |
|                |                           |
|                |                           |
|                |                           |



**BOD Direct - Respirometric BOD** laboratory measurement station, control unit and stirrer with electronic pressure sensors for up to 6 bottles.

| Article number  | Product description                    |
|-----------------|--|
| LQV158.98.00001 | BOD Direct Respirometric BOD apparatus |
|                 |  |
|                 |  |
|                 |  |



**BOD TRAK II - Manometric** BOD(n) laboratory measurement station, control unit and stirrer with pressure sensors for 6 bottles.

| Article number | Product description                                      |
|----------------|--|
| 2952400        | BOD TRAK II Respirometric BOD apparatus with accessories |
|                |  |
|                |  |

### Instrument quality assurance and documentation

| Product description                             | Article number |
|---|----------------|
| Quality assurance                               |                |
| Test solution set for spectrophotometers        | LZV810         |
| Validation filter kit for spectrophotometer     | LZV537         |
| Pipette validation kit                          | LCA722         |
| Documentation                                   |                |
| USB-A4-Printer for spectrophotometer            | LYW368         |
| USB barcode hand-scanner for spectrophotometers | LZV566         |
| USB keyboard QUERTY                             | LZV582         |



42 TURBIDIMETERS www.hach-lange.com

### **Turbidity instruments**

The HACH LANGE turbidity range comprises portable and benchtop instruments in accordance with DIN EN ISO. The available models for specific requirements can be found on our website.



### **Benchtop 2100 series**

2100 series laboratory turbidimeters are engineered to provide superior accuracy and sensitivity in any application. Since the first laboratory turbidimeter was introduced more than 40 years ago, the system has evolved to include advances in optics, signal processing, and software.

### Portable 2100Q series

The 2100Q portable turbidimeter offers a unique combination of advanced features, such as easy calibration and simplified data transfer, and measurement innovation for rapidly settling samples, giving you confidence you are getting accurate results every time.

### **Further information**

Find out more on our website, keywords:

- 2100N
- 2100AN
- 2100Q

Technical specifications, data sheets, manuals, recommended accessories as well as optional service contracts can be retrieved online!

### Electrochemistry

The HACH LANGE electrochemistry portfolio provides the right solution for your testing needs, backed by years of innovation and technical support. Whether you require a simple, dedicated pH meter and electrode or an advanced, expandable, multi-parameter system, HACH LANGE has your answer.



### **Accurate**

You don't just need an answer, but the right one, and fast. Optimised stabilisation algorithms in HACH LANGE meter platforms eliminate operator guesswork and reduce inaccurate measurements. Many electrodes work specifically to deliver accurate results and the quickest response, even in challenging environments.

### Easy to use

Now, more than ever, you are looking for new ways to stretch your budget and resources. Skip studying complicated user manuals and spend your time measuring, thanks to HACH LANGE innovative menu designs and simple keypads for truly intuitive operation.

### LDO - the best method to measure oxygen

Luminescence based LD0 technology is an established HACH LANGE innovation launched in 2003. The INTELLICAL LD0 is a drift-free sensor providing error-free and accurate results at high and low  $\rm O_2$  concentrations with minimum effort. There is no calibration and no replacement of electrolyte. Since it was launched, LD0 has proved itself everywhere where  $\rm O_2$  is measured!



### **HQD Meters and INTELLICAL Probes**



### Easy to use - no matter where you put them to work

HQD's user-friendly meters allow even new operators to produce accurate measurements while reducing possible errors — saving you time and hassle. Additionally, the durable design is built to handle the toughest conditions — in both the field and the lab.

- ► Automatically detect the testing parameter and calibration history
- ► Large graphic display makes the results easy to read, even in difficult light conditions
- ▶ Mix and match of electrodes and meters ensures reliability and flexibility
- ► HQD meters communicate clearly in 13 languages



#### **INTELLICAL** smart probes

INTELLICAL digital probes provide ultimate traceability of calibration history. The probes can be moved between meters without the need to recalibrate or re-enter measurement settings.

- ► Minimises errors
- ► Minimises setup time



#### **Break-through technology**

Breakthrough HACH LANGE Luminescent Dissolved Oxygen (LDO) technology, now the standard for measuring dissolved oxygen, eliminates the numerous reliability and maintenance concerns inherent in older DO membrane instruments.

- ► No membranes
- ► No electrolyte to replace or anode to polish
- ► Infrequent calibrations
- ► 1-year sensor cap life



### SENSION+ Meters and Probes



### Comes complete with everything you need to start testing

Each portable and benchtop kit comes complete with everything you need to start testing.

- ► SENSION+ Field Kits work as complete mobile measurement stations, which include meter, probe, robust carrying case and all supporting chemistries.
- ▶ Portable instruments offer real one-hand operation with ergonomic and lightweight design and IP67 protection. Screw-on calibration tubes simplify on-site calibration and minimize buffer consumption.
- ► SENSION+ Laboratory Stations include meter, probe, integrated probe stand, instrument controlled magnetic stirrer and supporting chemistries.
- ► The SENSION+ laboratory portfolio covers all application requirements from a basic pH meter up to the multi-channel GLP system for pH, ORP, Conductivity and ISE with data management and full PC driven options.



### Simple and fast measurements

A guided navigation menu allows you to follow simple prompts to set up your measurements, reducing the time necessary to start testing, and improving the accuracy of measurements.

- ► Guided menu makes navigation easy
- ► Simple prompts to ease set-up
- ► Keypads designed for truly intuitive operation
- ► Accurate results with fast response, even in challenging environments



### Available for a wide variety of applications

The portfolio of high quality probes built for SENSION+ means that you have the right option for nearly every testing environment.

- **▶** pH
- ► Conductivity
- ► Dissolved Oxygen
- ► Multi sensors for pH, ORP and Conductivity
- ▶ and much, much more!



For laboratory instruments you may also choose red rod pH electrodes or other high end probes from our large Radiometer probe portfolio - ask for Radiometer probes dedicated to SENSION+.

### **NEW: SENSION+ Portable Data Logger versions**

All-in-one systems - now making testing and data management fast and simple: interval measurements, 500 data point memory and wireless data transfer to your PC. Data management is as fast and simple as can be.



# **HQD Benchtop Meters**



All the benefits of the digital HQD system, with simplified data transfer and easy-to-read results on a large, backlit screen.

| Parameter  | HQ411D<br>DEDICATED<br>pH/Mv<br>BENCHTOP METER | HQ430D SINGLE INPUT<br>MULTI-PARAMETER<br>BENCHTOP METER WITH<br>PROBE STAND | HQ440D DUAL INPUT<br>Multi-Parameter<br>Benchtop Meter With<br>Probe Stand |
|--|--|--|--|
| Temperature  |  |  |  |
| pH<br>Glass<br>Non-glass (ISFET)                         |  |  | •  |
| mV   |  |  |  |
| Conductivity   |  |  |  |
| TDS  |  |  |  |
| Salinity   |  |  |  |
| Resistivity  |  |  |  |
| Dissolved Oxygen Luminescent (LDO) BOD Sensor (with LDO) |  | :  | :  |
| ORP/Redox  |  |  |  |
| Ammonia  |  |  |  |
| Ammonium   |  |  |  |
| Chloride   |  |  |  |
| Fluoride   |  |  |  |
| Nitrate  |  |  |  |
| Sodium   |  |  |  |

| Specifications                   | HQ411D<br>DEDICATED<br>pH/Mv<br>BENCHTOP METER | HQ430D SINGLE INPUT<br>MULTI-PARAMETER<br>BENCHTOP METER WITH<br>PROBE STAND | HQ440D DUAL INPUT<br>Multi-Parameter<br>Benchtop Meter With<br>Probe Stand |
|----------------------------------|--|--|--|
| Casing IP Rating                 | IP54   | IP54   | IP54   |
| Internal result storage capacity | 500*   | 500*   | 500*   |
| Inputs                           | M12 digital (1) for INTELLICAL probes          | M12 digital (1) for INTELLICAL probes  | M12 digital (2) for INTELLICAL probes                                      |
| Outputs                          | USB to PC / flash stick                        | USB to PC / flash stick  | USB to PC / flash stick  |
| Resolution                       | 0.1/ 0.01/ 0.001                               | 0.1/ 0.01/ 0.001   | 0.1/ 0.01/ 0.001   |
| Interface languages              | 13**   | 13**   | 13**   |
| Warranty                         | 3 years  | 3 years  | 3 years  |
| Compliance                       | CE.WEEE  | CE.WEEE  | CE.WEEE  |
| GLP features                     |  |  |  |
| PC data transfer software        | included                                       | included   | included   |
| Backlight                        |  |  |  |
| Battery requirements             | 4, AA  | 4, AA  | 4, AA  |
| AC and USB operation             |  |  |  |
| Article number                   | HQ411D.98.00002                                | HQ430D.98.00002  | HQ440D.98.00002  |

<sup>\*</sup>Expanded storage with external USB storage device. \*\*English, German, Spanish, Portuguese, French, Italian, Dutch, Turkish, Polish, Danish, Swedish, Russian, Czech

### **HQD Portable Meters**



Now offering a complete water analysis portfolio of testing parameters with standard and rugged options. HACH LANGE'S HQD system gives maximum measurement flexibility and ease of operation with interchangeable probes and automatic parameter recognition.

| Parameter  | HQ11D<br>Dedicated<br>pH/mV<br>Portable Meter | HQ14D<br>Dedicated<br>Conductivity<br>Portable Meter | HQ30D<br>Single input<br>Multi-parameter<br>Portable meter | HQ40D<br>DUAL INPUT<br>MULTI-PARAMETER<br>PORTABLE METER |
|--|---|--|--|--|
| Temperature  |   |  |  |  |
| pH<br>Glass<br>Non-glass (ISFET)                         |   |  |  | •  |
| mV   |   |  |  |  |
| Conductivity   |   |  |  |  |
| TDS  |   |  |  |  |
| Salinity   |   |  |  |  |
| Resistivity  |   |  |  |  |
| Dissolved Oxygen Luminescent (LDO) BOD Sensor (with LDO) |   |  | :  | :  |
| ORP/Redox  |   |  |  | 100  |
| Ammonia  |   |  |  |  |
| Ammonium   |   |  |  |  |
| Chloride   |   |  |  |  |
| Fluoride   |   |  |  | 100  |
| Nitrate  |   |  |  | 100  |
| Sodium   |   |  |  |  |

| Specifications                   | HQ11D<br>DEDICATED<br>pH/mV<br>PORTABLE METER | HQ14D<br>DEDIGATED<br>CONDUCTIVITY<br>PORTABLE METER | HQ30D<br>SINGLE INPUT<br>MULTI-PARAMETER<br>PORTABLE METER | HQ40D<br>DUAL INPUT<br>MULTI-PARAMETER<br>PORTABLE METER |
|----------------------------------|---|--|--|--|
| Casing IP Rating                 | IP67  | IP67   | IP67   | IP67   |
| Internal result storage capacity | 500*  | 500*   | 500*   | 500*   |
| Inputs                           | M12 digital (1) for INTELLICAL probes         | M12 digital (1) for INTELLICAL probes                | M12 digital (1) for INTELLICAL probes                      | M12 digital (2) for INTELLICAL probes                    |
| Outputs                          | USB to PC / flash stick                       | USB to PC / flash stick                              | USB to PC / flash stick                                    | USB to PC / flash stick                                  |
| Resolution                       | 0.1/ 0.01/ 0.001                              | 0.1/ 0.01/ 0.001                                     | 0.1/ 0.01/ 0.001   | 0.1/ 0.01/ 0.001   |
| Interface languages              | 13**  | 13**   | 13**   | 13**   |
| Warranty                         | 3 years                                       | 3 years  | 3 years  | 3 years  |
| Compliance                       | CE.WEEE                                       | CE.WEEE  | CE marking   | CE marking   |
| Benchtop                         | with stand                                    | with stand   | with stand   | with stand   |
| GLP features                     |   |  |  |  |
| PC data transfer software        | included                                      | included   | included   | included   |
| Backlight                        |   |  |  |  |
| Battery requirements             | 4, AA   | 4, AA  | 4, AA  | 4, AA  |
| AC and USB operation             | Optional                                      | Optional   | Optional   | Included   |
| Article number                   | HQ11D.99.000000                               | HQ14D.99.000000                                      | HQ30D.99.000000  | HQ40D.99.000000  |

<sup>\*</sup>Expanded storage with external USB storage device. \*\*English, German, Spanish, Portuguese, French, Italian, Dutch, Turkish, Polish, Danish, Swedish, Russian, Czech



# **HQD INTELLICAL Probes**



| Specifications        | DO Luminescent<br>Dissolved Oxygen                   | DO Luminescent<br>Dissolved Oxygen                           | BOD Luminescent<br>Dissolved Oxygen                          | Conductivity   | Conductivity  |
|-----------------------|--|--|--|--|---|
| Special feature       | LDO technology.<br>No membranes.                     | LDO technology.<br>No membranes.                             | LDO technology. Designed for BOD applications. No membranes. |  |   |
| Electrode type        | Laboratory Luminescent<br>Dissolved Oxygen           | Rugged Luminescent<br>Dissolved Oxygen                       | Laboratory   | Laboratory   | Rugged Outdoor  |
| Measuring range       | 0.05 - 20.00 mg/L                                    | 0.05 - 20.00 mg/L  | 0.05 - 20.00 mg/L  | Cond.: 0.01 µS/cm - 200 mS/cm<br>TDS: 0 - 50000 mg/L as NaCl<br>Salinity: 0 - 42 g/kg or<br>Resist.: 2.5 Ωcm - 49 ΜΩcm | Cond.: 0.01 $\mu$ S/cm - 200 mS/cm<br>TDS: 0 - 50000 mg/L as NaCl<br>Salinity: 0 - 42 g/kg or<br>Resist.: 2.5 $\Omega$ cm - 49 $M\Omega$ cm |
| Accuracy              | ±0.1 from 0 - 8 mg/L<br>±0.2 for greater than 8 mg/L | $\pm 0.1$ from 0 - 8 mg/L $\pm 0.1$ for greater than 10 mg/L | ±0.05 from 0 - 10 mg/L<br>±0.1 for greater than 10 mg/L      | Cond: $\pm 0.5\%$ of range TDS: $\pm 0.5\% \pm 1$ digit Salinity: $\pm 0.1$ , $\pm 1$ digit                            | Cond: $\pm 0.5\%$ of range<br>TDS: $\pm 0.5\%$ $\pm 1$ digit<br>Salinity: $\pm 0.1$ , $\pm 1$ digit   |
| Temperature range     | 0 - 50 °C  | 0 - 50 °C  | 0 - 50 °C  | -10 - 110 °C   | -10 - 110 °C  |
| Dimensions<br>(D x L) | 15 mm x 200 mm                                       | 45 mm x 250 mm   | 15.875 mm x 200 mm   | 14 mm x 200 mm   | 45 mm x 250 mm  |
| Sensor type           | Lumiphore  | Lumiphore  | Lumiphore  | 4-pole graphite, $k = 0.40 \text{ cm}^{-1}$  | 4-pole graphite, k = 0.40 cm <sup>-1</sup>  |
| Material              | Sensor Body:<br>Polycarbonate / Abs                  | Sensor Body: Polycarbonate /<br>Abs with Stainless Steel     | Sensor Body:<br>Polycarbonate / Abs                          | Sensor Body:<br>Noryl  | Sensor Body:<br>Noryl with Stainless Steel  |
| Article number        | LD010101   | LD010105   | LB0D10101  | CDC40101   | CDC40105  |



| Specifications        | pH Ultra<br>Refillable Combination | pH<br>Gel-filled Combination     | pH<br>Gel-filled Combination             | pH<br>Refillable Combination | pH<br>Gel-filled Combination                |
|-----------------------|------------------------------------|----------------------------------|--|------------------------------|---|
| Special feature       | Fast response time                 | Wastewater and difficult samples | General purpose.<br>Clean water samples. |                              | Low maintenance                             |
| Electrode type        | Laboratory                         | Laboratory                       | Laboratory                               | Laboratory                   | Rugged Outdoor                              |
| Measuring range       | 0 - 14 pH                          | 2 - 14 pH                        | 0 - 14 pH                                | 0 - 14 pH                    | 2 - 14 pH                                   |
| Accuracy              | ±0.02 pH                           | ±0.02 pH                         | ±0.02 pH                                 | ±0.02 pH                     | ±0.02 pH                                    |
| Temperature range     | 0 - 80 °C                          | 0 - 50 °C                        | 0 - 80 °C                                | 0 - 80 °C                    | 0 - 50 °C                                   |
| Dimensions<br>(D x L) | 12 mm x 200 mm                     | 12 mm x 200 mm                   | 12 mm x 200 mm                           | 12 mm x 200 mm               | 45 mm x 250 mm                              |
| Sensor type           | Glass                              | Glass                            | Glass                                    | Glass                        | Glass                                       |
| Reference             | Ag/AgCl (double junction)          | Ag/AgCl (double junction)        | Ag/AgCl (double junction)                | Ag/AgCl (double junction)    | Ag/AgCl (double junction)                   |
| Electrode junction    | Open                               | Open                             | Ceramic pin                              | Ceramic pin (x2)             | Open  |
| Material              | Sensor Body:<br>Zeonor             | Sensor Body:<br>Epoxy            | Sensor Body:<br>Epoxy                    | Sensor Body:<br>Zeonor       | Sensor Body:<br>Zeonor with Stainless Steel |
| Filling solution      | 2965026                            | Non-refillable solid gel         | Non-refillable gel                       | 2841700                      | Non-refillable solid gel                    |
| Article number        | PHC28101                           | PHC10101                         | PHC20101                                 | PHC30101                     | PHC10105                                    |

# **HQD INTELLICAL Probes**











|                       | ORP/Redox  | ORP/Redox                              | ORP/Redox                                   | Ammonia  | Ammonium  |
|-----------------------|--|--|---|--|---|
| Specifications        | Combination Gel-filled                               | Combination Refillable                 | Combination Gel-filled                      | Combination ISE  | Combination ISE   |
| Special feature       | Flat disc sensor for easy cleaning. Low maintenance. | Flat disc sensor for easy cleaning.    | Flat disc sensor for easy cleaning.         | Easy-to-replace<br>membrane modules.                                       | Dry storage & fast<br>response time.<br>No replacement membranes.           |
| Electrode type        | Laboratory   | Laboratory                             | Rugged Outdoor                              | Laboratory Combination ISE   | Laboratory Combination ISE  |
| Measuring range       | ± 1200 mV  | ± 1200 mV                              | ± 1200 mV                                   | 0.01 mg/L (5x10 <sup>7</sup> M) -<br>14,000 mg/L (1 M) NH <sub>3</sub> -N  | 0.018 mg/L (10 <sup>-6</sup> M) -<br>9,000 mg/L (0.5 M) NH <sub>4</sub> +-N |
| Accuracy              | ±0.02mV or 0.05%, whichever is greater               | ±0.02mV or 0.05%, whichever is greater | ±0.02mV or 0.05%, whichever is greater      | ±0.02mV or 0.05%, whichever is greater                                     | ±0.02mV or 0.05%, whichever is greater                                      |
| Temperature range     | 0 - 80 °C  | 0 - 80 °C                              | 0 - 80 °C                                   | 0 - 50 °C  | 0 - 50 °C   |
| Dimensions<br>(D x L) | 12 mm x 200 mm                                       | 12 mm x 200 mm                         | 45 mm x 250 mm                              | 12 mm x 220 mm   | 12 mm x 220 mm  |
| Sensor type           | Platinum disc  | Platinum disc                          | Platinum disc                               | Gas sensing:<br>glass w/ replaceable NH <sub>3</sub><br>sensitive membrane | Solid-state<br>PVC membrane   |
| Reference             | Ag/AgCl  | Ag/AgCI                                | Ag/AgCI                                     | Ag/AgCl  | Ag/AgCl   |
| Electrode junction    | Open   | Ceramic pin                            | Open  | Porous Teflon<br>Annular Ring  | Porous Teflon<br>Annular Ring   |
| Material              | Sensor Body:<br>Epoxy                                | Sensor Body:<br>Epoxy                  | Sensor Body:<br>Zeonor with Stainless Steel | Sensor Body:<br>Epoxy  | Sensor Body:<br>Epoxy   |
| Filling solution      | Non-refillable gel                                   | 2841700                                | Non-refillable gel                          | 4447226  | Non-refillable Dritek gel   |
| ISA required          |  |  |   | 4447169  | 2980699   |
| Article number        | MTC10101   | MTC30101                               | MTC10105                                    | ISENH318101  | ISENH418101   |









| Specifications        | Chloride<br>Combination ISE  | Fluoride<br>Combination ISE  | Nitrate<br>Combination ISE   | Sodium<br>Combination ISE  |
|-----------------------|--|--|--|--|
| Special feature       | Dry storage & fast response time.<br>No replacement membranes.         | Dry storage & fast response time.  No replacement membranes.           | Dry storage & fast response time. No replacement membranes.                            | Fast response time   |
| Electrode type        | Laboratory Combination ISE   | Laboratory Combination ISE   | Laboratory Combination ISE   | Laboratory Combination ISE   |
| Measuring range       | 0.1 mg/L (3x10 <sup>-6</sup> M) -<br>35,500 mg/L (1 M) Cl <sup>-</sup> | 0.01 mg/L (5x10 <sup>-7</sup> M) -<br>19,000 mg/L (1 M) F <sup>-</sup> | 0.1 mg/L (7x10 <sup>-6</sup> M) -<br>14,000 mg/L (1 M) NO <sub>3</sub> <sup>-</sup> -N | 0.023 mg/L (1x10 <sup>-6</sup> M) -<br>23,000 mg/L (1 M) Na <sup>+</sup> |
| Accuracy              | ±0.02mV or 0.05%, whichever is greater                                 | ±0.02mV or 0.05%, whichever is greater                                 | ±0.02mV or 0.05%,<br>whichever is greater  | ±0.02mV or 0.05%, whichever is greater                                   |
| Temperature range     | 5 - 50 °C  | 5 - 50 °C  | 0 - 50°C   | 0 - 50 °C  |
| Dimensions<br>(D x L) | 12 mm x 220 mm   | 12 mm x 220 mm   | 12 mm x 220 mm   | 12 mm x 220 mm   |
| Sensor type           | Solid-state<br>crystal membrane  | Solid-state<br>crystal membrane  | Solid-state<br>PVC membrane  | Glass  |
| Reference             | Ag/AgCl  | Ag/AgCl  | Ag/AgCl  | Ag/AgCl  |
| Material              | Sensor Body: Epoxy   | Sensor Body: Epoxy   | Sensor Body: Epoxy   | Sensor Body: Zeonor  |
| Filling solution      | Non-refillable Dritek gel  | Non-refillable Dritek gel  | Non-refillable Dritek gel  | 2965126  |
| ISA required          | 2318069  | 258999   | 2984799  | 4451569  |
| Article number        | ISECL18101   | ISEF12101  | ISEN0318101  | ISENA38101   |



### **HQD Meter & Probe Accessories**

### Meter Stand



Article number Description

4754900 Perfect for use in the laboratory or other applications where

hands-free meter operation is desired. Black molded plastic

is durable and easy to clean.

### **Universal Probe Stand** for Standard Electrode Sizes



Article number

Description

8508850

Also supports HACH LANGE INTELLICAL and SENSION+ probes.

### Protective Glove



Article number

Description

5828700

Provides added impact protection that field use requires. Wrist and neck straps keep the meter secure. Two probe holders (not included in kit) can slide onto glove.

### **HQD Field Meter Protective Caps**



Article number

9345200

Further protect your HQD meter from harsh (corrosive) field environments. For use with HQD Meter Protective Glove kit (5828700).

### **Probe Holder** (Standard INTELLICAL Probes Only)



Article number

Description

5829400 For use with the Protective Glove. Simply wrap 1- or 3-meter

cables around the holder and slide probe into protective sheath. Protective Glove can hold one or two Probe Holders.

### Field Kit



Article number

5825800

Includes protective glove and (5) 120 mL sample cups. Carry your HQD system anywhere. Custom inserts organise and securely position your meter and probes. NOTE: probes and

meter not included.

### Rugged Field Case for Two Probes with 5 m Cables



Article number

Description

8505500

Includes: Empty case and insert for meter and probe storage; (4) Containers for sample collection; 500 mL Wash Bottle.

### **Rugged Field Case for Three Probes** with 5 m Cables



Article number

8505501

Includes same accessories as case 8505500.

### **HQD Meter & Probe Accessories**

### **Colour Coded Probe Clips**



Article number

Description

5818400

Clip a colour-coded band to each end of your probe cable for identification. Ten clips of five different colors per package. NOTE: One package of clips included with each probe.

### Probe Cable Depth Markers



Article number

Description

5828610

Taking readings at specific depths with rugged probes is a snap with these Depth Markers. Visually identify the depth of your probe by attaching Depth Marker securely on cable at points to meet your needs. Five Depth Markers per package. NOTE: For use with rugged probes only.

### **Replacement Stirrer Assembly for INTELLICAL LBOD Probe**



Article number 5850800

### Replacement Shroud Kit for Rugged Probes



Article number

Description

5825900

Includes protective bell and locking ring. The rugged shroud absorbs the impact from rough wear and tear. NOTE: probe not included.

### Replacement LDO Sensor Cap



Article number

5811200

Description

Includes i-button and sensor cap.

### Replacement LBOD Sensor Cap



Article number 5838000

### Replacement Membranes for INTELLICAL ISENH3181 Ammonia ISE



Article number

Description

5812711

Pack of 3 replacement membrane modules designed for the INTELLICAL ISENH3181 Ammonia ISE.

### HQD Meter USB and AC Power Adapter Kit



Article number

Description

5834100

Provides A/C power in the lab. The USB adapter can be used to transfer data from the meter to a PC, memory stick, or printer. Included with HQ40D meter systems.







### **SENSION+ Portable Meters**



SENSION+ Portable instruments offer real one-hand operation with ergonomic and lightweight design and IP67 protection. Screw-on calibration tubes simplify on-site calibration and minimize buffer consumption.

| Parameter        | PH1<br>Basic<br>Portable<br>PH Meter | MM110<br>Portable<br>pH/ORP<br>Multi-Sensor Meter' | EC5<br>Portable<br>Conductivity<br>Meter | DOG<br>Portable<br>Dissolved<br>Oxygen Meter | MM150<br>Portable<br>Multi-Meter² |
|------------------|--------------------------------------|--|--|--|-----------------------------------|
| Temperature      |                                      |  |  |  |                                   |
| pH               |                                      |  |  |  |                                   |
| Conductivity     |                                      |  |  |  |                                   |
| Salinity         |                                      |  |  |  |                                   |
| TDS              |                                      |  |  |  |                                   |
| Dissolved Oxygen |                                      |  |  |  |                                   |
| ORP/Redox        |                                      |  |  |  |                                   |

| Specifications                   | PH1<br>Basic<br>Portable<br>PH Meter | MM110<br>Portable<br>pH/ORP<br>Multi-Sensor Meter' | EC5<br>Portable<br>Conductivity<br>Meter  | DO6<br>Portable<br>Dissolved<br>Oxygen meter | MM150<br>Portable<br>Multi-Meter <sup>2</sup>   |
|----------------------------------|--------------------------------------|--|---|--|---|
| Casing IP Rating                 | IP67                                 | IP67   | IP67  | IP67   | IP67  |
| Internal result storage capacity | _                                    | _  | _   | _  | _   |
| Resolution                       | 0.01 pH<br>0.1 mV (±199.9 mV range)  | 0.01 pH<br>1 mV                                    | EC: 0.01 µS/cm - 1 mS/cm<br>TDS: 1mg/L - 1g/L<br>depending on range<br>Salinity: 0.1 mg/L - 0.1 g/L<br>depending on range | 0.01 mg/L                                    | pH: 0.01<br>ORP: 1 mV<br>EC: 0.01 μS/cm - 1 mS/cm<br>depending on range<br>1 mg/L - 1 g/L<br>depending on range |
| Inputs                           | MP5                                  | MP8  | MP5   | MP5  | MP8   |
| Output                           | _                                    | _  | _   | _  | _   |
| Interface languages              | Icon based,<br>no language needed    | Icon based,<br>no language needed                  | lcon based,<br>no language needed   | Icon based,<br>no language needed            | lcon based,<br>no language needed   |
| Warranty                         | 2 years                              | 2 years  | 2 years   | 2 years                                      | 2 years   |
| Compliance                       | CE.WEEE                              | CE.WEEE  | CE.WEEE   | CE.WEEE                                      | CE.WEEE   |
| Measurement method               | Auto stabilisation, manual           | Auto stabilisation, manual                         | Auto stabilisation, manual  | Auto stabilisation, manual                   | Auto stabilisation, manual  |
| PC data transfer software        | _                                    | _  | _   | _  | _   |
| Backlight                        |                                      |  |   |  |   |
| Battery requirements             | 3, AA                                | 3, AA  | 3, AA   | 3, AA  | 3, AA   |
| Article number                   | LPV2500.98.0002                      | LPV2600.98.0002                                    | LPV3500.98.0002   | LPV4500.98.0002                              | LPV4000.98.0002   |

¹pH and ORP/Redox simultaneously with one multi electrode. ² Parameter measured depends on electrode selected.

# SENSION+ Portable Data Logger Meters



HACH LANGE'S SENSION+ Data Logger, an all-in-one system - now making testing and data management fast and simple. Interval measurements, 500 data point memory and wireless data transfer to your PC. Data management as fast and simple as it can be.

| Parameter        | pH1 DL PORTABLE<br>Data logger<br>ph meter | MM110 DL PORTABLE<br>Data logger<br>Multi-Meter¹ | EC5 DL<br>Portable<br>Data logger<br>Conductivity meter | DOG DL<br>PORTABLE<br>DATA LOGGER DISSOLVED<br>OXYGEN METER | MM150 DL<br>Portable data logger<br>Multi-Meter² |
|------------------|--|--|---|---|--|
| Temperature      |  |  |   |   |  |
| рН               |  |  |   |   |  |
| Conductivity     |  |  |   |   |  |
| Salinity         |  |  |   |   |  |
| TDS              |  |  |   |   |  |
| Dissolved Oxygen |  |  |   |   |  |
| ORP/Redox        |  |  |   |   |  |

| Specifications                   | pH1 DL PORTABLE<br>Data logger<br>pH Meter                      | MM110 DL PORTABLE<br>Data logger<br>Multi-Meter'                | EC5 DL<br>Portable<br>Data logger<br>Conductivity meter   | DO6 DL<br>Portable<br>Data logger dissolved<br>Oxygen meter | MM150 DL<br>Portable data logger<br>Multi-Meter²  |
|----------------------------------|---|---|---|---|---|
| Casing IP Rating                 | IP67  | IP67  | IP67  | IP67  | IP67  |
| Internal result storage capacity | 500 results   | 500 results   | 500 results   | 500 results   | 500 results   |
| Resolution                       | 0.01 pH<br>0.1 mV (±199.9 mV range)                             | 0.01 pH<br>1 mV   | EC: 0.01 µS/cm - 1 mS/cm 0.01 mg/L TDS: 1mg/L - 1g/L depending on range Salinity: 0.1 mg/L - 0.1 g/L depending on range |   | pH: 0.01<br>ORP: 1 mV<br>EC: 0.01 μS/cm - 1 mS/cm<br>depending on range<br>1 mg/L - 1 g/L<br>depending on range |
| Inputs                           | MP-5  | MP-8  | MP-5  | MP-5  | MP-8  |
| Output                           | Wireless to USB   | Wireless to USB   | Wireless to USB   | Wireless to USB   | Wireless to USB   |
| Interface languages              | lcon based,<br>no language needed                               | Icon based,<br>no language needed                               | Icon based,<br>no language needed   | lcon based,<br>no language needed                           | Icon based,<br>no language needed   |
| Warranty                         | 2 years   | 2 years   | 2 years   | 2 years   | 2 years   |
| Compliance                       | CE; RTTE directive see<br>NB opinion statement;<br>FCC approval | CE; RTTE directive see<br>NB opinion statement;<br>FCC approval | CE; RTTE directive see NB opinion statement; FCC approval   | CE; RTTE directive see NB opinion statement; FCC approval   | CE; RTTE directive see NB opinion statement; FCC approval   |
| Measurement method               | Interval, auto stabilisation,<br>manual                         | Interval, auto stabilisation,<br>manual                         | Interval, auto stabilisation, manual  | Interval, auto stabilisation,<br>manual                     | Interval, auto stabilisation, manual  |
| PC data transfer software        | Yes   | Yes   | Yes   | Yes   | Yes   |
| Backlight                        |   |   |   |   |   |
| Battery requirements             | 3, AA   | 3, AA   | 3, AA   | 3, AA   | 3, AA   |
| Article number                   | LPV2500DL.98.02   | LPV2600DL.98.02   | LPV3500DL.98.02   | LPV4500DL.98.02   | LPV4000DL.98.02   |

 ${}^1\!P\!H\ and\ ORP/Redox\ simultaneously\ with\ one\ multi\ electrode.\ {}^2\!P\!arameter\ measured\ depends\ on\ electrode\ selected.$ 



SENSION+ Benchtop Meters



SENSION+ benchtop meters are developed as complete work stations, delivered with integrated magnetic stirrer, probe holder and consumables.

| Parameter    | PH3<br>Basic<br>Benchtop<br>ph Meter | PH31<br>GLP<br>Benchtop<br>PH Meter | EC7<br>Benchtop<br>Conductivity<br>Meter | EC71<br>GLP<br>Benchtop<br>Conductivity<br>Meter | MM340<br>BENCHTOP<br>MULTI-<br>Parameter<br>Meter | MM374<br>BENCHTOP<br>MULTI-<br>PARAMETER<br>METER |
|--------------|--------------------------------------|-------------------------------------|--|--|---|---|
| Temperature  |                                      |                                     |  |  |   |   |
| рН           |                                      |                                     |  |  |   |   |
| mV           |                                      |                                     |  |  |   |   |
| Conductivity |                                      |                                     |  |  |   |   |
| Salinity     |                                      |                                     |  |  |   |   |
| TDS          |                                      |                                     |  |  |   |   |
| ORP/Redox    |                                      |                                     |  |  |   |   |
| Ammonia      |                                      |                                     |  |  |   |   |
| Nitrate      |                                      |                                     |  |  |   |   |
| Fluoride     |                                      |                                     |  |  |   |   |
| Sodium       |                                      |                                     |  |  |   |   |
| Chloride     |                                      |                                     |  |  |   |   |
| Ammonium     |                                      |                                     |  |  |   |   |

| Specifications                   | PH3<br>BASIC<br>BENCHTOP<br>pH METER | PH31<br>GLP<br>BENCHTOP<br>pH METER | EC7<br>BENCHTOP<br>CONDUCTIVITY<br>METER                                 | EC71<br>GLP<br>BENCHTOP<br>CONDUCTIVITY<br>METER  | MM340<br>BENCHTOP<br>MULTI-<br>PARAMETER<br>METER | MM374<br>BENCHTOP<br>MULTI-<br>PARAMETER<br>METER  |
|----------------------------------|--------------------------------------|-------------------------------------|--|---|---|--|
| Internal result storage capacity | — PH WEIEN                           | 330                                 |  | 400   | 330   | 330  |
| Resolution                       | pH: 0.01<br>ORP: 1 mV                | pH: 0.001<br>ORP: 0.1 mV            | EC: 0.01 µS/cm - 1 mS/cm depending on range Salinity: Dependent on range | EC: 0.01 µS/cm - 1 mS/<br>cm depending on range<br>TDS: 1 mg/L - 1 g/L<br>Salinity:<br>Dependent on range | pH: 0.001<br>ORP 0.1 mV                           | pH: 0.001 pH ORP: 0.1 mV EC: 0.01 µS/cm - 1 mS/ cm depending on range TDS: 1 mg/L - 1 g/L Salinity: Dependent on range |
| Inputs                           | BNC, ref, Pt1000                     | BNC, ref,<br>Pt1000                 | E.C. cell,<br>Pt1000   | E.C. cell,<br>Pt1000  | 2 BNC,<br>2 ref, Pt1000                           | 2 BNC,<br>2 ref, 1 E.C.  |
| Outputs                          |                                      | RS232*<br>(two-way), USB            |  | RS-232*<br>(two-way), USB   | RS-232*<br>(two-way), USB                         | RS-232*<br>(two-way), USB  |
| Interface languages              | 6**                                  | 6**                                 | 6**  | 6**   | 6**   | 6**  |
| Warranty                         | 2 years                              | 2 years                             | 2 years  | 2 years   | 2 years   | 2 years  |
| Compliance                       | CE.WEEE                              | CE.WEEE                             | CE.WEEE  | CE.WEEE   | CE.WEEE   | CE.WEEE  |
| GLP features                     | _                                    |                                     | _  |   |   |  |
| PC data transfer software        |                                      |                                     | •  | •   | •   |  |
| Backlight                        |                                      |                                     |  |   |   |  |
| AC and USB operation             |                                      |                                     |  |   |   |  |
| Article number                   | LPV2000.98.0002                      | LPV2100.98.0002                     | LPV3010.98.0002  | LPV3110.98.0002   | LPV2200.98.0002                                   | LPV4110.98.0002  |

 $<sup>{}^\</sup>star \text{USB compatibility with optional adapter.storage device.} \\ {}^{\star\star} \text{English, German, Spanish, French, Italian, and Portuguese} \\$ 

## **SENSION+ Probes for Portable Meters**



| Specifications        | pH Combination<br>Gel-filled  | pH Combination<br>Solid Gel<br>Wastewater                     | pH Combination<br>Low Conductivity<br>High Temperature                  | ORP/Redox<br>Combination<br>Gel-filled | Dissolved Oxygen<br>Polarographic                  |
|-----------------------|-------------------------------|---|---|--|--|
| Special feature       |                               | With Working Protector.<br>Solid gel for high solids content. | With Working Protector. For low ionic strength and/or high temperature. |  |  |
| Measuring range       | 0 - 14 pH                     | 2 - 14 pH   | 0 - 14 pH   | ± 2000 mV                              | 0.03 mg/L to saturation                            |
| Accuracy              | 0.02 pH with pH1 Meter        | 0.02 pH with pH1 Meter  | 0.02 pH with pH1 Meter  | 1 mV with pH1 Meter                    | ±0.5 mg/L  |
| Temperature range     | 0 - 80 °C                     | 0 - 80 °C   | 0 - 100 °C  | 0 - 80 °C                              | 0 - 50 °C  |
| Thermistor            | Pt1000                        | Pt1000  | Pt1000  | -                                      | Integrated 30kohm NTC                              |
| Dimensions<br>(D x L) | 12 mm x 85 mm                 | 12 mm x 85 mm   | 12 mm x 85 mm   | 12 mm x 85 mm                          | 12 mm x 120 mm                                     |
| Sensor type           | Glass                         | Glass   | Glass   | Platinum annular ring                  | Replaceable PTFE membrane;<br>Pt cathode; Ag anode |
| Electrode junction    | Ceramic pin                   | Open  | Porous annular PTFE   | Ceramic pin                            |  |
| Material              | Sensor Body:<br>Polycarbonate | Sensor Body:<br>Glass   | Sensor Body:<br>Glass   | Sensor Body:<br>Polycarbonate          | Probe body:<br>ABS & Delrin (nylon)                |
| Filling solution      | Non-refillable gel            | Non-refillable solid polymer                                  | Non-refillable gel  | Non-refillable gel                     | 2759123  |
| Article number        | LZW5050T.97.002               | LZW5051T.97.002   | LZW5052T.97.002   | LZW5055.97.0002                        | LZW5130.97.0002                                    |



| Specifications        | Conductivity  | Conductivity                               | Multi-<br>Combination**  | Multi-<br>Combination*                           | Multi-<br>Combination**                               |
|-----------------------|---|--|--|--|---|
| Special feature       |   | For harsh samples                          | pH, Conductivity, ORP  | pH, ORP  | pH, Conductivity                                      |
| Measuring range       | 0.2 μS/cm - 200 mS/cm                                 | 5 μS/cm - 50 mS/cm                         | pH: 0 - 14 pH  | pH: 0 - 14 pH                                    | pH: 0 - 14 pH   |
|                       |   |  | Conductivity:<br>20 µS/cm - 200 mS/cm  |  | Conductivity:<br>20 µS/cm - 200 mS/cm                 |
|                       |   |  | ORP:± 2000 mV  | ORP: ± 2000 mV                                   |   |
| Accuracy              | 0.50% with EC5 Meter                                  | 0.50% with EC5 Meter                       | 0.02 pH with MM150 Meter<br>1 mV with MM150 Meter<br>0.50% EC with MM150 Meter | 0.02 pH with MM110 Meter<br>1mV with MM110 Meter | 0.02 pH with MM150 Meter<br>0.50% EC with MM150 Meter |
| Temperature range     | 0 - 80 °C   | 0 - 80 °C                                  | 0 - 80 °C  | 0 - 80 °C  | 0 - 80 °C   |
| Thermistor            | Pt1000  | Pt1000                                     | Pt1000   | Pt1000   | Pt1000  |
| Dimensions<br>(D x L) | 12 mm x 85 mm   | 12 mm x 85 mm                              | 12 mm x 85 mm  | 12 mm x 85 mm                                    | 12 mm x 85 mm   |
| Sensor type           | 3-pole Platinum, $k = 1.0 \text{ cm}^{-1}$            | 2-pole Titanium, $k = 0.3 \text{ cm}^{-1}$ | pH: Glass<br>ORP: Platinum<br>Conductivity: Platinum                           | pH: Glass<br>ORP: Platinum                       | pH: Glass<br>Conductivity: Platinum                   |
| Electrode junction    |   |  | Ceramic  | Ceramic  | Ceramic   |
| Material              | Sensor Body: Outside:<br>Polycarbonate; Inside: Glass | Sensor Body: Titanium                      | Sensor Body: Polycarbonate   | Sensor Body: Polycarbonate                       | Sensor Body: Polycarbonate                            |
| Filling solution      |   |  | Non-refillable gel   | Non-refillable gel                               | Non-refillable gel                                    |
| Article number        | LZW5060.97.0002                                       | LZW5062.97.0002                            | LZW5048.97.0002  | LZW5045.97.0002                                  | LZW5059.97.0002                                       |

\*Only for use on MM110 Meter. \*\*Only for use with MM150 Meter.



# SENSION+ Probes for Benchtop Meters



| Specifications        | pH Combination<br>Gel-filled   | pH Combination<br>Refillable   | pH Combination Gel-filled<br>Wastewater  | pH Combination<br>Refillable   |
|-----------------------|--|--|--|--|
| Special feature       |  | TRIS compatible  | Solid gel for high solids content  | Clog-free Sleeve Junction for low ionic strength.  |
| Measuring range       | 0 - 14 pH  | 0 - 14 pH  | 2 - 14 pH  | 0 - 14 pH  |
| Accuracy              | 0.02 pH with pH3 Meter<br>0.002 pH with pH31 Meter<br>0.002 pH with MM340 Meter<br>0.002 pH with MM374 Meter | 0.02 pH with pH3 Meter<br>0.002 pH with pH31 Meter<br>0.002 pH with MM340 Meter<br>0.002 pH with MM374 Meter | 0.02 pH with pH3 Meter<br>0.002 pH with pH31 Meter<br>0.002 pH with MM340 Meter<br>0.002 pH with MM374 Meter | 0.02 pH with pH3 Meter<br>0.002 pH with pH31 Meter<br>0.002 pH with MM340 Meter<br>0.002 pH with MM374 Meter |
| Temperature range     | 0 - 80 °C  | -10 - 100 °C   | 0 - 80 °C  | 0 - 60 °C  |
| Thermistor            | Pt1000   | Pt1000   | Pt1000   | Pt1000   |
| Dimensions<br>(D x L) | 12 mm x 130 mm   |
| Sensor type           | Glass  | Glass  | Glass  | Glass  |
| Electrode junction    | Ceramic pin  | 2 x Ceramic pin  | Open   | Open with sleeve   |
| Material              | Sensor Body: Polycarbonate   | Sensor Body: Glass   | Sensor Body: Glass   | Sensor Body: Glass   |
| Filling solution      | Non-refillable gel   | LZW9500.99   | Non-refillable solid polymer   | LZW9500.99   |
| Prod. No.             | LZW5010T.97.002  | LZW5014T.97.002  | LZW5011T.97.002  | LZW5021T.97.002  |







| Specifications        | Conductivity                                       | ORP<br>Gel-filled Combination   | ORP Refillable Combination  |
|-----------------------|--|---|---|
| Special feature       |  |   |   |
| Measuring range       | 0.2 μS/cm - 200 mS/cm                              | ± 2000 mV   | ± 2000 mV   |
| Accuracy              | 0.50% EC & TDS                                     | 1mV with pH3 Meter<br>0.2mV with pH31 Meter<br>0.2mV with MM340 Meter<br>0.2mV with MM374 Meter | 1mV with pH3 Meter<br>0.2mV with pH31 Meter<br>0.2mV with MM340 Meter<br>0.2mV with MM374 Meter |
| Temperature range     | 0 - 80 °C  | 0 - 80 °C   | 0 - 80 °C   |
| Thermistor            | Pt1000   | _   | _   |
| Dimensions<br>(D x L) | 12 mm x 130 mm                                     | 12 mm x 130 mm  | 12 mm x 130 mm  |
| Sensor type           | 3-pole Platinum, $k = 0.7 \text{ cm}^{-1}$         | Platinum  | Platinum  |
| Electrode junction    | -  | Ceramic pin   | Ceramic pin   |
| Material              | Sensor Body: Outside: Polycarbonate; Inside: Glass | Sensor Body: Glass  | Sensor Body: Glass  |
| Filling solution      | _  | Non-refillable gel  | LZW9500.99  |
| Prod. No.             | LZW5070.97.0002                                    | LZW5056.97.0002   | LZW5057.97.0002   |

# SENSION+ ISE Probes



| Specifications        | Sodium ISE  | Chloride Combination ISE   | Fluoride Combination ISE   | Nitrate Combination ISE  |
|-----------------------|---|--|--|--|
| Special feature       | Indicator only; requires a reference electrode.                 | Combination electrode<br>Maintenance free: No electrolyte<br>or spare membranes needed | Combination electrode<br>Maintenance free: No electrolyte<br>or spare membranes needed | Combination electrode<br>Maintenance free: No electrolyte<br>or spare membranes needed |
| Measuring range       | 0.05 - 23,000 mg/L Na   | 0.1 mg/L - 35,500 mg/L Cl  | 0.01 mg/L - 19,000 mg/L F  | 0.1 mg/L - 14,000 mg/L NO <sub>3</sub>   |
| Accuracy              | ±0.02 mV or 0.05%, whichever is greater (application dependent) | ±0.02 mV or 0.05%, whichever is greater (application dependent)                        | ±0.02 mV or 0.05%, whichever is greater (application dependent)                        | ±0.02 mV or 0.05%, whichever is greater (application dependent)                        |
| Temperature range     | 0 - 60 °C   | 5 - 50 °C  | 5 - 50 °C  | 0 - 50 °C  |
| Dimensions<br>(D x L) | 12 mm x 120 mm  | 12 mm x 120 mm   | 12 mm x 120 mm   | 12 mm x 120 mm   |
| Sensor type           | Sodium Selective Glass Membrane                                 | Solid-state crystal membrane   | Solid-state crystal membrane   | Solid-state PVC membrane   |
| Electrode junction    | _   | Porous Teflon Annular Ring   | Porous Teflon Annular Ring   | Porous Teflon Annular Ring   |
| Material              | Sensor Body: Glass  | Sensor Body: Epoxy   | Sensor Body: Epoxy   | Sensor Body: Epoxy   |
| Filling solution      | _   | Non-refillable Dritek gel  | Non-refillable Dritek gel  | Non-refillable Dritek gel  |
| ISA required          | 4451569   | 2318069  | 258999   | 2984799  |
| Prod. No.             | LZW9650.97.0002   | LZW9652C.97.002  | LZW9655C.97.002  | LZW9662C.97.002  |







| Specifications        | Ammonium<br>Combination ISE  | Ammonia<br>Combination ISE   | Reference<br>Electrode   |
|-----------------------|--|--|--|
| Special feature       | Combination electrode. Maintenance free:<br>No electrolyte or spare membranes needed | Gas Sensing Electrode  | Double junction reference electrode for ISEs,<br>large ceramic diaphragm |
| Measuring range       | 0.018 mg/L - 9,000 mg/L NH <sub>4</sub>  | 0.06 mg/L - 17,000 mg/L NH <sub>3</sub>                                |  |
| Accuracy              | ±0.02 mV or 0.05%, whichever is greater (application dependent)                      | ±0.02 mV or 0.05%, whichever is greater (application dependent)        |  |
| Temperature range     | 5 - 50 °C  | 0 - 50 °C  | 0 - 60 °C  |
| Dimensions<br>(D x L) | 12 mm x 120 mm   | 12 mm x 149 mm   | 12 mm x 120 mm   |
| Sensor type           | Solid-state PVC membrane   | Gas sensing: glass with replaceable NH <sub>3</sub> sensitive membrane |  |
| Electrode junction    | Porous Teflon Annular Ring   | Porous Teflon Annular Ring   | Ceramic  |
| Material              | Sensor Body: Epoxy   | Sensor Body: ABS   | Sensor Body: Glass   |
| Filling solution      | Non-refillable Dritek gel  | 4447226  | LZW9901.00   |
| ISA required          | 2980699  | 4447169  |  |
| Prod. No.             | LZW9663C.97.002  | 5192700  | LZW5044.97.0002  |



# SENSION+ Accessories



LZW2598.99.0002

| Article number   | Description   |  |
|--|---|--|
| Accessories for SENSION+ Portable Instruments & Electrodes |   |  |
| LZW9137.98   | 3 x 10 mL printed tubes for portable pH calibration                                 |  |
| LZW9161.99   | Electrode Storage Protector, Polypropylene  |  |
| LZW9162.99   | Electrode Measurement Protector, Polypropylene                                      |  |
| LZW5123.99   | Protector-calibration flask for DO probe LZW5130                                    |  |
| 5196800  | Service Kit for 5130 DO Probe (contains 2 Membrane Modules and DO Filling Solution) |  |
| LZW2598.99.0002  | PortCom Kit for SENSION+ Data Logger (USB dongle and CD-ROM with PortCom software)  |  |



| 5192711 |
|---------|

|            | 1 |
|------------|---|
| L7W9325.99 |   |

| Accessories for SE | Accessories for SENSION+ Benchtop Instruments                       |  |  |
|--------------------|---|--|--|
| LZW8997.99         | LabCom Easy PC Software for SENSION+ GLP instruments                |  |  |
| LZW8999.99         | LabCom PC Software for SENSION+ GLP instruments                     |  |  |
| LZW9008.99         | Power Supply for SENSION+ Benchtop instruments, 230-115VAC          |  |  |
| LZW9110.98         | 3x50 mL printed flasks for pH calibration, Benchtop instruments     |  |  |
| LZW9111.99         | 3x50 mL printed flasks for benchtop conductivity calibration        |  |  |
| LZW9118.99         | Pyrex glass chamber, continuous flow measurements                   |  |  |
| LZW9321.99         | Three-probe holder for SENSION+ Benchtop instruments                |  |  |
| LZW9325.99         | Radiometer probe holder for SENSION+ Benchtop instruments           |  |  |
| LZW9325.99.T014    | Radiometer probe holder kit for SENSION+ with X31T014 adapter       |  |  |
| LZW9325.99.T031    | Radiometer probe holder kit for SENSION+ with X31T031 adapter       |  |  |
| LZW9319.99         | Second Magnetic stirrer with probe holder for SENSION+ Multi Meters |  |  |
| 5192711            | Ammonia Membrane Replacement Kit, 5 pieces                          |  |  |

| Standards  |   |
|------------|---|
| LZW9463.99 | pH buffer solution 4.01, 250 mL               |
| LZW9464.97 | pH buffer solution 7.00, 250 mL               |
| LZW9471.99 | pH buffer solution 10.00, 250 mL              |
| LZW9700.99 | Conductivity Standard 147 µS/cm, 250 mL       |
| LZW9710.99 | Conductivity Standard 1413 µS/cm, 250 mL      |
| LZW9720.99 | Conductivity Standard 12.88 mS/cm, 250 mL     |
| LZW9500.99 | Electrolytic solution, KCl 3M, 250 mL (GHS07) |

## **POCKET PRO Testers**



Take the guesswork out of your measurements. The large LCD display, intuitive user interface, and standard AAA batteries make POCKET PRO the easiest to use tester for your application. POCKET PRO's superior probe and calibration diagnostics deliver confidence in your results.

### NEW

| Specifications              | POCKET PRO pH | POCKET PRO ORP  | POCKET PRO TDS <sub>LR</sub> | POCKET PRO TDS <sub>HR</sub> |
|-----------------------------|---------------|-----------------|------------------------------|------------------------------|
| Parameter                   | pH, Temp      | ORP, Temp       | TDS, Temp                    | TDS, Temp                    |
| Operating temperature range | 0 - 50 °C     | 0 - 50 °C       | 0 - 50 °C                    | 0 - 50 °C                    |
| Range                       | 0.0 - 14.0 pH | -999 to +999 mV | 0 - 1999 ppm                 | 0 - 10.00 ppt                |
| Accuracy                    | 0.1 pH        | ± 2 mV          | 1% FS                        | 2% FS                        |
| Resolution                  | 0.1 pH        | 1 mV            | 1 ppm                        | 0.01 ppt                     |
| TDS factor                  |               |                 | adjustable; 0.71 default     | adjustable; 0.71 default     |
| Battery requirements        | 4, AAA        | 4, AAA          | 4, AAA                       | 4, AAA                       |
| Casing IP rating            | IP67          | IP67            | IP67                         | IP67                         |
| Backlight                   | No            | No              | No                           | No                           |
| Article number              | 9531000       | 9531100         | 9531200                      | 9531300                      |

| Specifications              | POCKET PRO Conductivity LR | POCKET PRO Conductivity HR | POCKET PRO Salt | POCKET PRO Temperature   |
|-----------------------------|----------------------------|----------------------------|-----------------|--|
| Parameter                   | Conductivity, Temp         | Conductivity, Temp         | Salinity, Temp  | Temperature  |
| Operating temperature range | 0 - 50 °C                  | 0 - 50 °C                  | 0 - 50 °C       | 0 - 50 °C  |
| Range                       | 0 - 1990 μS/cm             | 0.0 - 19.99 mS/cm          | 0 - 10.00 ppt   | -15 - 170°C  |
| Accuracy                    | 1% FS                      | 2% FS                      | 1% FS           | 1°C  |
| Resolution                  | 1 μS/cm                    | 0.01 mS/cm                 | 0.01 ppt        | 0.1°C  |
| TDS factor                  |                            |                            |                 |  |
| Battery requirements        | 4, AAA                     | 4, AAA                     | 4, AAA          | 4, AAA   |
| IP rating                   | IP67                       | IP67                       | IP67            | IP67   |
| Backlight                   | No                         | No                         | No              | No   |
| Article number              | 9531400                    | 9531500                    | 9531600         | <b>9531700</b> 9531701 is replacement sensor for POCKET PRO Temp |

### **POCKET PRO+ Testers**



The POCKET PRO+ series offers all the benefits of the POCKET PRO standard testers, plus convenient backlight and replaceable sensors. Have confidence in your results with superior probe and calibration diagnostics.

NEW

| Specifications                    | POCKET PRO+ pH   | POCKET PRO+ ORP  | POCKET PRO+ Multi 1  | POCKET PRO+ Multi 2  |
|-----------------------------------|--|--|--|--|
| Parameter                         | pH, Temp   | ORP, Temp  | Conductivity, TDS, Salinity, Temp  | pH, Conductivity, TDS, Salinity, Temp  |
| Operating<br>temperature<br>range | 0 - 50 °C  | 0 - 50 °C  | 0 - 50 °C  | 0 - 50 °C  |
| Range                             | 0.00 - 14.00 pH  | -999 to +999 mV  | Cond: Auto-ranging (0.0 - 199.9 µS/cm; 200 - 1999 µS/cm; 2.00 - 19.99 mS/cm) TDS: Auto-ranging (0.0 - 99.9 ppm; 100 - 999 ppm; 1.00 - 10.00 ppt) Sal: Auto-ranging (0.00 - 10.00 ppt; 0.00 - 1.00%) Temp: 0.0 - 50°C   | pH: 0.00 - 14.00 Cond: Auto-ranging (0.0 - 199.9 μS/cm; 200 - 1999 μS/cm) TDS: Auto-ranging (0.0 - 99.9 ppm; 100 - 999 ppm; 1.00 - 10.00 ppt) Sal: Auto-ranging (0.00 - 10.00 ppt; 0.00 - 1.00%) Temp: 0.0 - 50°C  |
| Accuracy                          | ± 0.01 pH  | 2 mV   | Cond: ± 1%<br>TDS: ± 1%<br>Sal: ± 1%<br>Temp: ± 0.5°C  | pH: ±0.01 pH<br>Cond: ±1%<br>TDS: ±1%<br>Sal: ±1%<br>Temp: ±0.5 °C   |
| Resolution                        | 0.01 рН  | 1 mV   | Cond: 0.1 µS/cm from 0.0 - 199.9 µS/cm; 1 µS/cm from 200 - 1999 µS/cm; 0.01 mS/cm from 2.00 - 19.99 mS/cm TDS: 0.1 ppm from 0.0 - 99.9 ppm; 1 ppm from 100 - 999 ppm; 0.01 ppt from 0.00 - 10.00 ppt Sal: 0.01 ppt from 0.00 - 10.00 ppt; 0.01% from 0.00 - 1% Temp: 0.1°C | pH: 0.01 pH  Cond: 0.1 µS/cm from 0.0 - 199.9 µS/cm; 1 µS/cm from 200 - 1999 µS/cm; 0.01 mS/cm from 2.00 - 19.99 mS/cm  TDS: 0.1 ppm from 0.0 - 99.9 ppm; 1 ppm from 100 - 999 ppm; 0.01 ppt from 0.00 - 10.00 ppt  Sal: 0.01 ppt from 0.00 - 10.00 ppt; 0.01% from 0.00 - 1%  Temp: 0.1°C |
| TDS factor                        |  |  | adjustable; 0.71 default   | adjustable; 0.71 default   |
| Battery requirements              | 4, AAA   | 4, AAA   | 4, AAA   | 4, AAA   |
| Casing IP rating                  | IP67   | IP67   | IP67   | IP67   |
| Backlight                         | Yes  | Yes  | Yes  | Yes  |
| Article number                    | 9532000  | 9532100  | 9532700  | 9532800  |
|                                   | 9532001 is replacement sensor for POCKET PRO+pH Tester | 9532101 is replacement sensor for POCKET PRO+ ORP Tester | 9532701 is replacement sensor for POCKET PRO+Multi 1 Tester  | 9532801 is replacement sensor for POCKET PRO+Multi 2 Tester  |

# pH Buffer and Conductivity Standard Solutions

Article number

### **pH Standard Solutions**

Description

| can, guaranteed shelf life, with DKD certificate, and traceable to standard reference materials with given tolerances. |         |  |
|--|---------|--|
| pH 1.679 ±0.010 at 25°C, 500 mL  | S11M001 |  |
| pH 4.005 ±0.010 at 25°C, 500 mL  | S11M002 |  |
| pH 6.865 ±0.010 at 25°C, 500 mL  | S11M003 |  |
| pH 7.000 ±0.010 at 25°C, 500 mL  | S11M004 |  |
| pH 7.413 ±0.010 at 25°C, 500 mL  | S11M005 |  |
| pH 9.180 ±0.010 at 25°C, 500 mL  | S11M006 |  |
| pH 10.012 ±0.010 at 25°C, 500 mL S11M007   |         |  |
| pH 12.45 ±0.05 at 25°C, 500 mL S11M008   |         |  |
| Quality huffer solutions Ready-to-use huffer solutions in hottles  |         |  |

| Quality buffer solutions. Ready-to-use buffer solutions | in bottles, |  |
|---|-------------|--|
| with and without colour coding                          |             |  |
| mII 4 04 Park F00 ml                                    | 0000440     |  |

| pH 4.01 Red, 500 mL                           | 2283449 |
|---|---------|
| pH 7.00 Yellow, 500 mL                        | 2283549 |
| pH 10.01 Blue, 500 mL                         | 2283649 |
| pH 4.01 No colour code, 500 mL                | 1222349 |
| pH 7.00 No colour code, 500 mL                | 1222249 |
| pH 10.01 No colour code, 500 mL               | 1222149 |
| pH 1.09 Technical buffer solution (DIN 19267) | S11M009 |
| pH 4.65 Technical buffer solution (DIN 19267) | S11M010 |
| pH 9.23 Technical buffer solution (DIN 19267) | S11M011 |

### SINGLET buffer solutions. Buffer solutions in individually sealed airtight pouches, colour coded, 25 mL/pouch

| SINGLET Single use pH buffer solution, pH 4.01, 20 pcs  | 2770020 |
|---|---------|
| SINGLET Single use pH buffer solution, pH 7.00, 20 pcs  | 2770120 |
| SINGLET Single use pH buffer solution, pH 10.01, 20 pcs | 2770220 |

### **Conductivity Standard Solutions**

| Description   | Article number |  |
|---|----------------|--|
| SINGLET buffer solutions. Buffer solutions in individually sealed airtight pouches, colour coded, 25 mL/pouch |                |  |
| SINGLET Single use conductivity standard,<br>147 µS/cm, 20 pcs  | 2771320        |  |
| SINGLET Single use conductivity standard, 2771420 1413 µS/cm, 20 pcs  |                |  |
| SINGLET Single use conductivity standard,<br>12.88 mS/cm, 20 pcs  | 2771520        |  |

| Description   |                    | Volume /<br>Package size | Article number |
|---|--------------------|--------------------------|----------------|
| Certified conductivity standard solutions. Supplied in airtight sealed can, guaranteed shelf life, with certificate and traceable to standard reference materials |                    |                          |                |
| KCI 1 D   | 111.3 mS/cm ±0.5%  | 500 mL                   | S51M001        |
| KCI 0.1 D   | 12.85 mS/cm ±0.35% | 500 mL                   | S51M002        |
| KCI 0.01 D  | 1408 μS/cm ±0.5%   | 500 mL                   | S51M003        |
| NaCI 0.05%  | 1015 μS/cm ±0.5%   | 500 mL                   | S51M004        |
| NaCl solutions  |                    |                          |                |
| 491 mg/L<br>as NaCl   | 1000 ±10 μS/cm     | 100 mL                   | 1440042        |
| 85.47 mg/L<br>as NaCl   | 180 ±10 μS/cm      | 100 mL                   | 2307542        |
| 1 000 mg/L<br>as NaCl   | 1990 ±20 μS/cm     | 100 mL                   | 210542         |
| 10 246 mg/L<br>as NaCl  | 18000 ±50 μS/cm    | 100 mL                   | 2307442        |
| Molar KCI solut   | ions               |                          |                |
| KS 910 KCI<br>0.1 M   | 12.88 mS/cm        | 500 mL                   | C20C250        |
| KS 920 KCI<br>0.01 M  | 1.413 mS/cm        | 500 mL                   | C20C270        |
| KS930 KCI<br>0.001 M  | 146.9 μS/cm        | 500 mL                   | C20C280        |
|   | 1                  |                          |                |








62 SERVICES www.hach-lange.com

### **HACH LANGE Services**

With HACH LANGE, you receive instruments, reagents and services from a single source - offering expertise from development to consulting. Our experienced team can provide you with expert advice to ensure your application satisfies all your requirements.



### From beginning to end

At HACH LANGE, we support the use of our instruments throughout their life-cycle, from initial start-up to disposal. Combined with our extensive range of services, we take care of you directly – with locally based professionals.

### Various service packages

HACH LANGE offers comprehensive service packages to suit customer needs. Our flexible service packages provide the necessary assurance for smooth operation of the instrument to deliver top performance for both routine laboratory tasks to demanding photometry applications.

### Area-wide customer support

In over 20 countries in Europe, HACH LANGE is represented by its own subsidiaries with a qualified field service team, well-trained service technicians, and an experienced support team on the hotline. Available by phone, fax and e-mail, we offer experienced employees who speak your language and will be glad to help.

# HACH LANGE Services – What is behind it? More than you think.



#### On-Site

- ► Customer service training
- ► Intensive training
- ► Symposia
- ► Customer Advisory Board
- **►** Exhibitions



#### **Device**

- ► European on-site
- ► Service hotline
- ► Device instruction
- ► Installation Qualification / Operational Qualification (IQ / OQ)
- ▶ Service contracts



#### Information

- ► Hotlines
- ► Website
- ► E-shop
- ► Order tracking
- ► Configurators
- ► Newsletters
- ► Practice reports



#### **Environmental**

- ► Minimal chemical use
- ► Free redemption of used reagents
- ► Detailed product labeling
- ► Used reagents and instruments are professionally recycled
- ► Comprehensive services: collection, treatment and proper disposal



### **Quality and security**

- ► Standard solutions
- ► Test equipment
- ► Proficiency test
- ► Detailed, understandable documents
- ► Hazard information



### Consulting

- ► Expert advice from the first contact until long after purchase
- ► Detailed needs analysis for your situation and operation requirements
- ► Complete knowledge of regional specific needs
- ➤ Solutions to meet the requirements of your industry
- ► Results that will ensure quality and save time



64 SERVICES www.hach-lange.com

## Everything you need is just one click away

The HACH LANGE web portal helps you find what you need, quickly and easily. You benefit from: quick response to all your questions plus the latest analysis and monitoring news; tailored solutions for your everyday tasks; innovative products and a whole host of services designed with the customer in mind!



#### **Get started immediately**

You will feel completely at home: The uncluttered, comprehensive layout of all elements plus the intuitive user navigation concept of our web portal makes finding information and ordering products a breeze!

### Get things done more easily

Quick order lists for items you require on a regular basis ensure that the ordering processing is as convenient for you as possible. Orders are completed quickly, saving your valuable time.

### Find what you need, faster

With the intuitive search function, you can navigate by industry, parameter, application and product group.

### Keep up-to-date, 24/7

In addition to catching up with the latest news, industry trends, product innovations and technical developments, you can download safety data sheets, datailed product brochures, application reports and much more...



# **Parameter Index**

| Parameter                      | Page(s)                                    |
|--------------------------------|--|
| Acid capacity                  | 7  |
| Alcohol                        | 7  |
| Alkalinity                     | 21   |
| Aluminum                       | 7, 12, 20, 37                              |
| Ammonia12, 1                   | 6, 18, 20, 21, 37, 46, 47, 49, 54, 57, 58  |
| Ammonium                       | 7, 20, 28, 29, 37, 46, 47, 49, 54, 57      |
| Ammonium compounds, quaternary | 12   |
| AOX                            | 7, 21, 23                                  |
| Barium                         | 12   |
| Benzotriazole                  | 12   |
| Bitter units                   | 7  |
| BOD <sub>5</sub>               | 7, 21, 23, 41, 46, 47                      |
| Boron                          | 8, 12                                      |
| Bromine                        | 12, 15, 37                                 |
| Cadmium                        |  |
| Carbonate/carbon dioxide       | 8  |
| Chloramine, mono               | 12, 37                                     |
| Chloride                       |  |
| Chlorine                       | 8, 16, 18, 21, 37                          |
| Chlorine dioxide               | 8, 12, 15, 16, 37                          |
| Chlorine, free                 |  |
| Chlorine, total                | 12, 13, 14, 15, 37                         |
| Chromic acid baths             | 8  |
| Chromium                       | 8, 12, 20, 22, 37                          |
| Cobalt                         | 12, 16, 37                                 |
| COD                            |  |
| Colour                         | 21   |
| Conductivity21, 46, 4          | 17, 48, 52, 53, 54, 55, 56, 58, 59, 60, 61 |
| Copper                         | 9, 12, 15, 20, 37                          |
| Copper baths, acidic           | 9  |
| Cyanide                        | 9, 12, 22, 23                              |
| Cyanuric acid                  | 12   |
| Fluoride                       | 9, 15, 16, 37, 46, 47, 49, 54, 57          |
| Formaldehyde                   | 9, 16, 37                                  |
| Hardness, Ca and Mg            | 16   |
| Hydrazine                      | 16   |
| Hydrogen peroxide              | 16   |
| lodine                         | 12, 15                                     |
| Iron                           | 9, 12, 13, 15, 16, 20, 21, 37              |
| Lead                           | 9, 16, 20, 21, 37                          |
| Magnesium                      | 9  |
| Manganese                      | 13, 17, 20, 37                             |
| Menthol                        | 9  |
|                                | 17   |
| Molybdate                      | 13, 37                                     |
|                                | 9, 13, 37                                  |
|                                |  |

| Parameter                       | Page(s)                                       |
|---------------------------------|---|
| Nickel                          | 9, 13, 16, 20, 37                             |
| Nickel baths, acidic            | 9   |
| Nitrate                         | 3, 15, 18, 20, 28, 29, 37, 46, 47, 49, 54, 57 |
| Nitrite                         | 10, 13, 15, 18, 20, 21, 28, 29                |
| Nitrogen, total                 | 10, 13, 19, 20, 28, 29, 37                    |
| Nitrogen, total, Kjeldahl       | 22  |
| Organic acids                   | 10, 23  |
| ORP/Redox                       | 46, 47, 49, 52, 53, 54, 55, 56, 59, 60        |
| Oxygen scavengers               | 13  |
| Oxygen, dissolved               | 15, 37, 46, 47, 48, 52, 53, 55, 58            |
| Ozone                           | 15, 16, 37                                    |
| pH17, 37, 46                    | 6, 47, 48, 52, 53, 54, 55, 56, 58, 59, 60, 61 |
| Phenols                         | 10, 13, 23                                    |
| Phosphate                       | 10, 13, 15, 19, 20, 21, 29, 37                |
| Phosphonates                    | 13, 37  |
| Phosphorus, ortho               | 10, 15, 19, 20, 37                            |
| Phosphorus, total               | 17, 19, 28                                    |
| Photometric lodine sample (PIS) | 10  |
| Potassium                       | 10, 13, 20                                    |
| Reducing agent                  | 17  |
| Resistivity                     | 46, 47  |
| Salinity                        | 46, 47, 52, 53, 54, 60                        |
| Silica                          | 13, 17, 21, 37                                |
| Silver                          | 10, 13  |
| Sludge activity                 | 11, 23  |
| Sodium                          |   |
| Starch                          | 11  |
| Sulphate                        | 11, 13, 20, 21, 37                            |
|                                 | 11, 17  |
| Sulphite                        | 11, 17  |
|                                 | 11, 21, 23                                    |
| Surfactants, cationic           | 11, 21, 23                                    |
|                                 | 11, 21, 23                                    |
|                                 | 17  |
|                                 | 46, 47, 52, 53, 54, 59, 60                    |
|                                 | 11  |
|                                 | 11, 19, 20, 23                                |
|                                 | 12  |
| -                               | 17  |
|                                 | 42  |
|                                 | 11  |
|                                 | 17  |
|                                 | 11, 16  |
|                                 | 11  |
|                                 | 11, 13, 20, 37                                |
| Zirconium                       | 11  |



66 HOW TO ORDER www.hach-lange.com

### Order online

Register for the HACH LANGE web shop and you will have access to prices, brochures, application reports and much more... You can also view previous orders and track any current orders via www.hach-lange.co.uk or www.hach-lange.ie.



#### **Improved search**

Our improved search function will allow you to easily browse through our web shop by industries, applications, parameters, or product groups. Alternatively, find products via search box. Simply enter the appropriate keyword or article number. If you don't find what you are looking for, just email us.

### Intuitive ordering

Once you receive your personal login details, you can order online. Go to the product detail page and choose ADD TO BASKET. You will then be guided through the ordering process in five easy steps.

### Speed up the process

Do you know the article number of the product you want to order? Choose the quick order list at MY ACCOUNT. Insert the article number as well as the quantity you want to order. One further click and the product is saved in your basket and ready to order.

#### Save time

Do you frequently order identical or similar products? Make use of our wish list for regular ordering. To create your personal wish list, click on the product detail page and select ADD TO WISHLIST. The item is now saved and you can continue shopping. To buy the product, click on ADD TO BASKET and you will be guided through the purchasing process.









### Looking forward to hearing from you!

Contact us to place an order, request product information, technical support, arrange a service contract or service engineer visit. Enquire about our training courses and workshops. We can also provide the disposal of your used reagents.



### **United Kingdom:**

HACH LANGE LTD
Pacific way
Salford
Manchester M50 1DL
United Kingdom
Phone: +44 (0)161 872 1487

Fax: +44 (0)161 848 7235

Order/delivery enquiries:

Quotes/enquiries: sales@hach-lange.co.uk

orders@hach-lange.co.uk

Technical support: uktechnicalsupport@ ach-lange.co.uk

### Ireland:

HACH LANGE LTD Unit 1 Chestnut Road Western Industrial Estate Dublin 12 Ireland Phone: +353 (0)1 460 2522

Fax: +353 (0) 1 450 9337

Order/delivery enquiries: orders@hach-lange.ie

Quotes/enquiries: sales@hach-lange.ie

Technical support: ietechnicalsupport@hach-lange.ie

#### **Free UK Service Hotline:**

The UK Service Hotline provides essential contact with a trained engineer to identify and resolve instrument failures remotely. Understanding the needs of our customers has always been a priority for HACH LANGE. In our continuing commitment to improving quality, reliability and service of our instrument range.

### International: +49 (0)211 52 88 0

\* Charges will apply if calling from Republic of Ireland.

### Service Hotline:

United Kingdom: Phone: 0800 9555 075 Email: servicehotlineuk@ hach-lange.co.uk

Ireland:

Phone 0800 9555 075\* Email: servicehotlineie@ hach-lange.ie







HACH LANGE GMBH

Willstaetterstrasse 11 D-40549 Duesseldorf Tel. +49 (0)211 5288-0 Fax +49 (0)211 5288-143 info@hach-lange.de www.hach-lange.com

HACH LANGE LTD

Pacific Way
Salford
GB-Manchester, M50 1DL
Tel. +44 (0)161 8721487
Fax +44 (0)161 8487324 info@hach-lange.co.uk www.hach-lange.co.uk

HACH LANGE LTD
Unit 1, Chestnut Road
Western Industrial Estate
IRL-Dublin 12
Tel. +353(0)1 4602522
Fax +353(0)1 4509337
info@hach-lange.ie www.hach-lange.ie

