


Cross flow filtration method



Contents

01	Introduction	3	06	Cell and lysate clarification	47
	1.1 What is cross flow filtration?	5		6.1 Introduction	49
	1.2 Key features of CFF	6		6.2 Membrane and cartridge selection	51
	1.3 CFF application areas	7		6.3 Operating parameters	53
02	Cross flow filtration systems	8		6.4 Clarification strategy examples	53
	2.1 System configuration	10		6.4.1 Mammalian cells	54
	2.2 Filters for cross flow filtration	12		6.4.2 Bacterial cells	55
	2.3 System volumes and process capacity	14	07	6.4.3 Yeast	56
	2.4 Systems from Cytiva	15		Concentration and diafiltration	57
03	Process design and operation	16		7.1 Introduction	59
	3.1 Process considerations	18		7.2 Product and process considerations	61
	3.2 Filter selection	19		7.3 Membrane selection	63
	3.3 Filter preparation	21	08	Appendix	64
	3.4 Operating parameters	23		Appendix A	65
	3.5 Recovering product	26		Abbreviations and glossary	
	3.6 Cleaning and testing filters	27		Glossary	66
	3.7 Scaling up processes	29		Index	69
04	Optimizing CFF processes	30			
	4.1 Optimizing process parameters	32			
	4.2 Optimizing yield	36			
05	Cell harvesting	38			
	5.1 Introduction	40			
	5.2 Cell harvesting process	41			
	5.3 Membrane and cartridge selection	44			
	5.4 Operating parameters	45			

01

Introduction

About this handbook

The Cross flow filtration method handbook gives a general introduction to the principles and applications of cross flow filtration using systems and filters from Cytiva. Detailed instructions for preparing, using, cleaning and storing particular filters are provided with each filter.

About this chapter

This chapter introduces the principles, basic terminology and applications of cross flow filtration.

In this chapter

This chapter contains the following sections:

Section	See page
1.1 What is cross flow filtration?	5
1.2 Key features of CFF	6
1.3 CFF application areas	7

1.1 What is cross flow filtration?

Cross flow filtration (CFF, also known as **tangential flow filtration** TFF) is a filtration technique in which the starting solution passes tangentially along the surface of the filter. A pressure difference across the filter drives components that are smaller than the pores through the filter. Components larger than the filter pores are retained and pass along the membrane surface, flowing back to the feed reservoir (see Fig 1.1). CFF is simple in concept, but its proper execution requires detailed knowledge and good filtration technique.

Solution that is directed to the membrane surface is called the **feed**. Solution that passes along the membrane surface and back to the feed reservoir is the **retentate**. This solution is usually pumped back to the feed reservoir and recirculated. Solution that passes across the membrane is the **permeate**.

Other terms used in CFF are explained in *Appendix A Abbreviations and glossary*, on page 65 and 66.




Fig 1.1. Principle of cross flow filtration.

1.2 Key features of CFF

Features

A key feature of CFF is the flow of fluid along the membrane surface that sweeps away the build up of material on the filter surface and reduces fouling of the filter. In addition, retentate solution can easily be recirculated, allowing thorough processing of large volumes of solution.

Differences between CFF and normal flow filtration

CFF differs from normal flow filtration (where the starting material is simply passed through the filters, also known as **dead-end filtration** or **direct flow filtration**) in three main ways:

- CFF filters use membranes exclusively, while conventional filtration may use membranes, paper, or other materials such as glass fiber to separate components in a feed stream
- CFF supports recirculation of the retentate solution. In normal flow filtration, the feed usually passes once through the filter.
- In a CFF system, the retentate remains as a solution and may be recovered directly. Retentate recovery is relatively uncommon in normal flow filtration and requires resuspension of material collected on the filter.

1.3 CFF application areas

CFF is used in research, product development, and production in the biopharmaceutical and medical industries.

Table 1.1 summarizes the characteristics of typical CFF applications.

Table 1.1. Typical uses of cross flow filtration

Application	Comments	Used in
Cell harvesting	Separation of cells from fermentation broth Cells are recovered in the retentate	Upstream processing
Cell or lysate clarification	Separation of target molecules from intact cells, cell debris and molecular aggregates Target molecules are recovered in the permeate	
Product fractionation	Separation of components on the basis of molecular size	Downstream processing
Product concentration	Concentration of product solution prior to further processing by removal of solvent and small molecules Product is recovered in the retentate	Downstream processing Product formulation
Diafiltration	Buffer exchange Product is recovered in the retentate	

Documents describing some examples of CFF applications are available from Cytiva.

02

Cross flow filtration systems

About this chapter

This chapter describes the essential components of CFF systems and introduces systems from Cytiva.

In this chapter

This chapter contains the following sections:

Section	See page
2.1 System configuration	10
2.2 Filters for cross flow filtration	12
2.3 System volumes and process capacity	14
2.4 Systems from Cytiva	15

2.1 System configuration

The illustration to the right shows a generalized scheme for the basic configuration of a CFF system.

Pumps and valves

Liquid flow in the system is maintained by one or more pumps in the process lines (only the feed pump is shown in Fig 2.1):

- A **feed pump** maintains the flow of feed into the filter
- A **retentate pump** may be used to maintain and control flow of retentate back into the feed reservoir
- A **permeate pump** may be used to control flow of permeate from the filter. The permeate pump, if used, should not create negative pressure on the permeate side of the filter (i.e., the pump flow rate must be less than the spontaneous permeate flux).
- A **transfer pump** may be used in washing and diafiltration applications to add liquid (usually buffer) to the feed reservoir at a controlled rate

Flow in the various process lines may also be regulated by valves with flow restrictors. Together, the controlled pump rates and valve restriction create the pressure across the membrane that drives the filtration process.




Fig 2.1. Basic configuration of a CFF system.

Pressure sensors

Pressure sensors in the feed and retentate lines are essential for monitoring and controlling the process. A pressure sensor on the permeate side may also be used to monitor permeate pressure.

Flow sensors

Measurement of flow rates for feed solution, retentate and/or permeate and any addition of fluid to the feed reservoir is necessary for monitoring and controlling process conditions during filtration. Flow sensors are placed at strategic points in the system.

Reservoir level sensors

A level sensor in the feed reservoir monitors and controls the amount of liquid in the reservoir.

Air sensors

An air sensor located in the feed stream allows continuous monitoring for air bubbles in the feed prevents the introduction of air into the filter. An air sensor can also be used in the transfer line to detect when the transfer reservoir is empty.

Additional sensors

Temperature, pH, UV absorbance and conductivity sensors may be included in the system according to the requirements of the specific process.

2.2 Filters for cross flow filtration

Filters for CFF may be classified according to filter configuration or filter pore size.

Filter configuration

Two basic filter configurations are generally used (see Fig 2.2):

- In **cartridge filters** (often called **hollow fiber filters**), the membrane forms a set of parallel hollow fibers. The feed stream passes through the lumen of the fibers and the permeate is collected from outside the fibers. Cartridges are characterized in terms of fiber length, lumen diameter and number of fibers, as well as filter pore size. Cartridges from Cytiva have lengths of 30 to 110 cm and lumen diameters of 0.5 to 1.75 mm.
- In **cassette filters**, several flat sheets of membrane are held apart from each other and from the cassette housing by support screens. The feed stream passes into the space between two sheets and permeate is collected from the opposite side of the sheets. Cassettes are characterized in terms of flow path length and channel height, as well as membrane pore size. The channel height is determined by the thickness of the support screen.

Both cartridges and cassettes are constructed from materials chosen for mechanical strength, chemical and physical compatibility, and low levels of extractable and/or toxic compounds.

The total membrane surface area is one of the factors determining how much feed can be handled in a process run. Typical values for the selection of filters by membrane surface area are as follows:

- For microfiltration, 30 to 100 liters of feed per m² of membrane surface area
- For ultrafiltration, 100 to 200 liters of feed per m² of membrane surface area




Fig 2.2. Configuration of hollow fiber cartridges and flat sheet cassettes.

Filter pore size

Two classes of filters are distinguished on the basis of filter pore size, which also distinguishes between application areas. In both cases, the pore size is expressed as an average value. The range of actual pore sizes in a given filter determines the selectivity of the filter.

Microfiltration filters

Filters with pore sizes of 0.1 µm and larger are classified as microfilters. In CFF applications the pore size is usually in the range 0.1 to 1 µm. These membranes are used for separation of cultured cells from the growth medium (broth), as well as for removal of particulate material in numerous biopharmaceutical processes.

Ultrafiltration filters

Ultrafiltration membranes have pore sizes in the range 20 to 100 nm, and are generally characterized in terms of the nominal molecular weight cutoff (NMWC), which is the molecular weight of the largest globular protein that can pass through the membrane. NMWC values range from 1 to 100 kD (kiloDalton). These filters are used for concentrating and fractionating protein streams, virus concentration, desalting and buffer exchange. The objective of most ultrafiltration processes is to retain soluble macromolecules such as proteins above a certain size, while allowing smaller molecules such as salts, amino acids, and mono- or disaccharides to pass through the membrane.

Filters from Cytiva

Cytiva supplies cartridges for both microfiltration and ultrafiltration, and cassettes for ultrafiltration. The choice of filter type to use for a given application is made in the first place on the filtration requirements: thus cartridges may be used for all applications, while cassettes are suitable for handling proteins.

See *Section 3.2 Filter selection*, on page 19 for more detailed guidelines concerning filter selection. Additional information is available from Cytiva.

2.3 System volumes and process capacity

Minimum working volume

The minimum working volume of a CFF system represents the amount of feed/retentate fluid required to operate the system at the desired cross flow rate without drawing air into the feed pump. The minimum working volume is determined by the design of the system (feed and retentate tubing volume, reservoir bottom design), the filter hold-up volume, and the cross flow rate. It is important to consider the minimum working volume of a system in the design of a CFF process and in particular to confirm that the final target retentate volume is not less than the system's minimum working volume.

Hold-up volume

The term **hold-up volume** refers to the volume of liquid in the filtration system. For some purposes it is useful to distinguish between the **filter hold-up volume** (the volume in the filter itself) and the **system hold-up volume** (the volume in the system tubing and pumps). The filter and system should be chosen with the smallest hold-up volume that is compatible with other performance requirements in the process.

Process capacity

The system process capacity (the volume of starting material that can be processed in one run) should be chosen in relation to the planned volume of starting material. Process capacity is partly a function of system size and design, but also varies according to the tendency of starting material to foul the filter. Using a high capacity system for a small sample volume will lead to unnecessary loss of material in the system dead volume. For processes that will be scaled up for production, it will be necessary to switch between different systems one or more times during process development (see Table 2.1). CFF processes using systems from Cytiva are scalable, making transitions between systems straightforward.

2.4 Systems from Cytiva

CFF systems from Cytiva cover a full range of filtration capacity from laboratory research to production, and offer different levels of monitoring and control from manual to fully automated systems. Customized systems can also be provided.

Table 2.1 summarizes the characteristics of the available systems. Additional information including system application matrix and selection guides is available from Cytiva.

Table 2.1. Cross flow filtration systems from Cytiva

System	Scale	Applications
ÄKTACrossflow™	Laboratory	Research
		Process development
ÄKTA™ flux	Laboratory to pilot	Research
		Process development
		Small-scale production
UniFlux™	Pilot to production	Manufacturing

03

Process design and operation

About this chapter

This chapter describes design and operation of a CFF process.

In this chapter

This chapter contains the following sections:

Section	See page
3.1 Process considerations	18
3.2 Filter selection	19
3.3 Filter preparation	21
3.4 Operating parameters	23
3.5 Recovering product	26
3.6 Cleaning and testing filters	27
3.7 Scaling up processes	29

3.2 Filter selection

The choice of filter between microfiltration and ultrafiltration and between cartridges and cassettes is governed by the nature of the application.

Application	Separation principle	Technique and filter
Cell harvesting	Cells are separated from soluble molecules	Microfiltration (cartridge)
Cell or lysate clarification	Cells and cell debris are separated from soluble molecules	
Protein fractionation	Macromolecules are separated on the basis of size	Ultrafiltration (cartridge or cassette)
Concentration and diafiltration	Macromolecules are separated from low molecular weight buffer components	

Within that framework, filters are chosen from consideration of selectivity, filter pore size and protein binding.

Filter selectivity

Filter selectivity describes the ability of a filter to separate particles or molecular species on the basis of size. Filters with a narrow pore size distribution will be highly selective, while a broader pore size distribution will give a less selective filter.

Harvesting and cell clarification applications involve separation of relatively large particles (cells and/or cell debris) from macromolecules, so that high selectivity is generally not required. Lysate clarification on the other hand may make more stringent demands, since the lysate will contain a wide range of proteins and other macromolecules. The most important factor is that the target protein can pass freely through the filter so that yields are not compromised.

Filter selectivity is also not critical for diafiltration, where a macromolecular product is separated from buffer components. The controlling factor here is that the target protein is completely retained by the filter, so that product is not lost in the permeate. An ultrafiltration filter with a nominal molecular weight cut-off (NMWC) that is 3× to 5× less than the molecular weight of the target molecule is generally recommended.

Filter pore size

All filters will tend to become fouled (blocked by particulate material accumulating in the filter pores), especially during applications such as lysate clarification involving particulate starting material. Fouling will shorten the time for which a filter can be used before it must be cleaned, and therefore restrict the maximum processing capacity per run. The choice of filter pore size is important for minimizing filter fouling. In general, filters with smaller pores will show less tendency to fouling, because particulate material cannot penetrate and block the pores.

Use the guidelines in the table below for selecting filters for harvesting and clarification applications:

Application	Filter pore size
Mammalian cell harvesting	0.2 to 0.65 µm
Yeast and bacterial cell harvesting	0.1 µm
Cell and lysate clarification	NMWC about 10× target molecule size

Protein binding

The level of protein binding depends upon the material in the filter and the protein characteristics, and increases with increasing hydrophobicity. Protein binding is seen as a reduction in yield that is not accounted for by product remaining in the retentate. Normally, protein binding remains insignificant at the laboratory scale, but for low NMWC ultrafiltration membranes it can be an indicator of a propensity towards filter fouling.

For both clarification and diafiltration applications, protein binding to the filter usually becomes an issue when attempting to process small amounts of protein. In such cases it is important to select a filter with low binding tendency for the target protein, and to choose buffer conditions that minimize binding.

3.3 Filter preparation

Preparation of filters for a process involves rinsing to remove storage solution and conditioning the filter with process buffer. For some processes the filter may need to be sanitized and depyrogenated before use. This section gives a brief overview of preparation procedures. More details are provided in the instructions accompanying each filter.

Sanitization and depyrogenation

Follow the steps below to sanitize and depyrogenate the filter if necessary.

Step	Action
1	Clean and rinse the filter thoroughly
2	Recirculate a solution of 0.1 to 0.5 M NaOH (pH 13) for 30 to 60 minutes at 30°C to 50°C
3	Drain the system thoroughly
4	Rinse the filter with clean water for 30 minutes

3.4 Operating parameters

A CFF process may be controlled by an interplay of several operating parameters according to the specific process requirements. The most important parameters are:

- Pressure at various points in the system
- Flow rate at various points in the system
- Process time

Pressure and flow rates

Liquid pressure and flow rates are essential factors for controlling and monitoring a CFF process.

Pressure

Pressure may be monitored in the feed stream, the retentate stream or the permeate stream. Two differential pressure measurements are generally used, ΔP and transmembrane pressure (TMP).

- ΔP is the difference in pressure between the feed and retentate streams, and can be used to control cross flow
- TMP represents the driving force for transfer of material across the filter, and is calculated as shown below.

TMP can be used to control flux.

$$\text{TMP} = \frac{\text{Feed pressure} + \text{Retentate pressure}}{2} - \text{Permeate pressure}$$

3.7 Scaling up processes

The ability to scale a process from the laboratory to manufacturing is a key factor in process development. Normally, the scale-up sequence is completed in multiple steps: lab scale to pilot scale, and pilot scale to production scale. Reasonable scale-up increments are typically 5 to 20 times.

Scaling up a process involves increasing the filter area in order to handle larger volumes of starting material without significantly changing process conditions. The following parameters should be kept constant where possible:

- Ratio of filter area to feed volume
- Fiber or cassette path length
- Channel height (cassettes) or lumen size (hollow fiber cartridges)
- Membrane characteristics (pore size, selectivity, materials)
- Cross flow rate per unit filter area
- TMP
- Temperature
- Feed concentration
- Process steps and sequence

04

Optimizing CFF processes

About this chapter

This chapter considers optimization of CFF processes.

In this chapter

This chapter contains the following sections:

Section	See page
4.1 Optimizing process parameters	32
4.2 Optimizing yield	36

Process parameters may be optimized with respect to several factors, such as capacity, total process time, product yield and purity and so on.

Temperature

For heat sensitive proteins the process solution temperature can be modified in a number of ways during processing to optimize product yield:

- Precondition the system with cooled buffer before starting
- Lower the feed temperature before beginning the filtration process
- Use chilled buffer during diafiltration
- Use low pressure and low pump speed to reduce heat generation in the flow path
- Use a low ratio of feed volume to filter surface area to shorten process time
- Place the system in a cold room
- Use a heat exchanger or tank jacket

Enzymatic action

Proteolytic enzymes released during cell culture or lysis may follow the target protein through the filtration process and cause degradation and loss of yield. The effect may be more pronounced in ultrafiltration or diafiltration applications, where both enzymes and target proteins may be concentrated together in the concentration layer at the filter surface or in the retentate.

Enzymatic activity may be reduced by:

- Including enzyme inhibitors in the sample
- Using a low ratio of feed volume to filter surface area to shorten process time
- Lowering the temperature of the process
- Adjusting buffer conditions (ionic strength, pH, metal ions, etc.) to minimize enzyme activity

05

Cell harvesting

About this chapter

This chapter considers the use of CFF for harvesting cultured cells.

In this chapter

This chapter contains the following sections:

Section	See page
5.1 Introduction	40
5.2 Cell harvesting process	41
5.3 Membrane and cartridge selection	44
5.4 Operating parameters	45

5.1 Introduction

Cell harvesting is the process of separating cells from the fermentation broth in which the cells are grown. Hollow fiber microfiltration or higher NMWC ultrafiltration cartridges may be used effectively for cell harvesting. The harvested cells are recovered in the retentate.

Successful cell harvesting relies on knowledge of parameters such as:

- Robustness of the cultured cells
- Starting volume and concentration of cells
- Desired finished concentration and volume
- Desired yield and integrity of the cells

5.2 Cell harvesting process

The cell harvesting process involves both concentration and washing of the cells.

Concentration

Cells are concentrated as part of the harvesting process (Fig 5.1).

The concentration factor that can be achieved depends on the starting concentration, and is usually limited by the fluidity of the concentrated cell suspension. Typical concentration factors are shown in Table 5.1.

Table 5.1. Typical concentration factors for different cell types

Cell type	Typical concentration factor
<i>E.coli</i>	5×
Yeast	2×
Mammalian	7 to 10×




Fig 5.1. Concentration is the first step in cell harvesting.

Washing

Cell harvesting usually includes a washing step to ensure effective removal of broth components from the cells (Fig 5.2).

After washing, the ideal end product would consist of the concentrated cells suspended in the buffer used to wash the cells. However, in practice the harvested cells in buffer can contain varying levels of unwanted elements such as precipitated proteins, enzymes, and cell debris.

The washing process is commonly a constant volume diafiltration process, in which buffer is added to the cell suspension at the same rate as the permeate flows across the membrane. Unlike centrifugal techniques where cells are packed in a dense cake or pellet, washing the cells in a buoyant state enables effective removal of contaminants.




Fig 5.2. Washing harvested cells after concentration.

Typical process steps

Typical steps in a cell harvesting method are outlined in Table 5.2:

Table 5.2. Typical steps in a cell harvesting method

Step	Operation	Purpose
Preparation	Rinsing	Rinse storage solution from filter
	Water flush	Flush cleaning solution from system
	Water flux test	Determine performance of filter before processing
	Buffer conditioning	Condition filter and system components before adding product
Harvesting	Cell harvesting (concentration)	Concentrate cells
	Cell washing (diafiltration)	Remove unwanted components
	Product recovery	Recover harvested cells
Finishing	Buffer flush	Flush residual product from system
	CIP	Recirculate cleaning solution
	Water flush	Flush cleaning solution from system
	Water flux test	Determine the performance of filter after use and cleaning
	Storage	Flush filter with storage solution to prevent bacterial growth in storage

06

Cell and lysate clarification

About this chapter

This chapter considers the use of CFF for clarification of cell cultures and lysates.

In this chapter

This chapter contains the following sections:

Section	See page
6.1 Introduction	49
6.2 Membrane and cartridge selection	51
6.3 Operating parameters	53
6.4 Clarification strategy examples	53

07

Concentration and diafiltration

About this chapter

This chapter considers the use of CFF for product concentration and diafiltration.

In this chapter

This chapter contains the following sections:

Section	See page
7.1 Introduction	59
7.2 Product and process considerations	61
7.3 Membrane selection	63

08

Appendix

