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(A) Determination of Physical Parameters

Quality Control:

The following general guidelines may be applied to many of the methods in this part :

- (1) Evaluate analyst performance for each method.
- (2) Determine competence by analyses of samples containing known concentrations.
- (3) Calibrate instruments and ensure that instrument measurements do not drift.
- (4) Assess the precision of analytical procedures by analyzing at least 10% of samples in duplicate. Analyze a minimum of one duplicate with each set of samples.
- (5) Determine bias of an analytical procedure in each sample batch by analysis of blanks, known additions with a frequency of at least 5% of samples, and, if possible, an externally provided standard.

pH Value (Method A-1, Electrometry)

Principle:

The pH of a water sample is determined by the measurement of a voltage produced between an electrode responsive to hydrogen ions (glass electrode) and a reference electrode. (usually calomel electrode) when both are immersed in the sample. A difference of 1 pH unit produces a potential change of 58.16 mV at 25°C. The pH-meter is calibrated to read a pH-value directly.

Equipment:

- (1) pH-meter (accuracy 0.1 pH units).
- (2) Combined glass/calomel, or glass/silver- silver chloride electrode.
- (3) Magnetic stirrer.
- (4) Teflon-coated stirring bar.

Reagents:

- (1) Standard pH buffer solutions (pH 4, 7 and 9).
- (2) Bidistilled water.
- (3) Soft tissue paper.

Procedure:

Rinse the electrode(s) with bidistilled water and dry by gentle wiping with a soft tissue paper. This should be done each time before applying the electrode to a new buffer solution or sample.

- (1) If manual temperature compensation is provided, adjust it to the temperature of the buffer or sample to be measured.
- (2) Calibrate the instrument by immersing the rinsed electrode(s) in a buffer solution of pH 7 and set the meter to read 7.00 by adjusting the asymmetry control.
- (3) Rinse the electrode and repeat using a buffer solution of pH 4 or pH 9. Select the possibility which is coming closest to the pH of the sample to be measured.
- (4) Adjust the reading to the appropriate number by the slope control and discard the portions of the buffer used in the calibration.
- (5) Immerse the rinsed combination glass electrode and temperature corrected electrode in the unknown sample solution and take a reading after 30 seconds.
- (6) Express the results to the nearest 0.1 pH unit. The temperature of the sample at the time of pH measurement should be indicated.

Precautions:

- (1) Before making any measurement, rinse the electrodes with bidistilled water and dry with tissue paper.
- (2) Always adjust the instrument to the temperature of the sample.
- (3) Ensure an adequate level of filling solution in the reference electrode.
- (4) For storage, always keep the electrode(s) in bidistilled water, never allow them to dry out.
- (5) Preferable the pH of any sample should be established at the sampling site, or as soon as possible after taking the sample.

- (6) No chemical preservation of the sample is permissible.
- (7) In case salinity is high, a sodium ion error at pH 10 and above has to be compensated.

References:

- (1) USEPA "Methods For Chemical Analysis of Water and Wastes", Ohio, 1983 (Method 150.1).
- (2) "Standard Methods for Examination of Water and Wastewater", 19th Edition, American Public Health Association, Washington 1995 (Method 4500-H+).

Acidity (Method A-4, Potentiometry)

Principle:

The pH of the sample is determined and a measured amount of standard acid is added, as needed , to lower the pH to 4 or less . Hydrogen peroxide is added, the solution boiled for several minutes, cooled, and titrated electrometrically with standard alkali to pH 8.2. The method covers the range from approximately 10 mg/L acidity to approximately 1000 mg/L as CaCO_3 , using a 50 ml sample.

Equipment:

- (1) pH meter.
- (2) Titration equipment including pipettes and burettes.

Reagents:

- (1) Hydrogen peroxide (H_2O_2 , 30% solution).
- (2) Standard sodium hydroxide , 0.02 N.
- (3) Standard sulfuric acid , 0.02 N.

Procedure:

- (1) Pipette 50 ml of the sample into a 250 ml beaker.
- (2) Measure the pH of the sample. If the pH is above 4.0, add standard 0.02 N sulfuric acid in 5.0 ml increments to lower the pH to 4.0 or less. If the initial pH of the sample is less than 4.0 , the incremental addition of sulfuric acid is not required.
- (3) Add 5 drops of 30% hydrogen peroxide.
- (4) Heat the sample to boiling and continue boiling for 2 to 4 minutes.
- (5) In some instances, the concentration of ferrous iron in a sample is such that an additional amount of hydrogen peroxide and slightly longer boiling time may be required.
- (6) Cool the sample to room temperature and titrate electrometrically with standard 0.02 N sodium hydroxide to reach a pH value of 8.2.

Calculations:

$$(1) \text{ Acidity, as CaCO}_3, (\text{mg/L}) = \frac{[(A \times B) - (C \times D)] \times 50000}{\text{ml sample}}$$

Where:

A = volume of standard sodium hydroxide used in titration (ml).

B = normality of standard sodium hydroxide.

C = volume of standard sulfuric acid (ml) used to reduce pH to 4 or less.

D = normality of standard sulfuric acid.

Precautions:

- (1) Suspended matter present in the sample, or precipitates formed during the titration may cause a sluggish electrode response.
- (2) This effect may be offset by allowing a 15-20 second pause between additions of titrant or by slow dropwise addition of titrant as the endpoint pH is approached.

References:

- (1) USEPA "Methods For Chemical Analysis of Water and Wastes", Ohio, 1983 (Method 305.1).
- (2) "Standard Methods for Examination of Water and Wastewater", 19th Edition, American Public Health Association, Washington 1995 (Method 2310 B).

Alkalinity (Method A-5, Potentiometry)

Principle:

An unaltered sample is potentiometrically titrated to pH 4.5. The sample must not be filtered, diluted, concentrated, or altered in any way. The method is suitable for all ranges of alkalinity.

Equipment:

- (1) pH meter (sensitivity 0.05 pH units). If automatic temperature compensation is not provided , make titration at $25 \pm 2^\circ\text{C}$.
- (2) Use appropriate sized vessel to keep the air space above the solution at a minimum.
- (3) Use a rubber stopper fitted with holes for the glass electrode, reference electrode (or combination electrode) and burette.
- (4) Magnetic stirrer, pipettes, flasks and other standard laboratory equipment.
- (5) Burettes, Pyrex 50, 25 and 10 ml.

Reagents:

- (1) Sodium carbonate solution, approximately 0.05 N: Place 2.5 ± 0.2 g Na_2CO_3 (dried at 250°C for 4 hours and cooled in desiccator) into a 1 liter volumetric flask and dilute to the mark.
- (2) Standard acid (sulfuric or hydrochloric), 0.1 N: Dilute 3.0 ml of concentrated H_2SO_4 or 8.3 ml of concentrated HCl to 1 liter with distilled water. Standardize versus 40.0 ml of 0.05 N Na_2CO_3 solution diluted with about 60 ml distilled water by titrating potentiometrically to pH of about 5. Lift electrode and rinse into beaker. Boil solution gently for 3-5 minutes under a watch glass cover. Cool to room temperature. Rinse cover glass into beaker. Continue titration to the pH inflection point . Calculate normality using:

$$N = \frac{A \times B}{53.00 \times C}$$

Where:

A = weight of Na_2CO_3 (g) dissolved in 1 liter.

B = volume of Na_2CO_3 solution used for titration (ml).

C = volume of acid used to reach inflection point(ml).

- (1) Standard acid (sulfuric or hydrochloric) , 0.02N : Dilute 200.0 ml of 0.1000 N standard acid to 1 liter with distilled water. Standardize by potentiometric titration of 15.0 ml 0.05 N Na_2CO_3 solution as above.

Procedure:

- (1) Sample size:

A. Use a sufficiently large volume of titrant (> 20 ml in a 50 ml burette) to obtain good precision while keeping

volume low enough to permit sharp end point .

- B. For $<1000 \text{ mg CaCO}_3 / \text{L}$, use a standard 0.02 N acid titrant.
- C. For $>1000 \text{ mg CaCO}_3 / \text{L}$, use a standard 0.1 N acid titrant.
- D. A preliminary titration is helpful.

(2) Potentiometric titration of normal alkalinity:

- A. Place 25- 50 ml of the sample solution in flask by pipetting wpipette tip near bottom of flask.
- B. Measure pH of sample.
- C. Add standard acid (0.02 or 0.1 N), being careful to stir thoroughly but gently to allow needle to obtain equilibrium.
- D. Titrate to pH 4.5 . Record volume of titrant.

(3) Potentiomettitration of low alkalinity:

- A. For alkalinity of $< 20 \text{ mg/L}$,titrate 100 - 200 ml as above using a 10 ml microburette and 0.02 N acid solution.
- B. Stop titration at pH in range of 4.3 - 4.7, record titrant volume and exact pH. Carefully add titrant to lower pH by exactly 0.3 pH units and record volume of titrant.

Calculations:

(1) Potentiometric titration to pH 4.5:

$$\text{Alkalinity as CaCO}_3, (\text{mg/L}) = \frac{A \times N \times 50,000}{\text{ml of sample}}$$

Where:

A = volume of standard acid consumed in the titration (ml).
B = normality standard acid.

(2) Potentiometric titration of low alkalinity :

$$\text{Total alkalinity, as CaCO}_3 (\text{mg/L}) = \frac{(2B - C) \times N \times 50,000}{\text{ml of sample}}$$

where :

B = volume (ml) of titrant consumed to first recorded pH.
C = total volume (ml) of titrant used to reach pH 0.3 units lower.
N = normality of acid.

Precautions:

- (1) The sample should be refrigerated at 4°C and run as soon as practical. Do not open sample bottle before analysis.
- (2) Substances, such as salts of weak organic and inorganic acids present in large amounts, may cause interference in the electrometric pH measurements.
- (3) For samples having high concentrations of mineral acids, such as mine wastes and associated receiving waters, titrate to an electrometric endpoint of pH 3.9.
- (4) Oil and grease, by coating the pH electrode, may also interfere, causing sluggish response.

References:

- (1) USEPA "Methods For Chemical Analysis of Water and Wastes", Ohio, 1983 (Method 301.1).
- (2) "Standard Methods for Examination of Water and Wastewater", 19th Edition, American Public Health Association, Washington 1995 (Method 2320 B).

Hardness

(Method A-6, Visual Titrimetry)

Principle:

Calcium and magnesium ions in the sample are sequestered upon the addition of disodium ethylenediamine tetraacetate (Na_2EDTA). The end point of the reaction is detected by means of Eriochrome Black T indicator, which has a red color in the presence of calcium and magnesium and a blue color when the cations are sequestered. The method is suitable for all concentration ranges of hardness.

Equipment:

Standard laboratory titrimetric equipment.

Reagents:

- (1) Buffer magnesium EDTA solution: Dissolve 1.179 g disodium EDTA (analytical reagent grade) and 780 mg $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ (or 644 mg $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$) in 50 ml distilled water. Add this solution to a 250 ml volumetric flask containing 16.9 g NH_4Cl and 143 ml of concentrated aqueous ammonia solution with mixing and dilute to the mark with distilled water. Store in a tightly stoppered plastic bottle; stable for approximately one month. Dispense with bulb operated pipette. Discard when 1 or 2 ml added to sample fails to produce a pH of 10.0 ± 0.1 at end point of titration.
 - (2) Indicator: Mix together 0.5 g Eriochrome Black T and 100 g NaCl .
 - (3) Standard EDTA titrant, (0.02) N: Place 3.723 g analytical reagent grade disodium ethylenediamine tetraacetate dihydrate, $\text{Na}_2\text{H}_2\text{C}_{10}\text{H}_{12}\text{O}_8\text{N}_2 \cdot 2\text{H}_2\text{O}$ in a 1 liter volumetric flask and dilute to the mark with distilled water. Check with standard calcium solution by titration. Store in polyethylene. Check periodically because of gradual deterioration.
- A. Standard calcium solution 0.02 N: Place 1.000 g anhydrous calcium carbonate (primary standard low in metals) in a 500 ml flask. Add , 200 ml distilled water. Boil for a few minutes to expel CO_2 . Cool , add a few drops of methyl red indicator and adjust to intermediate orange color by adding 3 N aqueous ammonia solution or (1 + 1) HCl as required. Quantitatively transfer to a 1 liter volumetric flask and dilute to mark with distilled water.
- B. Hydrochloric acid solution, (1 + 1).
- C. Methyl red indicator : Dissolve 0.10 g methyl red in distilled water in a 100 ml volumetric flask and dilute to the mark .

- D. Aqueous ammonia solution , (3 N): Dilute 210 ml of concentrated aqueous ammonia to 1 liter with distilled water.
- E. Standardization titration procedure : Place 10.0 ml standard calcium solution in vessel containing about 50 ml distilled water . Add 1 ml buffer solution. Add 1-2 drops indicator or small scoop of dry indicator. Titrate slowly with continuous stirring until the last reddish tinge disappears; adding last few drops at 3-5 second intervals. At end point the color is blue. Total titration duration should be 5 minutes from the time of buffer addition.

$$\text{Normality of of EDTA} = \frac{0.2}{\text{ml of EDTA}}$$

- (1) Aqueous ammonia (1 N): Dilute 70 ml of concentrated aqueous ammonia solution to 1 liter with distilled water.
- (2) Inhibitor: Sodium cyanide (extremely poisonous).

Procedure:

(1) Pretreatment:

- A. For drinking waters , surface waters , saline waters , and dilutions thereof , no pretreatment steps are necessary.
- B. For most wastewaters , and highly polluted waters, the sample must be digested as given in the Atomic Absorption Methods section of this manual.

(2) Titration of sample-normal to high hardness:

- A. Sample should require < 15 ml EDTA titrant, titration should be completed within 5 minutes after buffer addition.
- B. Place 25.0 ml sample in titration vessels, neutralize with 1 N aqueous ammonia and dilute to about 50 ml.
- C. Add 1 to 2 ml buffer solution to adjust to pH 10 ± 0.1 .
- D. If end point is not sharp (as determined by practice run) add about 250 mg of sodium cyanide inhibitor at this point.
- E. Add small scoop of dried powder indicator formulation.
- F. Titrate slowly with continuous stirring with standard EDTA titrant until last reddish tint disappears. Solution is normally blue at end point.

(3) Titration of sample of low hardness (less than 5 mg/L):

- A. Use a larger sample volume (100 ml).
- B. Use proportionately larger amounts of buffer, inhibitor and indicator.
- C. Use a microburette and run a blank using bidistilled, distilled or deionized water.

Calculations:

Hardness (EDTA) expressed as

$$\text{CaCO}_3(\text{mg/L}) = \frac{A \times N \times 50,000}{\text{ml sample}}$$

where :

A = volume (ml) of EDTA titrant .

N = normality of EDTA titrant.

Precautions:

Excessive amounts of heavy metals can interfere. This is usually overcome by complexing the metals with cyanide. Routine addition of sodium cyanide solution (*caution: deadly poison*) to prevent potential metallic interference is recommended.

References:

- (1) USEPA "Methods For Chemical Analysis of Water and Wastes", Ohio, 1983 (Method 130.2).
- (2) "Standard Methods for Examination of Water and Wastewater", 19th Edition, American Public Health Association, Washington 1995 (Method 2340 C).

Total Suspended Solids (Method A-8, Gravimetry)

Principle:

A well mixed sample is filtered through a glass fiber filter, and the residue retained on the filter is dried to constant weight at 103 - 105°C. The filtrate from this method may be used for dissolved solids. The practical range of the determination is 4 - 20000 mg/L.

Equipment:

- (1) Glass fiber filter discs, without organic binder.
- (2) Filter support: Filtering apparatus with reservoir and a coarse (40-60 microns) fritted disc as a filter support.
- (3) Suction flask.
- (4) Drying oven, 103 -105 °C.
- (5) Desiccator.
- (6) Analytical balance (0.1 mg sensitivity).

Reagents:

Deionized bidistilled water.

Procedure:

- (1) Place the glass fiber filter on the membrane filter apparatus or insert into bottom of a suitable Gooch crucible with wrinkled surface up.
- (2) While vacuum is applied, wash the disc with three successive 20 ml volumes of distilled water. all traces of water by continuing to apply vacuum after water has passed through.
- (3) Remove filter from membrane filter apparatus or both crucible and filter if Gooch crucible is used, and dry in an oven at 103 -105 °C for one hour.
- (4) Transfer to desiccator and store until needed.
- (5) Repeat the drying cycle until a constant weight is obtained (weight loss is less than 0.5 mg) . Weigh immediately before use . After weighing , handle the filter ocrucible / filter with forceps or tongs only.
- (6) For a 4.7 cm diameter filter, filter 100 ml sample. If weight of captured residue is less than 1.0 mg, the sample volume must be increased to provide at least 1.0 mg of residue. If other filter diameters are used, start with a sample volume equal to 7 ml/cm³ of filter area and collect at least a weight of residue proportional to the 1.0 mg stated above.
- (7) Assemble the filtering apparatus and begin suction. Wet the filter with a small volume of distilled water to seat it against the fritted support.
- (8) Shake the sample vigorously and quantitatively transfer the sample volume to the filter using a graduated cylinder. Remove all traces of water by continuing to apply vacuum after sample has passed through.
- (9) With suction on, wash the graduated cylinder, filter, non-filterable residue and filter funnel wall with three portions of distilled water allowing complete drainage between washing. Remove all traces of water by continuing to apply vacuum after water has passed through.
- (10) Carefully remove the filter from the filter support. Alternatively, remove crucible and filter from crucible adapter. Dry at least one hour at 103 - 105°C. Cool in a desiccator and weigh. Repeat the drying cycle until a constant weight is obtained (weight loss is less than 0.5 mg).

Calculations:

$$\text{Total suspended solids, (mg/L)} = \frac{(A-B) \times 1,000}{C}$$

where :

A = weight (mg) of filter (or filter and crucible) + residue.

B = weight (mg) of filter (or filter and crucible).

C = volume (ml) of sample filtered.

Precautions:

- (1) Non-representative particulates such as leaves , sticks , fish and lumps of fecal matter should be excluded from the sample if it is determined that their inclusion is not desired in the final result.
- (2) Preservation of the sample is not practical; analysis should begin as soon as possible.
- (3) Refrigeration or icing to 4°C, to minimize micro- biological decomposition of solids, is recommended.
- (4) Filtration apparatus, filter material, pre-washing, post-washing, and drying temperature are specified because these variables have been shown to affect the results.
- (5) Samples high in filterable residue (dissolved solids), such as saline waters, brines and some wastes, may be subject to a positive interference. Care must be taken in selecting the filtering apparatus so that washing of the filter and any dissolved solids in the filter minimizes this potential interference.

References:

- (1) USEPA “Methods For Chemical Analysis of Water and Wastes”, Ohio, 1983 (Method 160.2).
- (2) “ Standard Methods for Examination of Water and Wastewater”, 19th Edition, American Public Health Association, Washington 1995 (Method 2540 D).

Fixed and Volatile Solids (Method A-11, Gravimetry)

Principle:

The residue from the experiments: total solids, dissolved solids and suspended solids is ignited to constant weight at 550°C. The remaining solids represent the fixed total, dissolved or suspended solids while the weight lost on ignition is the volatile solids. The determination is useful in control of wastewater treatment plant operation because it offers a rough approximation of the amount of organic matter present in the solid fraction of wastewater, activated sludge, and industrial wastes.

Equipment:

See total dissolved and suspended solids.

Procedure:

- (1) Warm up the muffle furnace to 550°C before inserting samples.
- (2) Ignite residual produced by the methods of total, dissolved or suspended solids to constant weight in a muffle furnace at a temperature of 550°C (usually, 15 - 20 min. ignition are required for 200 mg residue). However, more than one sample and/or heavier residues may overtax the furnace and necessitate longer ignition times.
- (3) Let dish or filter disk cool partially in air until most of the heat has been dissipated.
- (4) Transfer to a desiccator for final cooling in a dry atmosphere. Do not overload desiccator.
- (5) Weigh dish or disk as soon as it has cooled to balance temperature. Repeat cycle of igniting, cooling, desiccating, and weighing until a constant weight is obtained or until weight change is less than 4% or 0.5 mg, whichever is less.
- (6) Duplicate determinations should agree within 5% of their average.

Calculations:

$$\text{Volatile solids, (mg/L)} = \frac{(A - B) \times 1000}{\text{Sample volume, ml}}$$

$$\text{Fixed solids, (mg/L)} = \frac{(B - C) \times 1000}{\text{Sample volume, ml}}$$

where :

A = weight of residue + dish before ignition (mg).

B = weight of residue + dish or filter after ignition (mg).

C = weight of dish or filter (mg) .

Precautions:

- (1) Negative errors in the volatile solids may be produced by loss of volatile matter during drying .
- (2) Determination of low concentrations of volatile solids in the presence of high fixed solids concentrations may be subjected to considerable error. In such cases, measure for suspect volatile components by another test , for example , total organic carbon .

References:

- (1) USEPA “Methods For Chemical Analysis of Water and Wastes”, Ohio, 1983 (Method 160.4).
- (2) “Standard Methods for Examination of Water and Wastewater”, 19th Edition, American Public Health Association, Washington 1995 (Method 2540 D).

Turbidity (Method A-12, Nephelometry)

Principle:

Turbidity in water is caused by suspended matter and is an expression of the optical property that causes light to be scattered and absorbed rather than transmitted in straight lines through the sample. The most commonly used method is a nephelometric method. This method is based on a comparison of the amount of light scattered by a standard reference suspension under the same conditions. The more light is scattered the higher is the turbidity. The method is applicable in the turbidity range 0 - 40 NTU.

Equipment:

- (1) Turbidity meter.
- (2) Sample tubes.

Reagents:

(1) Stock turbidity suspension:

Solution (I):

Dissolve 1.00 g of hydrazine sulfate $(\text{H}_2\text{N})_2\text{H}_2\text{SO}_4$ in distilled water and dilute to 100 ml.

Solution (II):

Dissolve 10 g hexamethylene-tetramine, $(\text{CH}_2)_6\text{N}_4$, in distilled water and dilute to 100 ml.

In a 100 ml volumetric flask, mix 5.0 ml of solution (I) and 5.0 of solution (II). Allow to stand for 24 hrs. Dilute to mark and mix. The turbidity of this suspension is 400 NTU (Nephelometric Turbidity Units). These suspensions can be stored for a month.

(2) Standard turbidity suspension:

Dilute 25 ml of stock turbidity suspension to 100 ml with distilled water. Prepare daily. This suspension has a turbidity equivalent to 100 NTU. Prepare another dilution by diluting the 100 NTU suspension 10 times to give a 10 NTU suspension.

Procedure:

- (1) Calibrate the nephelometer with the 100 NTU suspension and verify the setting with a reading on the 10 NTU suspension.
- (2) Thoroughly shake sample, wait until air bubbles disappear and pour sample into turbidity meter tube. Read turbidity directly from the instrument scale.
- (3) In case the turbidity is going beyond the scale dilute the sample and measure the turbidity. The dilution factor should be taken into account for the final calculation of the result.

Calculations:

$$\text{Sample NTU} = \frac{\text{NTU measured}}{\text{proportion of sample in dilution}}$$

Precautions:

Rapidly settling sand, the presence of air bubbles, and dirty glassware are the main sources of interference and should be removed before analysis.

Referenc:

1. USEPA "Methods For Chemical Analysis of Water and Wastes", Ohio, 1983 (Method 180.1).
2. "Standard Methods for Examination of Water and Wastewater", 19th Edition, American Public Health Association, Washington 1995 (Method 2130 B).

(B) Determination of Inorganic Nonmetallic Constituents

Quality Control:

1. Detailed recommendations and general information concerning laboratory quality assurance, as well as specific details of applicable quality-control measures, should be followed.
2. Keep in mind the overall purpose of the measurements. Can duplicates (replicates) adequately establish precision and can recovery of known additions determine bias. Are suitable standards, control charts blanks, calibrations, and other necessary measures incorporated to assure the quality of data produced. Furthermore, can adequate documentation be presented to support these contentions.

The same degree of substantiation may not be necessary in all cases, e.g., requirements for an operational control determination may differ from those for data submitted for regulatory purposes. Therefore, the intended use of the data is an important factor in determining overall quality control as well as the accompanying quality assurance measures.

Chlorine, Total Residue (Method B-3, Spectrophotometry)

Principle:

Chlorine (hypochlorite ion, hypochlorous acid) and chloramines stoichiometrically liberate iodine from potassium iodide at pH 4 or less. The liberated iodine reacts with N,N-diethyl-*p*-phenylenediamine (DPD) to produce a red colored solution with maximum absorbance at 515 nm. The method is applicable to the concentration range of 0.2-4 mg/L chlorine.

Equipment:

Spectrophotometer for use at 515 nm and cells of light path 1 cm or longer.

Reagents:

- (1) Phosphate buffer solution: Dissolve 24 g anhydrous disodium hydrogen phosphate, Na_2HPO_4 , and 46 g anhydrous potassium dihydrogen phosphate, KH_2PO_4 , in distilled water. Dissolve 800 mg disodium ethylenediamine tetraacetate dihydrate in 100 ml distilled water. Combine these two solutions and dilute to 1 liter with distilled water. Add 20 mg HgCl_2 as a preservative.
- (2) N,N-Diethyl-*p*-phenylenediamine (DPD) indicator solution : Dissolve 1 g DPD oxalate or 1.5 g *p*-amino-N,N-diethylaniline sulfate in chlorine free distilled water containing 8 ml of 1 + 3 H_2SO_4 and 200 mg disodium ethylenediamine tetraacetate dihydrate. Dilute to 1 liter, store in a colored, glass-stoppered bottle. Discard when discolored. The buffer and indicator sulfate are available as a combined reagent in stable powder form.
- (3) Sulfuric acid solution (1+3) : Slowly add one part of H_2SO_4 (sp.gr.1.84) to three parts of distilled water.
- (4) Potassium iodide, KI crystals.
- (5) Stock potassium permanganate solution : Place 0.891 g KMnO_4 in a volumetric flask and dilute to 1 liter.
- (6) Standard potassium permanganate solution: Dilute 10.00 ml of stock potassium permanganate solution to 100 ml with distilled water in a volumetric flask. One milliliter of this solution diluted to 100 ml with distilled water is equivalent to 1.00 mg/L chlorine.

Procedure:

- (1) Calibration:
 - A. Prepare a series of permanganate standards covering the chlorine equivalent range of 0.05 to 4 mg/L.
 - B. Place 5 ml phosphate buffer in a flask.
 - C. Add 5 ml DPD reagent.
 - D. Add 100 ml permanganate standard.
 - E. Read at 515 nm on a spectrophotometer and record the absorbance.
 - F. Return the contents of the cell to the flask.
 - G. Titrate the contents of the flask with standard ferrous ammonium sulfate until the red color is discharged. Record the result.
- (2) Sample analysis:
 - A. Place 0.5 ml phosphate buffer in flask.
 - B. Add 0.5 ml DPD reagent.

- C. Add approximately 0.1 g KI.
- D. Add 10 ml of sample.
- E. Let stand 2 minutes.
- F. Read at 515 nm spectrophotometer , and record the absorbance.

Calculations:

- (1) Plot the absorbance of the standard permanganate solutions on the vertical axis versus the titrated concentration on the horizontal axis.
- (2) Draw the line of best fit through the points.
- (3) Locate the absorbance of the sample on the vertical axis.
- (4) Read the concentration on the horizontal axis at the intersect of the absorbance and the calibration line .

Precautions:

- (1) Interferences are caused by any oxidizing agent but these are usually present at insignificant concentrations compared to the residual chlorine concentrations .
- (2) Turbidity and color will essentially prevent the colorimetric analysis .

References:

- (1) USEPA “Methods For Chemical Analysis of Water and Wastes”, Ohio, 1983 (Method 330.5).
- (2) “Standard Methods for Examination of Water and Wastewater”, 19th Edition, American Public Health Association, Washington 1995 (Method [4500-Cl (C)]).

Ammonia-Nitrogen (Method B-6, Potentiometry)

Principle:

The ammonia-selective electrode uses a hydrophobic gas-permeable membrane to separate the sample solution from an electrode internal solution of ammonium chloride. Dissolved ammonia is converted to ammonia gas by raising the pH to above 11 with a strong base. Diffusion of ammonia through the membrane causes change of the pH of the internal solution that is sensed by a pH electrode inside the ammonia sensor. Sample distillation is unnecessary. The method covers the range from 0.03 to 1400 mg NH₃/L.

Equipment:

- (1) Electrometer: A pH meter with expanded millivolt scale capable of 0.1 mV resolution between - 700 mV and + 700 mV or a specific ion meter.
- (2) Ammonia-selective electrode.
- (3) Magnetic stirrer, thermally insulated, with TFE-coated stirring bar.

Reagents:

- (1) Ammonia free water.
- (2) Sodium hydroxide, 10 N.
- (3) NaOH/EDTA solution, 10 N : Dissolve 400 g NaOH in 800 ml water. Add 45.2 ethylenediaminetetraacetic acid, tetrasodium salt, tetrahydrate (Na₄EDTA.4H₂O) and stir to dissolve. Cool and dilute to 1000 ml.
- (4) Stock ammonium chloride solution: Dissolve 3.819 g anhydrous NH₄Cl (dried at 100°C) in water, and dilute to 1000 ml ; 1.00 ml = 1.00 mg N = 1.22 mg NH₃.
- (5) Standard ammonium chloride solutions.

Procedure:

- (1) Preparation of standards: Prepare a series of standard solutions covering the concentration of 1000 , 100, 10, 1, and 0.1 mg NH₃-N/L by making decimal dilutions of stock NH₄Cl solution with water.
- (2) Electrometer calibration: Place 100 ml of each standard solution in a 150 ml beaker. Immerse electrode in standard of lowest concentration and mix with a magnetic stirrer. Limit stirring speed to minimize possible loss of ammonia from the solution. Maintain the same stirring rate and a temperature of about 25°C throughout calibration and testing procedures. Add a sufficient volume of 10 N NaOH solution (1 ml is usually sufficient) to raise pH above 11. If the presence of silver or mercury is possible, use NaOH / EDTA solution in place of NaOH solution. If it is necessary to add more than 1 ml of either NaOH or NaOH / EDTA solution, note volume used, because it is required for subsequent calculations. Keep electrode in solution until stable millivolt reading is obtained. Do not add NaOH solution before immersing electrode, because ammonia may be lost from a basic solution repeat procedures with remaining standards, proceeding from lowest to highest concentration. Wait until the reading has stabilized (at least 2 to 3 minutes) before recording millivolts for standards and samples containing less than 1 mg NH₃-N/L.
- (3) Preparation of standard curve: Using semi-logarithmic graph paper, plot ammonia concentration in milligrams NH₃-N per liter on the log axis vs. potential in millivolts on the linear axis starting with

the lowest concentration at the bottom of the scale. If the electrode is functioning properly a tenfold change of $\text{NH}_3\text{-N}$ concentration produces a potential change of about 59 mV.

- (4) Calibration of specific ion meter: Refer to manufacturer's instructions.
- (5) Measurement of samples: Dilute if necessary to bring $\text{NH}_3\text{-N}$ concentration to within calibration curve range. Place 100 ml sample in 150 ml beaker and follow procedure in the electrometer calibration above. Record volume of 10 N NaOH added. Read $\text{NH}_3\text{-N}$ concentration from standard curve.

Calculations:

$$\text{NH}_3\text{-N, (mg/L)} = A \times B \times \frac{100 + D}{100 + C}$$

Where:

A = dilution factor.

B = concentration of $\text{NH}_3\text{-N}$ (mg/L) from calibration curve.

C = volume (ml) of 10 N NaOH added to calibration standards.

D = volume (ml) of 10 N NaOH added to sample.

Precaut:

- (1) Amines are a positive interference.
- (2) Mercury and silver interfere by complexing with ammonia, unless the NaOH /EDTA solution is used.

References:

- (1) USEPA "Methods For Chemical Analysis of Water and Wastes", Ohio, 1983 (Method 350.3).
- (2) "Standard Methods for Examination of Water and Wastewater", 19th Edition, American Public Health Association, Washington 1995 (Method [4500- $\text{NH}_3\text{(D)}$]).

Nitrite-Nitrogen (NO₂-N) (Method B-7, Spectrophotometry)

Principle:

Nitrite reacts in strongly acid medium with sulphanilamide. The resulting diazo compound is coupled with N-(1-naphthyl)-ethylenediamine-dihydrochloride to form an intensely red-coloured azo-compound. The absorption of the dye is proportional to the concentration of nitrite present. The method is applicable in the range of 0.01 to 1.0 mg nitrite-nitrogen/L. Samples containing higher concentration must be diluted before analysis.

Equipment:

- (1) Spectrophotometer.
- (2) Nessler tubes (50 ml) and 1.0 cm glass cuvettes.

Reagents:

1. Buffer-colour reagent: Put 105 ml of HCl in glass stoppered dark bottle. Add 5.0 g of sulphanilamide and 0.5 g of N(1-naphthyl)-ethylenediamine dihydrochloride. Stir until dissolved. Add 136 g of sodium acetate and again stir until dissolved. Dilute to 500 ml with distilled water. This solution is stable for several weeks in the dark (wrap in aluminum foil).
2. Stock nitrite solution: Dissolve 0.4926 g of dried anhydrous sodium nitrite in distilled water and dilute to 1 liter. Preserve with 2 ml chloroform (1.0 ml = 100 µg nitrite-N).
3. Working standard solution: Dilute 10.0 ml of the stock solution to 1 L of distilled water (1.0 ml = 1.0 µg nitrite-N).

Procedure:

- (1) Filter the sample, if necessary, through a 0.45 µm pore-size membrane filter.
- (2) Put 10-50 ml of the sample in a 50 ml Nessler tube.
- (3) Prepare a series of standards in 50 ml Nessler tubes (1.0-20.0 µg of nitrite-N per tube).
- (4) Add 2 ml of buffer color reagent to each standard and sample, mix and allow color to develop for at least 15 min.
- (5) Check the pH of the solution to ensure that the value is between pH 1.5 and 2.0.
- (6) Measure the absorbance of standard and sample in the spectrophotometer at 540 nm against a blank.
- (7) Plot the amount of nitrite-nitrogen per standard tube against the absorbance.
- (8) Read the concentration of nitrite-nitrogen per sample tube from the calibration curve.

Calculations:

Nitrite-nitrogen, (µg/L) = 1000 x reading graph (ug nitrite-nitrogen)

volume of sample (ml)

Precautions:

- (1) The nitrite concentration in collected samples can change very rapidly due to bacterial oxidative or reductive conversions and the analysis should be carried out as early as possible.
- (2) Short time preservation for one or two days is possible by deep freezing (-20°C) or the addition of 40 mg mercuric ion/ liter of sample with storage at 4°C .

References:

1. USEPA “Methods For Chemical Analysis of Water and Wastes”, Ohio, 1983 (Method 354.1).
2. “Standard Methods for Examination of Water and Wastewater”, 19th Edition, American Public Health Association, Washington 1995 (Method [4500- NO_2^- (B)]).

(Nitrate and Nitrite)-Nitrogen (Method B-8, Spectrophotometry)

Principle :

Nitrate is reduced to nitrite when a sample is passed over a column containing copper amalgamated cadmium granules or filings. The nitrite if originally present , plus reduced nitrate) is determined by the same method as described earlier for nitrite . The applicable range of this method is 0.01 to 1.0 mg/L nitrate-nitrite nitrogen . The range can be extended by sample dilution .

Equipment:

- (1) Reduction columns: Constructed from three pieces of glass tubing joined end to end (as shown in the following figure 1).
- (2) All glass ware needed for nitrite analysis.
- (3) Spectrophotometer with 1.0 cm glass cuvettes.

Reagents:

- (1) Prepare copper - cadmium granules by washing 25 g new or used 40 - 60 mesh cadmium metal (Cd) granules with 6 N HCl and rinse with water . Swirl Cd with 100 ml 2 % CuSO₄ solution for 5 minutes or until the blue color partially fades . Decant and repeat with fresh CuSO₄ until a brown colloidal precipitate begins to develop. Gently flush with water to remove all precipitated copper metal (Cu).
- (2) Buffer - color reagent : Put 105 ml of hydrochloric acid in glass stoppered dark bottle. Add 5.0 g sulphanilamide and 0.5 g of N- (1 -naphthyl) ethylenediamine dihydrochloride. Stir until dissolved . Add 136 g of sodium acetate and again stir until dissolved. Dilute to 500 ml with distilled water. This solution is stable for several weeks if stored in the dark (wrap in aluminum foil).
- (3) Ammonium chloride - EDTA solution: Dissolve 13 g of NH₄Cl and 1.7 g of (EDTA) in 900 ml distilled water. Adjust the pH to 8.5 with concentrated aqueous ammonia and dilute to one litre.
- (4) Diluted ammonium chloride-EDTA solution: Dilute 300 ml NH₄Cl-EDTA solution to 500 ml with water.
- (5) Hydrochloric acid: HCl , 6 N.
- (6) Copper sulphate solution, 2%.
- (7) Stock nitrite solution: Dissolve 0.4926 g of dried anhydrous sodium nitrite (24 hrs. In desiccator) in distilled water and dilute to one litre. Preserve with 2 ml chloroform per litre(1.0 ml = 100 µg nitrite-N).
- (8) Test nitrite solution: Working standard solution is prepared by diluting: Dilute 10.0 ml of the stock solution to 1 litre of distilled water (1.0 ml = 1.0 µg nitrite - N).
- (9) Stock nitrate solution : Dissolve 7.218 g of KNO₃ in distilled water, and dilute to 1 litre . Addition of 2 ml of chloroform allows for storage up to 6 months (1.0 ml = 1.0 mg NO₃-N).
- (10) Standard nitrate solution: Dilute 10 ml of stock nitrate solution to 1 litre (1.0 ml = 10 µg NO₃-N)

Procedure:

- (1) Preparation of reduction column : Plug the bottom of the reduction column and fill with water. Add sufficient Cu-Cd granules to produce a column of 18.5 cm long. Maintain the water level above Cu-Cd

- granules to prevent entrapment of air. Wash the column with 200 ml of diluted NH_4Cl -EDTA solution. Activate the column by passing through it, at least 100 ml of a solution composed of 25 % 1.0 mg $\text{NO}_3\text{-N/L}$, at 7 to 10 ml/min.
- (2) Filter the sample through a 0.45 μm pore size membrane filter.
 - (3) Adjust pH to between 7 and 9, using NaOH or HCl as necessary. This ensures a pH of 8.5 after adding NH_4Cl -EDTA solution.
 - (4) To 25 ml sample or a portion diluted to 25 ml, add 75 ml NH_4Cl -EDTA solution and mix. Pour mixed sample into column and collect at a rate of 7 to 10 ml/min. Discard first 25 ml . Collect the rest in the original sample flask. There is no need to wash columns between samples , but if the columns are not to be used for several hours or longer , pour 50 ml dilute NH_4Cl -EDTA solution on , and let it pass through the system . Store Cu-Cd column in this solution and never let it dry.
 - (5) Carry out the reduction of $\text{NO}_3\text{-N}$ standards exactly as prescribed for the samples. At least one nitrite standard should be compared to a reduced nitrate standard at the same concentration to verify the efficiency of the reduction column.
 - (6) As soon as possible , and not more than 15 min after reduction, add 2.0 ml color reagent to 50 ml of sample or standard and mix. Between 10 min and 2 hours afterward, measure absorbance at 543 nm against a distilled water reagent blank. Note: if the $(\text{NO}_3+\text{NO}_2)\text{-N}$ concentration exceeds the standard curve range (about 1 mg N/L), use the remainder of reduced sample to make an appropriate dilution and analyse again.
 - (7) Make a calibration curve relating the absorbance of the nitrate standards against the amount of NO_3+NO_2 present . Read the amount of NO_3+NO_2 in a sample from the graph.

Calculations:

$$(\text{NO}_3 + \text{NO}_2) - \text{N} , (\text{mg/L}) = \frac{\text{reading from graph } (\mu\text{g})}{\text{volume. of sample (ml)}}$$

Precautions:

- (1) To prevent any change in the nitrate it should be determined immediately after sampling. If storage is necessary, the sample should be kept at 4°C with or without preservatives such as H_2SO_4 (0.8 ml of $\text{coH}_2\text{SO}_4\text{/L}$ of sample).
- (2) Build up of suspended matter in the reduction column will restrict sample flow. Since nitrate nitrogen occurs in a soluble state, the sample can be prefiltered.
- (3) Highly turbid sashould be treated with zinc sulfate before filtration to remove turbidity.
- (4) Low results might be obtained for samples that contain high concentrations of iron, copper, cadmium or other metals. EDTA is added to eliminate this interference.

References:

- (1) USEPA “Methods For Chemical Analysis of Water and Wastes”, Ohio, 1983 (Method 353.3).
- (2) “ Standard Methods for Examination of Water and Wastewater”, 19th Edition, American Public Health Association, Washington 1995 (Method [4500- $\text{NO}_3(\text{E})$]).

Phosphorus (Orthophosphate) (Method B-16, Spectrophotometry)

Principle:

Orthophosphate reacts under acid conditions to form a heteropoly acid, molybdophosphoric acid. In the presence of vanadium, yellow vanadomolybdophosphoric acid is formed. The intensity of the yellow color is proportional to phosphate concentration. The minimum detectable concentration is 0.2 mg phosphorus/L in 1 cm spectrophotometric cell.

Equipment:

- (1) Spectrophotometer for use at 400-490 nm.
- (2) Acid-washed glassware.
- (3) Filtration apparatus and filter paper.

Reagents:

- (1) Phenolphthalein indicator.
- (2) Hydrochloric acid (1:1).
- (3) Activated carbon.
- (4) Vanadate-molybdate reagent :
Solution (I) : Dissolve 25 g ammonium molybdate $(\text{NH}_4)_6\text{Mo}_7\text{O}_{24} \cdot 4\text{H}_2\text{O}$ in 300 ml distilled water.

Solution (II): Dissolve 1.25 g ammonium metavanadate NH_4VO_3 by heating to boiling in 300 ml distilled water .
Cool and add 330 ml of concentrated HCl .

Cool solution (II) to room temperature , pour solution (I) into solution (II) , mix and dilute to 1 liter.
- (5) Standard phosphate solution : Dissolve in distilled water 219.5 mg anhydrous KH_2PO_4 and dilute to 1000 ml (1.00 ml = 50 μg $\text{PO}_4\text{-P}$).

Procedure:

- (1) If sample pH is greater than 10 , add 0.05 ml (1 drop) phenolphthalein indicator to 50 ml sample and discharge the red color with 1:1 HCl before diluting to 100 ml .
- (2) Remove excessive color in sample by shaking about 50 ml with 200 mg activated carbon in an erlenmeyer flask for 5 min. and filter to remove carbon.
- (3) Place 35 ml or less of sample containing 0.05-1.0 mg P , in a 50 ml volumetric flask. Add 10 ml vanadate-molybdate reagent and dilute to the mark with distilled water.
- (4) Prepare a blank in which 35 ml distilled water is substituted for the sample.
- (5) After 10 min., measure absorbance of sample versus a blank at a wavelength of 400-490 nm. and compare with the calibration curve.
- (6) Prepare a calibration curve by using 1-20 µg/ml standard phosphate solution and proceed as above.

Calculations:

$$\text{Phosphate-P, (mg/L)} = \frac{\text{mg P (in 50 ml final volume)} \times 1000}{\text{ml sample}}$$

Precautions:

- (1) Positive interference is caused by silica and arsenate only if the sample is heated.
- (2) Negative interferences are caused by arsenate, fluoride, thorium, bismuth, sulfide, thiosulfate, thiocyanate, or excess molybdate.
- (3) Blue color is caused by ferrous iron concentration but this does not affect results if ferrous concentration is less than 100 mg/L.
- (4) Sulfide interference may be removed by oxidation with bromine water.

References:

“Standard Methods for Examination of Water and Wastewater”, 19th Edition, American Public Health Association, Washington 1995 (Method [4500-P (C)]).

Total Coliform (Method E-1 Microbiology)

Principle:

This bacteria develops red colonies with a golden metallic sheen within 24 hours at 35°C on an Endo-type medium containing lactose. All red, pink, blue, white or colorless colonies lacking sheen are considered non-coliforms.

Equipment:

- (1) Sterilized glass or polypropylene wide-mouthed bottles (300 ml to 500 ml).
- (2) Sterilized pipettes and graduated cylinders
- (3) Sterilized glass Erlenmeyer flasks with metal caps, metal foil covers or screw caps.
- (4) Sterilized borosilicate glass or disposable plastic petri dishes (60x15 mm, 50 x 12 mm).
- (5) Filtration unit consisting of a membrane filter funnel and membrane holder. (The unit should be made of seamless steel, glass, or autoclavable plastic that does not leak and is not subject to corrosion.
- (6) Membrane filter (pre-sterilized and certified membranes).
- (7) Sterile absorbent pads 48 mm in diameter and of sufficient thickness to absorb 1.8 to 2.2 ml of medium.
- (8) Sterilized smooth-tipped forceps without corrugations on the inner side of the tips (Ethyl or absolute methyl alcohol [95%] and flaming may be used for sterilization).
- (9) Incubator adjusted to 35±0.5 °C and a 90% relative humidity.
- (10) Microscope with a 10 to 15 magnification and a white fluorescent light source

Reagents and Materials

- (1) LES Endo agar : Dehydrated Difco M-Endo Agar LES (No. 0736), dehydrated BBL M-Endo Agar LES (No. 11203), or equivalent.
- (2) Dehydrate in IL water containing 20 ml 95% ethanol. Bring to a near boil to dissolve to agar, then remove from heat to cool to 45-50 °C. Final pH should be 7.2±0.2. Use 5 to 7 ml for each 60 mm petri dish. Do not expose plates to sunlight; store in the dark at 4 to 8 °C. Discard the medium after 2 weeks.
- (3) M-Endo medium : Dehydrated Difco M-Endo Broth Mf (No. 0749), dehydrated BBL M-coliform broth (No. 11119), or equivalent.
- (4) Dehydrate in IL water containing 20 ml 95% ethanol. Bring to a near boil to dissolve agar, then remove from heat to cool to 45-50 °C. Final pH should be 7.0 to 7.3. Use ml for each 60 mm petri dish. Do not expose plates to sunlight, store in the dark at 4 to 8 °C. Discard broth medium after 96 hours and discharge agar after 3 weeks.

- (5) Liquid medium (without agar) can be used with absorbent pads if the pads are certified to be free of sulfite or any other agent that could inhibit bacterial growth. Use 2 ml of the liquid medium per plate.

Procedure:

- (1) Step 1 : Sample Size
- (2) For drinking water, filter 100 to 1000ml volume which yields about 50 coliform colonies through a membrane filter (MF).
- (3) For other waters, filter three different volumes (diluted or undiluted), depending on the expected bacterial density.
- (4) Step 2: Sample Filtration
- (5) Place a sterile MF (grid side up) over porous plate of receptacle using sterile forceps.
- (6) Place matched funnel unit over receptacle carefully, lock it in place, and filter the sample under a partial vacuum.
- (7) Rinse funnel by filtering a 20-to-30 ml portion of sterile diluted water. Disengage vacuum, unlock and remove funnel, and remove MF using sterile forceps. Place MF on selected medium with a rolling motion to avoid entrapment of air.
- (8) Prepare a sterile rinse water sample (100 ml) as above for every 10 sample series.
- (9) At the beginning of each filtration series, use a sterile filtration unit. If a time interval of more than 30 minutes elapses between sample filtrations, begin a new series using a sterile filtration unit. Sterilize unit with flowing stream of boiling water or ultraviolet disaffection (UV).
- (10) Step 3: Culture Dish and Incubation
- (11) If agar-based Endo-type medium is used: Place prepared filter directly on agar by rolling it onto agar surface without entrapping air.
- (12) If liquid medium is used : Place a pad in the culture dish and saturate it with 1.8 to 2.0 ml M-Endo medium. Place prepared filter directly on pad. Invert dish.
- (13) Incubate the prepared dishes for 22 to 24 h at $35 \pm 0.5^{\circ}\text{C}$.
- (14) Step 4: Counting
- (15) Use a low-power (10 to 15 magnification) binocular dissecting microscope or other optical device with a white fluorescent light. Count the number of colonies that are pink-to-dark-red colored with a metallic surface sheen.
- (16) Atypical coliform colonies can be dark-red or nucleated without sheen. Pink, red, white and colorless colonies lacking sheen are considered to be non-coliforms.

Calculations:

$$TC = \frac{CCC \times 100}{V_s}$$

Where:

TC = total coliform colonies/100 ml.

CCC = counted coliform colonies.

V_s = volume of sample filtered (ml).

Precautions:

- (1) The presence of algae, particulate and turbidity in the sample solution can cause interference.
- (2) This procedure is not applicable to wastewater samples that have received primary treatment followed by chlorination.
- (3) Thoroughly clean and sterilize all equipment.

References:

- (1) "Standard Methods for Examination of Water and Wastewater", 19th Edition, American Public Health Association, Washington 1995 (Method 9222 B).

Fecal Coliform (Method E-2, Microbiology)

Principle:

Fecal coliform causes fermentation of lactose and grows on M-FC broth within 48h at 44-45 °C. It forms blue-colored colonies. Non-fecal coliform colonies are gray to cream colored.

Equipment :

- (1) Sterilized glass or polypropylene, wide-mouthed bottles (300 ml to 500 ml).
- (2) Sterilized pipettes and graduated cylinders.
- (3) Sterilized glass Erlenmeyer flasks with metal caps, metal foil covers or screw caps.
- (4) Sterilized borosilicate glass or disposable plastic petri dishes (60x15 mm, 50x12 mm).
- (5) Filtration unit consisting of a membrane filter funnel and membrane holder. (The unit should be made of seamless stainless steel, glass, or autoclavable plastic that does leak and is not subject to corrosion.)
- (6) Membrane filter (pre-sterilized and certified membranes).
- (7) Sterile absorbent pads 48 mm in diameter and of sufficient thickness to absorb 1.8 to 2.0 ml of medium.
- (8) Sterilized smooth-tipped forceps without corrugations on the inner side of the tips (ethyl or absolute methyl alcohol (95%) and flaming may be used for sterilization).
- (9) Incubator adjusted to 44.5±0.5 °C or a suitable water bath.
- (10) Microscope with magnification of 10 to 15 times and white fluorescent light source.

Reagents and Materials:

- (1) Rehydrate M-FC medium in water containing 10 ml 1% rosolic acid in 0.2 N NaOH. Heat to near boiling, remove from heat and cool to below 50 °C. Do not sterilize by autoclaving. The final pH should be 7.4. Store the finished medium at 4 to 8 °C in sealed plastic bags or other suitable container to reduce moisture loss. Discard unused broth after 96 h and the unused agar medium after 2 weeks.
- (2) Use 5 to 7 ml of the agar medium per 50x12 mm petri dish and let solidify.

Procedure:

Step 1: Sample Size

For drinking water, a 50 to 100 ml sample is suitable. When the bacterial density of the sample is unknown, filter a few milliliters to establish fecal coliform density. Estimate the volume expected to yield a countable membrane. Two other portions, representing both one tenth and ten times the estimated volume, can be used.

Step 2: Sample Filtration

- 1) Place a sterile MF (grid side up) over the porous plate of receptacle using sterile forceps.
- 2) Carefully place matched funnel unit over receptacle and lock it in place. Filter the sample under a partial vacuum.
- 3) Rinse funnel by filtering a 20 to 30 ml portion of sterile dilution water. Disengage vacuum, unlock and remove funnel.
- 4) Remove MF using sterile forceps and place it on selected medium with rolling motion to avoid entrapment of air.
- 5) Prepare a sterile rinse water sample (100 ml) as above after every 10 samples in a series.

At the beginning of each filtration series, use a sterile filtration unit. If a time interval of more than 30 minutes elapses between sample filtrations, begin a new series using a sterile filtration unit. Sterilize unit with a flowing stream of boiling water or UV sterilizer.

Step 3 : Preparation of Culture Dish

If using liquid M-FC broth with pads:

- (1) Place a sterile absorbent pad in each culture dish
- (2) Saturate pad with about 2.0 ml M-FC medium.
- (3) Remove any excess liquid from culture dish.
- (4) Place MF the pad and invert the dish

If using M-FC broth without pads:

- (1) Add 1.5% agar to the M-FC broth.
- (2) Place 5 to 7 ml of the medium in each culture dish and let solidify.
- (3) Roll the prepared MF onto agar surface, being careful to avoid entrapment of air.

Step 4: Incubation:

- (1) Place prepared dishes in waterproof plastic bags or sealed petri dishes with waterproof tape to prevent leakage during submersion.
- (2) Submerge in water bath and incubate for 24 ± 2 h at $44.5 \pm 0.2^\circ \text{C}$.

Step 5: Counting

- (1) Use wide-field dissecting microscope (low power, 10-15 magnification) or other optical device to count the blue colored colonies of E. coli on the M-FC.

- (2) Do not count the gray to cream colored coliform colonies.

Calculations :

$$FC = \frac{CCC \times 100_s}{V_s}$$

Where:

FC = fecal coliforms /100 ml
 CCC = fecal coliform colonies counted
 V_s = volume of sample filtered.

Precautions

- (1) The presence of algae, particulates and turbidity in the sample solution can cause interference.
- (2) This procedure is not applicable to wastewater samples that have received primary treatment followed by chlorination
- (3) Be sure to clean and sterilize of all equipment

References:

“Standard Methods for Examination of Water and Wastewater”, 19th Edition, American Public Health Association, Washington 1995 (Method 9222 D).